## Quantitative determination of Haptoglobin (HAPTO) IVD

Store 2 - 8°C.

### PRINCIPLE OF THE METHOD

Anti-haptoglobin antibodies when mixed with samples containing haptoglobin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the haptoglobin concentration of the patient sample, that can be quantified by comparison from a calibrator of known haptoglobin concentration.

### CLINICAL SIGNIFICANCE

The haptoglobin is an  $\alpha_2$ -glycoprotein synthesized in the liver that binds hemoglobin irreversibly. The hapto-hemoglobin complexes, as well as free haptoglobin itself, play significant roles in the iron storage and prevention of possible renal damage as a consequence of hemoglobin excretion. As an acute-phase protein, haptoglobin is increased in the presence of acute inflammatory process, tissue necrosis or malignancy.

Haptoglobin deficiency in plasma is a consequence of hemolysis "in vivo", presence of estrogens in pregnancy and oral contraceptive therapy, as well as most forms of acute or chronic hepatocellular disease, including viral hepatitis.

Haptoglobin test is mainly used for the determination and monitoring of the hemolytic disorders. Under normal circumstances, approximately 1% of circulating red blood cells are destroyed from circulation. An increase of only 2% destruction per day will completely deplete plasma haptoglobin in the absence of a stimulus to production, such as acute inflammation or corticosteroids therapy.

### REAGENTS

<b>Diluent (R1)</b>	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Preservative.
<b>Antiserum (R2)</b>	Goat serum, anti-human haptoglobin pH 7.5. Preservative.
<b>Optional</b>	Ref: 1102003 PROT CAL.

### CALIBRATION

The assay has been standardized against the Reference Material ERM-DA470k/IFCC. It must be used with the PROT CAL Calibrator to calibrate the reagent.

### PREPARATION

**Reagents:** Ready to use.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

**Reagent deterioration:** The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

### ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus autoanalyzer
- Laboratory equipment.

### SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

### REFERENCE VALUES<sup>2</sup>

Between 30 - 200 mg/dL. Each laboratory should establish its own reference range

### QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

### SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME		
Nº	**	Vol. R1	240	
Test	HAPTO	Vol. R2	60	
Full Name	Haptoglobin	Vol. R3		
Standard nº	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.1 Slope 1	
Vol. Sample Indec.		Unit	mg/dL Inter. 0	
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	340	Reagent Blank	20-21	
Sec. Wave.		React. Time	56-57	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
Nº CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.1	10.0	90.0	3.0
1	CAL *0.25	30.0	90.0	3.0
2	CAL *0.4951	30.0	90.0	6.0
3	CAL *0.7354	30.0	90.0	9.0
4	CAL	0.0	0.0	3.0

### SPIN640Plus APPLICATION

EDIT PARAMETERS						
Test	HAPTO	No.	**			
Full name	HAPTOGLOBIN	Print name	HAPTO			
Reac. Type	End Point	Direction	Increase			
Pri. Wave.	340	Sec. Wave.				
Unit	mg/dL	Decimal	0.1			
Reagent Blank	20-21	React. Time	63 - 64			
Vol. Sample	3 ul	R1	240 ul			
Increased		R2	60 ul			
Decreased		R3				
Sample blank		R4				
CALIBRATION	Spline	Water	0	Vol aspirated	Vol net	Vol dil
Calibrator	*0.1			3	13	115
Calibrator	*0.25			3	35	105
Calibrator	*0.489			6	35	105
Calibrator	*0.735			9	35	105
Calibrator	(*)Cal Value			3	--	--

The Calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

### PERFORMANCE CHARACTERISTICS

**1. Measurement range:** Up to 300 mg/dL under the described assay conditions.

Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

**2. Detection Limit:** Values less than 1,3 mg/dL give non-reproducible results.

**3. Prozone effect:** No prozone effect was detected upon 1200 mg/dL.

**4. Sensitivity:**  $\Delta$  4,96 mA / mg/dL (100 mg/dL).

**5. Precision:** The reagent has been tested for 20 days, using three levels of serum in an EP5-based study.

EP5	CV (%)		
	39,25 mg/dL	97,35 mg/dL	191,5 mg/dL
Total	8%	3,2%	2,3%
Within Run	1,5%	0,9%	1,2%
Between Run	6,7%	2,3%	1,2%
Between Day	4%	2%	1,5%

**6. Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the method from Beckman (System Array 360 CE). 35 samples ranging from 10 to 400 mg/dL of Haptoglobin were assayed. The correlation coefficient (r) was 0,95 and the regression equation  $y = 0,88x + 4,8$ .

The results of the performance characteristics depend on the used analyzer.

### BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Kluthe R et al. Nature 1965; 205: 93-94.
5. Tarukosky PH et al. Scand J Clin Lab Invest 1966; 18: 80-86.
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
7. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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TEST INFORMATION		REAGENT VOLUME		
N°	**	Vol. R1	240	
Test	HAPTO	Vol. R2	60	
Full Name	Haptoglobine	Vol. R3		
Standard n°	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.1 Slope 1	
Vol. Sample Increas.		Unit	mg/dL Inter. 0	
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	340	Reagent Blank	20-21	
Sec. Wave.		React. Time	56-57	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
N° CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.1	10.0	90.0	3.0
1	CAL *0.25	30.0	90.0	3.0
2	CAL *0.4951	30.0	90.0	6.0
3	CAL *0.7354	30.0	90.0	9.0
4	CAL	0.0	0.0	3.0

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**BIBLIOGRAPHIE**

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Kluthe R et al. Nature 1965; 205: 93-94.
- Tarukosky PH et al. Scand J Clin Lab Invest 1966; 18: 80-86.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACCPres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACCPres, 1997.

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