Quantitative determination of Immunoglobulin G (IgG) IVD

Store at 2 - 8°C.

PRINCIPLE OF THE METHOD

Anti-human IgG antibodies when mixed with samples containing IgG, form insoluble complexes. These complexes cause an absorbance change, dependent upon the IgG concentration of the patient sample, that can be quantified by comparison from a calibrator of known IgG concentration.

CLINICAL SIGNIFICANCE

IgG is the most important immunoglobulin produced by plasma cells, and represents about 75% of the total immunoglobulins. Its main function is to neutralize toxins in tissue spaces. IgG deficit may be due to a congenital primary disturbance (immunodeficiency congenital and acquired) and is a special risk in children. Polyclonal hyperimmunoglobulinemia is the normal response to infections, especially in hepatitis and cirrhosis as well as autoimmune diseases. Increases of monoclonal IgG are found in multiple myeloma, lymphocytic leukemia, and Waldenström macroglobulinemia.

REAGENTS

R 1	Diluent	Tris buffer 20 mmol/L, PEG 8000, pH 8,3. Preservative.
R 2	Antibody	Goat serum, anti-human IgG, pH 7,5. Preservative.
Optional		Ref: 1102003 PROT CAL.

CALIBRATION

The assay is calibrated to the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL Calibrator to calibrate the reagent.

PREPARATION

Reagents: Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Between 700 - 1600 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Cod.:1102004). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME		
Nº	**	Vol. R1	240	
Test	IgG	Vol. R2	60	
Full Name	ImmunoglobulinG	Vol. R3		
Standard nº	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	2	Decimal	0.1	Slope 1
Vol. Sample Increases.		Unit	mg/dL	Inter. 0
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	605	Reagent Blank	40-41	
Sec. Wave.		React. Time	65-66	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
Nº CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.1	10.0	90.0	2.0
1	CAL *0.25	30.0	90.0	2.0
2	CAL *0.4967	30.0	90.0	4.0
3	CAL *0.7402	30.0	90.0	6.0
4	CAL	0.0	0.0	2.0

SPIN640Plus APPLICATION

EDIT PARAMETERS																																							
Test	IgG	No.	**																																				
Full name	Immunoglobulin G	Print name	IgG																																				
Reac. Type	End Point	Direction	Increase																																				
Pri. Wave.	605	Sec. Wave.																																					
Unit	mg/dL	Decimal	0.1																																				
Reagent Blank	46-47	React. Time	72 - 73																																				
Vol. Sample	2 ul	R1	240 ul																																				
Increased		R2	60 ul																																				
Decreased		R3																																					
Sample blank		R4																																					
CALIBRATION		<table style="width: 100%; border-collapse: collapse;"> <tr> <td>Spline</td> <td>Water</td> <td>0</td> <td>Vol aspirated</td> <td>Vol net</td> <td>Vol dil</td> </tr> <tr> <td>Calibrator</td> <td>*0.1</td> <td></td> <td>2</td> <td>13</td> <td>115</td> </tr> <tr> <td>Calibrator</td> <td>*0.25</td> <td></td> <td>2</td> <td>35</td> <td>105</td> </tr> <tr> <td>Calibrator</td> <td>*0.489</td> <td></td> <td>4</td> <td>35</td> <td>105</td> </tr> <tr> <td>Calibrator</td> <td>*0.735</td> <td></td> <td>6</td> <td>35</td> <td>105</td> </tr> <tr> <td>Calibrator</td> <td>(*)Cal Value</td> <td></td> <td>2</td> <td>--</td> <td>--</td> </tr> </table>		Spline	Water	0	Vol aspirated	Vol net	Vol dil	Calibrator	*0.1		2	13	115	Calibrator	*0.25		2	35	105	Calibrator	*0.489		4	35	105	Calibrator	*0.735		6	35	105	Calibrator	(*)Cal Value		2	--	--
Spline	Water	0	Vol aspirated	Vol net	Vol dil																																		
Calibrator	*0.1		2	13	115																																		
Calibrator	*0.25		2	35	105																																		
Calibrator	*0.489		4	35	105																																		
Calibrator	*0.735		6	35	105																																		
Calibrator	(*)Cal Value		2	--	--																																		

The Calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

1. **Measurement range:** Up to 3000 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and re-tested again. The linearity limit and measurement range depends on the sample to reagent / ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. **Detection Limit:** Values less than 10,3 mg/dL give non-reproducible results.
3. **Prozone effect:** No prozone effect was detected upon 8000 mg/dL
4. **Sensitivity:** Δ 0,6 mA. mg/dL at 359 mg/dL:
5. **Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

EP5	CV (%)		
	340,3 mg/dL	801,96 mg/dL	1517.5 mg/dL
Total	2,1%	2,8%	4,8%
Within Run	0,9%	0,7%	1%
Between Run	1,5%	1,5%	1,8%
Between Day	1%	2,2%	4,4%

6. **Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the Elecsys method from Roche. 79 samples ranging from 450 to 2600 mg/dL of IgG were assayed. The correlation coefficient (r) was 0.94 and the regression equation $y = 0,957x + 105,67$.

The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Skoug Jonh W et al. Clin Chem 1988; 34/2: 309 – 315
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14: 401-406.
5. Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Skoug Jonh W et al. Clin Chem 1988; 34/2: 309 - 315
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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TEST INFORMATION		REAGENT VOLUME		
N°	**	Vol. R1	240	
Test	IgG	Vol. R2	60	
Full Name	Immunoglobulin G	Vol. R3		
Standard n°	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	2	Decimal	0.1 Slope 1	
Vol. Sample Increas.		Unit	mg/dL Inter. 0	
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	605	Reagent Blank	40-41	
Sec. Wave.		React. Time	65-66	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
N° CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.1	10.0	90.0	2.0
1	CAL *0.25	30.0	90.0	2.0
2	CAL *0.4967	30.0	90.0	4.0
3	CAL *0.7402	30.0	90.0	6.0
4	CAL	0.0	0.0	2.0

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BIBLIOGRAPHIE

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Skoug Jonh W et al. Clin Chem 1988; 34/2: 309 - 315
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAFIA

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Skoug Jonh W et al. Clin Chem 1988; 34/2: 309 - 315
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
5. Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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Ref: MD1103004

Cont.

R1.Solvente: 2 x 40 mL

R2.Anticorpo: 1 x 20 mL