## Quantitative determination of Antithrombin III (ATHROM-III) IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Anti-antithrombin III antibodies when mixed with samples containing antithrombin-III, form insoluble complexes. These complexes cause an absorbance change, dependent upon the antithrombin-III concentration of the patient sample, that can be quantified by comparison from a calibrator of known antithrombin-III concentration.

### CLINICAL SIGNIFICANCE

Antithrombin III is a protein synthesized in the liver, normally present in the human plasma. It is the major inhibitor of the thrombin, and inhibits coagulation and limits the forming of blood clots. Antithrombin-III is also capable of activating other components of the coagulation cascade (eg, factor Xa), as well as plasmin.

Antithrombin-III deficiency can cause or lead to thrombosis, a clot forming in a blood vessel. Clots forming in the legs and pulmonary embolism are most commonly reported. Antithrombin-III deficiency is usually inherited and affects males and females equally. All family members should be tested if there is history of the disease.

Acquired antithrombin-III deficiency can occur as a result of other conditions. It has been reported in patients with liver diseases, patients receiving certain kinds of chemotherapy, and patients using oral contraceptives.

### REAGENTS

<b>R 1</b> Diluent	Tris buffer 20 mmol/L, PEG 8000, pH, 8,3. Preservative.
<b>R 2</b> Antibody	Goat serum, anti-human antithrombin III, pH 7,5. Preservative.
<b>Optional</b>	Ref: 1102003 PROT CAL.

### CALIBRATION

It must be used the PROT CAL to calibrate the reagent.

### PREPARATION

**Reagents:** Ready to use.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

**Reagent deterioration:** The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

### ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus autoanalyzer
- Laboratory equipment.

### SAMPLES

Fresh serum or plasma. Sodium citrate should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

Do not use highly hemolyzed or lipemic samples.

### REFERENCE VALUES

Between 17 - 30 mg/dL. Each laboratory should establish its own reference range.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

### SPIN640 APPLICATION

<b>TEST INFORMATION</b>		<b>REAGENT VOLUME</b>	
Nº	**	Vol. R1	240
Test	ATIII	Vol. R2	60
Full Name	Antithrombin III	Vol. R3	
Standard nº	6	Vol. R4	
<b>SAMPLE VOLUME</b>		<b>RESULT SETUP</b>	
Vol. Sample Stand.	8	Decimal	0.1
Vol. Sample Increas.		Unit	mg/dL
Vol. Sample Dec		Slope	1
		Inter.	0
<b>REACTION PARAMETERS</b>			
Reac. Type	End Point	Direction	Increase
Pri. Wave.	340	Reagent Blank	15-16
Sec. Wave.		React. Time	75-76
<b>CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)</b>			
CALIBRATION TYPE: SPLINE			
<b>Nº CAL DIL</b>	<b>CONCENTRATION</b>	<b>SAMPLE DIL</b>	<b>VOL DIL</b>
0	CAL *0.1	10.0	90.0
1	CAL *0.25	30.0	90.0
2	CAL *0.4936	30.0	90.0
3	CAL *0.7358	30.0	90.0
4	CAL	0.0	0.0

### SPIN640Plus APPLICATION

<b>EDIT PARAMETERS</b>			
Test	ATIII	No.	**
Full name	Antitrombin III	Print name	ATIII
Reac. Type	End Point	Direction	Increase
Pri. Wave.	340	Sec. Wave.	
Unit	mg/dL	Decimal	0.1
Reagent Blank	15 - 16	React. Time	81 - 82
Vol. Sample	6 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	
<b>CALIBRATION</b>	Spline	Water	0
	Calibrator	*0.1	6
	Calibrator	*0.25	6
	Calibrator	*0.489	12
	Calibrator	*0.735	18
	Calibrator	(*1)Cal Value	6

<i>Vol aspirated</i>	<i>Vol net</i>	<i>Vol dil</i>
6	13	115
6	35	105
12	35	105
18	35	105
6	--	--

The Calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

### PERFORMANCE CHARACTERISTICS

- Measurement range:** Up to 70 mg/dL (Nota 1), under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Limit detection:** Values less than 7,4 mg/dL give non-reproducible results.
- Prozone effect:** No prozone effect was detected upon 200 mg/dL.
- Sensitivity:** Δ 7,5 mA / mg/dL.
- Precision:** The reagent has been tested for 20 days, using two levels of serum in a EP5-based study.

EP5	CV (%)	
	17,85 mg/dL	35,93 mg/dL
Total	3,2%	3,5%
Within Run	0,8%	1%
Between Run	2,4%	2,4%
Between Day	2%	2,3%

- Accuracy:** Results obtained using this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

The results of the performance characteristics depend on the used analyzer.

### BIBLIOGRAPHY

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Buller HR et al. Critical care Medicine 1982; 10: 311.
- Kauffman et al. Am J Med 1978; 65: 607.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

### PACKAGING

Ref: MD1102014	Cont.	R1.Diluent:	2 x 40 mL
		R2.Antibody:	1 x 20 mL

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**BIBLIOGRAFIA**

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
3. Buller HR et al. Critical care Medicine 1982; 10: 311.
4. Kauffman et al. Am J Med 1978; 65: 607.
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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SAMPLE VOLUME		RESULT SETUP	
Vol. Sample Stand.	8	Decimal	0.1 Slope 1
Vol. Sample Increas.		Unit	mg/dL Inter. 0
Vol. Sample Dec			

REACTION PARAMETERS			
Reac. Type	End Point	Direction	Increase
Pri. Wave.	340	Reagent Blank	15-16
Sec. Wave.		React. Time	75-76

### CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)

CALIBRATION TYPE: SPLINE

Nº CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.1	10.0	90.0	8.0
1	CAL *0.25	30.0	90.0	8.0
2	CAL *0.4936	30.0	90.0	16.0
3	CAL *0.7358	30.0	90.0	24.0
4	CAL	0.0	0.0	8.0

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EDIT PARAMETERS																																													
Test	ATIII	No.	**																																										
Full name	Antithrombin III	Print name	ATIII																																										
Reac. Type	End Point	Direction	Increase																																										
Pri. Wave.	340	Sec. Wave.																																											
Unit	mg/dL	Decimal	0.1																																										
Reagent Blank	15 - 16	React. Time	81 - 82																																										
Vol. Sample	6 ul	R1	240 ul																																										
Increased		R2	60 ul																																										
Decreased		R3																																											
Sample blank		R4																																											
<table border="1"> <thead> <tr> <th>CALIBRATION</th> <th>Spline</th> <th>Water</th> <th>0</th> <th>Vol aspirated</th> <th>Vol net</th> <th>Vol dil</th> </tr> </thead> <tbody> <tr> <td>Calibrator</td> <td>*0.1</td> <td></td> <td></td> <td>6</td> <td>13</td> <td>115</td> </tr> <tr> <td>Calibrator</td> <td>*0.25</td> <td></td> <td></td> <td>6</td> <td>35</td> <td>105</td> </tr> <tr> <td>Calibrator</td> <td>*0.489</td> <td></td> <td></td> <td>12</td> <td>35</td> <td>105</td> </tr> <tr> <td>Calibrator</td> <td>*0.735</td> <td></td> <td></td> <td>18</td> <td>35</td> <td>105</td> </tr> <tr> <td>Calibrator</td> <td>(* )Cal Value</td> <td></td> <td></td> <td>6</td> <td>--</td> <td>--</td> </tr> </tbody> </table>				CALIBRATION	Spline	Water	0	Vol aspirated	Vol net	Vol dil	Calibrator	*0.1			6	13	115	Calibrator	*0.25			6	35	105	Calibrator	*0.489			12	35	105	Calibrator	*0.735			18	35	105	Calibrator	(* )Cal Value			6	--	--
CALIBRATION	Spline	Water	0	Vol aspirated	Vol net	Vol dil																																							
Calibrator	*0.1			6	13	115																																							
Calibrator	*0.25			6	35	105																																							
Calibrator	*0.489			12	35	105																																							
Calibrator	*0.735			18	35	105																																							
Calibrator	(* )Cal Value			6	--	--																																							

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