**Quantitative determination of Transferrin (TRF) IVD**

Store 2 - 8°C.

**PRINCIPLE OF THE METHOD**

Anti-transferrin antibodies when mixed with samples containing TRF, form insoluble complexes. These complexes cause an absorbance change, dependent upon the TRF concentration of the patient sample, that can be quantified by comparison from a calibrator of known TRF concentration.

**CLINICAL SIGNIFICANCE**

Transferrin is a plasma protein that contains a single polypeptide chain with approximately 6% carbohydrate. It is synthesized in the liver and transfers iron through the serum.

Evaluation of plasma TRF levels is useful for the differential diagnosis of anemia and for monitoring its treatment. In the hypochromic anemia (iron deficiency), the TRF level is increased. On the other hand, if the anemia is due to a failure to incorporate iron into erythrocytes, the TRF level is normal or low but the protein is highly saturated with iron. In iron overload, the TRF concentration is normal but saturation exceeds 55% and may be as great as 90%.

TRF concentration may, in fact, be used for assessing nutritional status. In the congenital defect atransferrinemia, very low level of TRF is accompanied by iron overload and a severe hypochromic anemia that is resistant to iron therapy. High levels of TRF occur in pregnancy and during estrogen administration.

**REAGENTS**

<b>R 1</b> Diluent	Tris buffer 20 mmol/L, PEG 8000, pH 8,3. Preservative.
<b>R 2</b> Antibody	Goat serum, anti-human transferrin, pH 7,5. Preservative.
<b>Optional</b>	Ref: 1102003 PROT CAL.

**CALIBRATION**

The assay has been standardized against the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent.

**PREPARATION**

**Reagents:** Ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

**Reagent deterioration:** The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

**ADDITIONAL EQUIPMENT**

- SPIN640 / SPIN640Plus autoanalyzer
- Laboratory equipment.

**SAMPLES**

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C. The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

**REFERENCE VALUES**

Between 200 - 360 mg/dL. Each laboratory should establish its own reference range.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**NOTES**

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

**SPIN640 APPLICATION**

TEST INFORMATION		REAGENT VOLUME		
Nº	**	Vol. R1	240	
Test	TRF	Vol. R2	60	
Full Name	Transferrin	Vol. R3		
Standard nº	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.1 Slope 1	
Vol. Sample Increases.		Unit	mg/dL Inter. 0	
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	340	Reagent Blank	20-21	
Sec. Wave.		React. Time	56-57	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
Nº CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.0625	10.0	150.0	3.0
1	CAL *0.125	25.0	175.0	3.0
2	CAL *0.25	30.0	90.0	3.0
3	CAL *0.4967	30.0	90.0	6.0
4	CAL	0.0	0.0	3.0

**SPIN640Plus APPLICATION**

EDIT PARAMETERS						
Test	TRF	No.	**			
Full name	Transferrin	Print name	TRF			
Reac. Type	End Point	Direction	Increase			
Pri. Wave.	340	Sec. Wave.				
Unit	mg/dL	Decimal	0.1			
Reagent Blank	20-21	React. Time	63 - 64			
Vol. Sample	3 ul	R1	240 ul			
Increased		R2	60 ul			
Decreased		R3				
Sample blank		R4				
CALIBRATION Spline		Water		Vol aspirated Vol net Vol dil		
Calibrator	*0.0625	0		3	10	150
Calibrator	*0.125			3	25	175
Calibrator	*0.25			3	30	90
Calibrator	*0.495			6	30	90
Calibrator	Cal Value			3	--	--

The Calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

**PERFORMANCE CHARACTERISTICS**

- Measurement range:** Up to the calibrator value (aprox. 800 mg/dL) under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection Limit:** 0 mg/dL.
- Prozone effect:** No prozone effect was detected upon 2000 mg/dL.
- Sensitivity:** Δ 3,0 mA / mg/dL (94 mg/dL).
- Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

EP5	CV (%)		
	77,02 mg/dL	206,99 mg/dL	377 mg/dL
Total	5,4%	2,5%	5,4%
Within Run	1%	0,8%	1,2%
Between Run	1,7%	1,3%	2,1%
Between Day	5%	2%	4,9%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the Immage from Beckman. 100 samples ranging from 50 to 700 mg/dL of TRF were assayed. The correlation coefficient (r) was 0,95 and the regression equation  $y = 1,046x + 3,843$ .

The results of the performance characteristics depend on the used analyzer.

**BIBLIOGRAPHY**

- McNeely Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Kreutzer HJH. J Clin Chem Clin Biochem 1976; 14: 401-406
- Young DS.Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3tn ed. AACC Pres, 1997.

**PACKAGING**

Ref: MD1102134	Cont.	R1.Diluent:	2 x 40 mL
		R2.Antibody:	1 x 20 mL

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TEST INFORMATION		REAGENT VOLUME		
Nº	**	Vol. R1	240	
Test	TRF	Vol. R2	60	
Full Name	Transferrin	Vol. R3		
Standard nº	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.1 Slope 1	
Vol. Sample Increas.		Unit	mg/dL Inter. 0	
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	340	Reagent Blank	20-21	
Sec. Wave.		React. Time	56-57	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
Nº CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.0625	10.0	150.0	3.0
1	CAL *0.125	25.0	175.0	3.0
2	CAL *0.25	30.0	90.0	3.0
3	CAL *0.4967	30.0	90.0	6.0
4	CAL	0.0	0.0	3.0

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EDIT PARAMETERS						
Test	TRF	No.	**			
Full name	Transferrin	Print name	TRF			
Reac. Type	End Point	Direction	Increase			
Pri. Wave.	340	Sec. Wave.				
Unit	mg/dL	Decimal	0.1			
Reagent Blank	20-21	React. Time	63 - 64			
Vol. Sample	3 ul	R1	240 ul			
Increased		R2	60 ul			
Decreased		R3				
Sample blank		R4				
CALIBRATION	Spline	Water	0	Vol aspirated Vol net Vol dil		
Calibrator	*0.0625			3	10	150
Calibrator	*0.125			3	25	175
Calibrator	*0.25			3	30	90
Calibrator	*0.495			6	30	90
Calibrator	Cal Value			3	--	--

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**BIBLIOGRAFIA**

- McNeely Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Kreutzer HJJ. J Clin Chem Clin Biochem 1976; 14: 401-406
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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TEST INFORMATION		REAGENT VOLUME		
N°	**	Vol. R1	240	
Test	TRF	Vol. R2	60	
Full Name	Transferrin	Vol. R3		
Standard n°	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.1 Slope 1	
Vol. Sample Inceas.		Unit	mg/dL Inter. 0	
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	340	Reagent Blank	20-21	
Sec. Wave.		React. Time	56-57	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
N° CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.0625	10.0	150.0	3.0
1	CAL *0.125	25.0	175.0	3.0
2	CAL *0.25	30.0	90.0	3.0
3	CAL *0.4967	30.0	90.0	6.0
4	CAL	0.0	0.0	3.0

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**BIBLIOGRAPHIE**

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Kreutzer HJH. J Clin Chem Clin Biochem 1976; 14: 401-406
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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TEST INFORMATION		REAGENT VOLUME			
Nº	**	Vol. R1	240		
Test	TRF	Vol. R2	60		
Full Name	Transferrin	Vol. R3			
Standard nº	6	Vol. R4			
SAMPLE VOLUME		RESULT SETUP			
Vol. Sample Stand.	3	Decimal	0.1	Slope	1
Vol. Sample Increas.		Unit	mg/dL	Inter.	0
Vol. Sample Dec					
REACTION PARAMETERS					
Reac. Type	End Point	Direction	Increase		
Pri. Wave.	340	Reagent Blank	20-21		
Sec. Wave.		React. Time	56-57		
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)					
CALIBRATION TYPE: SPLINE					
Nº CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE	
0	CAL *0.0625	10.0	150.0	3.0	
1	CAL *0.125	25.0	175.0	3.0	
2	CAL *0.25	30.0	90.0	3.0	
3	CAL *0.4967	30.0	90.0	6.0	
4	CAL	0.0	0.0	3.0	

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EDIT PARAMETERS						
Test	TRF	No.	**			
Full name	Transferrin	Print name	TRF			
Reac. Type	End Point	Direction	Increase			
Pri. Wave.	340	Sec. Wave.				
Unit	mg/dL	Decimal	0.1			
Reagent Blank	20-21	React. Time	63 - 64			
Vol. Sample	3 ul	R1	240 ul			
Increased		R2	60 ul			
Decreased		R3				
Sample blank		R4				
CALIBRATION	Spline	Water	0	Vol aspirated	Vol net	Vol dil
	Calibrator	*0.0625		3	10	150
	Calibrator	*0.125		3	25	175
	Calibrator	*0.25		3	30	90
	Calibrator	*0.495		6	30	90
	Calibrator	Cal Value		3	--	--

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**BIBLIOGRAFIA**

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Kreutzer HJH. J Clin Chem Clin Biochem 1976; 14: 401-406
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACCPres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACCPres, 1997.

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Ref.:MD1102134

Cont.

 R1: Diluyente: 2 x 40 mL  
 R2: Anticorpo: 1 x 20 mL