IgM

IgM (Immunoglobulin M)

Turbidimetry

Quantitative determination of Immunoglobulin M (IgM) IVD

Store at 2 - 8°C.

PRINCIPLE OF THE METHOD

Anti-human IgM antibodies when mixed with samples containing IgM, form insoluble complexes. These complexes cause an absorbance change, dependent upon the IgM concentration of the patient sample, that can be quantified by comparison from a calibrator of known IgM concentration.

CLINICAL SIGNIFICANCE

IgM is the only immunoglobulin that a neonate normally synthesizes, and in adults, it represents the 5-10% of the total immunoglobulins. Its structure is a pentamer of five IgG molecules and its high molecular weight (900.000 daltons) prevents its passage into extravascular spaces.

IgM concentration is decreased in diseases related with hereditary or acquired deficiencies of the immunoglobulin production. Polyclonal increases in serum immunoglobulins are the normal response to infections. The IgM generally increases as a primary response to virus infections and blood stream infections such as malaria and primary biliary cirrhosis. In multiple myeloma, if the paraprotein proves to be IgM, the diagnosis is probably Waldenström macroglobulinemia.

REAGENTS

R 1 Diluent	Tris buffer 20 mmol/L, PEG 8000, pH 8,3. Preservative.
R 2 Antibody	Goat serum, anti-human IgM, pH 7,5. Preservative.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

The assay is calibrated to the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL Calibrator to calibrate the reagent.

PREPARATION

Reagents: Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Between 40 - 230 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Cod.:1102004). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

Clinical diagnosis should not be made on findings of a single test result but should integrate both clinical and laboratory data.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME		
Nº	**	Vol. R1	240	
Test	IgM	Vol. R2	60	
Full Name	ImmunoglobulinM	Vol. R3		
Standard nº	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.1	Slope 1
Vol. Sample Increas.		Unit	mg/dL	Inter. 0
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	340	Reagent Blank	20-21	
Sec. Wave.		React. Time	56-57	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
Nº CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.1	10.0	90.0	3.0
1	CAL *0.25	30.0	90.0	3.0
2	CAL *0.4951	30.0	90.0	6.0
3	CAL *0.7356	30.0	90.0	9.0
4	CAL	0.0	0.0	3.0

SPIN640Plus APPLICATION

EDIT PARAMETERS																																													
Test	IgM	No.	**																																										
Full name	Immunoglobulin M	Print name	IgM																																										
Reac. Type	End Point	Direction	Increase																																										
Pri. Wave.	340	Sec. Wave.																																											
Unit	mg/dL	Decimal	0.1																																										
Reagent Blank	20-21	React. Time	63 - 64																																										
Vol. Sample	3 ul	R1	240 ul																																										
Increased		R2	60 ul																																										
Decreased		R3																																											
Sample blank		R4																																											
<table border="1"> <tr> <td>CALIBRATION</td> <td>Spline</td> <td>Water</td> <td>0</td> <td>Vol aspirated</td> <td>Vol net</td> <td>Vol dil</td> </tr> <tr> <td></td> <td>Calibrator</td> <td>*0.1</td> <td></td> <td>3</td> <td>13</td> <td>115</td> </tr> <tr> <td></td> <td>Calibrator</td> <td>*0.25</td> <td></td> <td>3</td> <td>35</td> <td>105</td> </tr> <tr> <td></td> <td>Calibrator</td> <td>*0.489</td> <td></td> <td>6</td> <td>35</td> <td>105</td> </tr> <tr> <td></td> <td>Calibrator</td> <td>*0.735</td> <td></td> <td>9</td> <td>35</td> <td>105</td> </tr> <tr> <td></td> <td>Calibrator</td> <td>(*) Cal Value</td> <td></td> <td>3</td> <td>--</td> <td>--</td> </tr> </table>				CALIBRATION	Spline	Water	0	Vol aspirated	Vol net	Vol dil		Calibrator	*0.1		3	13	115		Calibrator	*0.25		3	35	105		Calibrator	*0.489		6	35	105		Calibrator	*0.735		9	35	105		Calibrator	(*) Cal Value		3	--	--
CALIBRATION	Spline	Water	0	Vol aspirated	Vol net	Vol dil																																							
	Calibrator	*0.1		3	13	115																																							
	Calibrator	*0.25		3	35	105																																							
	Calibrator	*0.489		6	35	105																																							
	Calibrator	*0.735		9	35	105																																							
	Calibrator	(*) Cal Value		3	--	--																																							

The Calibration is stable until 31 days. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

1. Measurement range: Up to 300 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and re-tested again. The linearity limit and measurement range depends on the sample to reagent / ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

2. Detection Limit: Values less than 1 mg/dL give non-reproducible results.

3. Prozone effect: No prozone effect was detected upon 2000 mg/dL

4. Sensitivity: Δ 2,4 mA. mg/dL at 30 mg/dL.

5. Precision: The reagent has been tested for 20 days, using two levels of serum in a EP5-based study.

EP5	CV (%)	
	68,67 mg/dL	143,3 mg/dL
Total	5,7%	2,8%
Within Run	1,1%	0,7%
Between Run	3,8%	2,3%
Between Day	4,2%	1,3%

6. Accuracy: Results obtained using this reagent (y) were compared to those obtained using the Elecsys system from Roche. 100 samples ranging from 50 to 210 mg/dL of IgM were assayed. The correlation coefficient (r) was 0.958 and the regression equation $y = 0,974x + 1,296$.

The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Skoug Jonh W et al. Clin Chem 1988; 34/2: 309 - 315
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
5. Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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EP5	CV (%)	
		68,67 mg/dL
Total	5,7%	2,8%
Within Run	1,1%	0,7%
Between Run	3,8%	2,3%
Between Day	4,2%	1,3%

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BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Skoug Jonh W et al. Clin Chem 1988; 34/2: 309 - 315
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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APPLICATION AU SPIN640

TEST INFORMATION		REAGENT VOLUME		
N°	**	Vol. R1	240	
Test	IgM	Vol. R2	60	
Full Name	Immunoglobulin M	Vol. R3		
Standard n°	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.1 Slope 1	
Vol. Sample Increas.		Unit	mg/dL Inter. 0	
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	End Point	Direction	Increase	
Pri. Wave.	340	Reagent Blank	20-21	
Sec. Wave.		React. Time	56-57	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
N° CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.1	10.0	90.0	3.0
1	CAL *0.25	30.0	90.0	3.0
2	CAL *0.4951	30.0	90.0	6.0
3	CAL *0.7356	30.0	90.0	9.0
4	CAL	0.0	0.0	3.0

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BIBLIOGRAPHIE

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Skoug Jonh W et al. Clin Chem 1988; 34/2: 309 - 315
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
5. Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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