Quantitative determination of low levels of C-Reactive Protein IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

The CRP-ultra turbilatex is a quantitative turbidimetric test for the measurement of low levels of C- reactive protein (CRP) in human serum or plasma. Latex particles coated with specific anti- human CRP are agglutinated when mixed with samples containing CRP. The agglutination causes an absorbance change, dependent upon the CRP contents of the patient sample that can be quantified by comparison from a calibrator of known CRP concentration.

CLINICAL SIGNIFICANCE

CRP is an acute-phase protein present in normal serum, which increases significantly after most forms of tissue injuries, bacterial and virus infections, inflammation and malignant neoplasia. CRP may be also useful in detecting atherosclerotic process and providing important prognostic information about patients with asymptomatic heart disease, unstable angina, and myocardial infarction. Recent studies in apparently healthy people show that CRP concentration in serum rise long before traditional symptoms of heart and vascular diseases are noticed.

REAGENTS

Diluent-ultra (R1)	Tris buffer 20 mmol/L, pH 8,2. Preservative.
Latex-ultra (R2)	Latex particles coated with goat IgG anti- human CRP, pH 7,3. Preservative.
U-CRP CAL	Liquid Calibrator. C-Reactive protein concentration is stated on the vial label.

PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

Use CRP Ultra Calibrator Reference 43035. The sensitivity of the assay and the target value of the calibrator have been standardized against the Reference Material ERM-DA 474/IFCC. Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

Ready for use.
CRP Calibrator: Ready for use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Reagent deterioration: Presence of particles (R1, R2) and turbidity (R1). Do not freeze; frozen Latex or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- SPIN640 / SPIN640Plus autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C. The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Below 3 mg/L is considered normal. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT CRP Ultra Control (Ref.:43036). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN640 APPLICATION

TEST INFORMATION		REAGENT VOLUME		
Nº	**	Vol. R1	240	
Test	CRP _u	Vol. R2	60	
Full Name	CRP-Ultrasetensive	Vol. R3		
Standard nº	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.01	Slope 1
Vol. Sample Inceas.		Unit	mg/L	Inter. 0
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	Fixed Time	Direction	Increase	
Pri. Wave.	546	Reagent Blank	20-21	
Sec. Wave.		React. Time	44-70	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
Nº CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.05	5.0	95.0	3.0
1	CAL *0.10	10.0	90.0	3.0
2	CAL *0.25	40.0	120.0	3.0
3	CAL *0.4967	40.0	120.0	6.0
4	CAL	0.0	0.0	3.0

SPIN640Plus APPLICATION

EDIT PARAMETERS			
Test	CRP _u	No.	**
Full name	CRP _u	Print name	CRP _u
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	546	Sec. Wave.	
Unit	mg/L	Decimal	0.1
Reagent Blank	20-21	React. Time	51 - 77
Vol. Sample	3 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	
CALIBRATION	Spline	Water	0
	Calibrator	*0.0625	3 10 150
	Calibrator	*0.125	3 25 175
	Calibrator	*0.25	3 30 90
	Calibrator	*0.4951	6 30 90
	Calibrator	*	3 -- --

The Calibration is stable until 31 days. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

- Linearity limit:** Up to 10 mg/L, under the described assay conditions. Samples with higher concentrations should be diluted 1/3 in NaCl 9 g/L and re-tested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection limit:** Values less than 0,05 mg/L give non-reproducible results.
- Prozone effect:** No prozone effect was detected upon 800 mg/L.
- Sensitivity:** Δ 44 mA.mg/L.
- Precision:** The reagent has been tested for 20 days, using three different CRP concentrations in a EP5-based study.

EP5	CV (%)		
	+/- 0,28 mg/L	+/- 3,09 mg/L	+/- 5,95 mg/L
Total	7,7%	2,7%	3,0%
Within Run	4,5%	1,7%	1,4%
Between Run	4,7%	1,9%	2,7%
Between Day	4,1%	0,7%	0,0%

6. Accuracy: Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 23 samples of different concentrations of CRP were assayed. The correlation coefficient (r) was 0,99 and the regression equation $y=1,0028x-0,0625$. The results of the performance characteristics depend on the analyzer used.

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APPLICATION AU SPIN640

TEST INFORMATION		REAGENT VOLUME		
N°	**	Vol. R1	240	
Test	CRPu	Vol. R2	60	
Full Name	CRP-Ultrasensitive	Vol. R3		
Standard n°	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.01 Slope 1	
Vol. Sample Increas.		Unit	mg/L Inter. 0	
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	Fixed Time	Direction	Increase	
Pri. Wave.	546	Reagent Blank	20-21	
Sec. Wave.		React. Time	44-70	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
N° CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.05	5.0	95.0	3.0
1	CAL *0.10	10.0	90.0	3.0
2	CAL *0.25	40.0	120.0	3.0
3	CAL *0.4967	40.0	120.0	6.0
4	CAL	0.0	0.0	3.0

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5. Precision: The reagent has been tested for 20 days, using three different CRP concentrations in a EP5-based study.

EP5	CV (%)		
	+/- 0,28 mg/L	+/- 3,09 mg/L	+/- 5,95 mg/L
Total	7,7%	2,7%	3,0%
Sur une goutte	4,5%	1,7%	1,4%
Entre la goutte	4,7%	1,9%	2,7%
Entre le jour	4,1%	0,7%	0,0%

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BIBLIOGRAPHIE

- Thomas A et al. IVD Technology 2000; March/April: 27-35.
- Macy E M et al. Clinical Chemistry 1997; 43: 52-58.
- Pearson TA et al. Circulation 2003;107:499-511.
- Haverkate F et al. Fibrinolysis and Proteolysis 1007; 11: 1331-134.
- Ronald D et al. Journal of Clinical Ligand Assay 1997; 313-315.
- Ridker PM et al. The New England Journal of Medicine 2000; 23: 836-843.
- Koenig W et al. Circulation 1999; 99: 237-242.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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TEST INFORMATION		REAGENT VOLUME		
Nº	**	Vol. R1	240	
Test	CRPu	Vol. R2	60	
Full Name	CRP-Ultrasensitive	Vol. R3		
Standard nº	6	Vol. R4		
SAMPLE VOLUME		RESULT SETUP		
Vol. Sample Stand.	3	Decimal	0.01	Slope 1
Vol. Sample Inceas.		Unit	mg/L	Inter. 0
Vol. Sample Dec				
REACTION PARAMETERS				
Reac. Type	Fixed Time	Direction	Increase	
Pri. Wave.	546	Reagent Blank	20-21	
Sec. Wave.		React. Time	44-70	
CALIBRATION DILUTION (5 CAL + 1 WATER BLANK)				
CALIBRATION TYPE: SPLINE				
Nº CAL DIL	CONCENTRATION	SAMPLE DIL	VOL DIL	VOL SAMPLE
0	CAL *0.05	5.0	95.0	3.0
1	CAL *0.10	10.0	90.0	3.0
2	CAL *0.25	40.0	120.0	3.0
3	CAL *0.4967	40.0	120.0	6.0
4	CAL	0.0	0.0	3.0

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EDIT PARAMETERS			
Test	CRPu	No.	**
Full name	CRPu	Print name	CRPu
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	546	Sec. Wave.	
Unit	mg/L	Decimal	0.1
Reagent Blank	20-21	React. Time	51 - 77
Vol. Sample	3 ul	R1	240 ul
Increased		R2	60 ul
Decreased		R3	
Sample blank		R4	
CALIBRATION			
Spline	Water	0	
Calibrator	*0.0625	3	10 150
Calibrator	*0.125	3	25 175
Calibrator	*0.25	3	30 90
Calibrator	*0.4951	6	30 90
Calibrator	*	3	-- --

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BIBLIOGRAFIA

- Thomas A et al. IVD Technology 2000; March/April: 27-35.
- Macy E M et al. Clinical Chemistry 1997; 43: 52-58.
- Pearson TA et al. Circulation 2003;107:499-511.
- Haverkate F et al. Fibrinolysis and Proteolysis 1007; 11: 1331-134.
- Ronald D et al. Journal of Clinical Ligand Assay 1997; 313-315.
- Ridker PM et al. The New England Journal of Medicine 2000; 23: 836-843.
- Koenig W et al. Circulation 1999; 99: 237-242.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AAC Press, 1995.

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Ref: MD43134

Cont.

 R1. Diluyente: 2 x 40 mL
 R2. Látex: 1 x 20 mL
 U-CRP CAL: 1 x 2 mL