
Quantitative determination of α-amylase (AMS) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

α-Amylase hydrolyzes the 2-chloro-4-nitrophenyl-α-D-maltotriose (CNPG₃) to release 2-chloro-4-nitrophenol (CNP) and form 2-chloro-4-nitrophenyl-α-D-maltoside (CNPG₂), maltotriose (G₃) y glucose (G) according to the following reaction:

The rate of 2-chloro-4-nitrophenol formation, measured photometrically, is proportional to the catalytic concentration of α-amylase present in the sample¹.

CLINICAL SIGNIFICANCE

α-Amylase (AMS) is an enzyme that helps to digest the glycogen and the starch. It is produced mainly by exocrine pancreas and salivary glands. This determination is made mainly in diagnosis or to control diseases of the pancreas as acute or chronic pancreatitis. It can also reflect biliary or gastrointestinal disease and other upheavals^{2,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	MES pH 6.0	100 mmol/L
	CNPG ₃	2,25 mmol/L
	Sodium chloride	350 mmol/L
	Calcium acetate	6 mmol/L
	Potassium thiocyanate	900 mmol/L
	Sodium azide	0.95 gr/L

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

After opening, the reagent is stable for 60 days when properly capped immediately after each opening and stored at 2-8°C.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm ≥ 0,40.

ADDITIONAL EQUIPMENT

- SPIN 800 Autoanalyzer.
- General laboratory equipment (Note 2).

SAMPLES

- Serum or plasma¹, remove from cells as soon as possible. It is recommended to use heparin as anticoagulant.

- Urine, adjust pH to approximately 7.0 prior to storage.

Stability: 1 month at 2-8°C.

REFERENCE VALUES¹

Serum or plasma	Up to 90 U/L of α-amylase
Urine	Up to 450 U/L of α-amylase

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN 800 APPLICATION

EDIT PARAMETERS			
Test	AMY	No.	**
Full name	AMYLASE	Print name	AMY
Reac. Type	Kinetic	Direction	Increase
Pri. Wave.	412	Sec. Wave.	700
Unit	U/L	Decimal	0
Reagent Blank	-	React. Time	6 - 16
Vol. Sample	6 ul	Increas.	Decrease
			Sample blank
R1	300 ul	R2	R3 R4
CALIBRATION			
Rule	Linear two points	Water	0
		Calibrator	*

The Calibration is stable until **36 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 0,2439 U/L to linearity limit of 2200 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	77	194	77
SD	1,12	2,22	1,08	2,96
CV (%)	1,45	1,15	1,39	1,50

Sensitivity: 1 U/L = 0,00025 ΔA / min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained were the following:

Correlation coefficient (r)²: 0,98628.

Regression equation: y=0,746x - 1,2697.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. α-Amylase enzyme activity is temperature dependent. Assays performed at temperatures <37°C or >37°C will show an apparent decrease or increase levels.
2. Saliva and sweat contain α-amylase. Avoid mouth pipetting and skin contact with the reagent or material used.
3. Contains potassium thiocyanate. Avoid inhalation, skin or eyes contact. If it happens, wash with plenty of water and consult a doctor.

BIBLIOGRAPHY

1. Ying Foo A et al. Amylase measurement with 2-chloro-4nitrophenyl maltotriose as substrate. Clin Chim 272, 1998; 137-147.
2. McNeely M. Amylase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1112-116.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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1. Ying Foo A et al. Amylase measurement with 2-chloro-4nitrophenyl maltotriósido as substrate. Clin Chim 272, 1998; 137-147.
2. McNeely M. Amylase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1112-116.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAFIA

1. Ying Foo A et al. Amylase measurement with 2-chloro-4-nitrophenyl maltotrioside as substrate. Clin Chim 272, 1998; 137-147.
2. McNeely M. Amylase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1112-116.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: MX41201	Cont.	R:6x 40 mL
--------------	-------	------------