**Quantitative determination of calcium IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Calcium with Arsenazo III (1,8-Dihydroxy-3,6-disulpho-2,7-naphthalene-bis (azo)-dibenzene-arsonic acid), at neutral pH, yields a blue colored complex.

The intensity of the colour formed is proportional to the calcium concentration in the sample<sup>1,2,3</sup>.

**CLINICAL SIGNIFICANCE**

Calcium is the most abundant and one of the most important minerals in the human body. Approximately 99% of body calcium is found in bones.

A decrease in albumin level causes a decrease in serum calcium. Low levels of calcium are found in hypoparathyroidism, pseudohypoparathyroidism, vitamin D deficiency, malnutrition and intestinal malabsorption.

Among causes of hypercalcemia are cancers, large intake of vitamin D, enhanced renal retention, osteoporosis, sarcoidosis, thyrotoxicosis, hyperparathyroidism<sup>1,6,7</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

R	Imidazol Buffer pH 6,5	100 mmol/L
Arsenazo III	Arsenazo III	120 mmol/L

**PRECAUTIONS**

R: H360- May damage fertility or the unborn child.

Follow the precautionary statements given in MSDS and label of the product.

**PREPARATION**

The reagent is ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations are prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 660 nm  $\geq$  0.50.

**ADDITIONAL EQUIPMENT**

- SPIN 800 Autoanalyzer.
- General laboratory equipment (Note 1, 2).

**SAMPLES**

- Serum or plasma<sup>1</sup>: Separated from cells as rapidly as possible. Blood anticoagulants with oxalate or EDTA are not acceptable since these chemicals will strongly chelate calcium.

- Urine<sup>1</sup>: Collect 24 hour urine specimen in calcium free containers. The collecting bottles should contain 10 ml of diluted Nitric acid (50% v/v). Record the volume.

Dilute a sample 1/2 in distilled water. Mix. Multiply results by 2 (dilution factor).

Stability of the samples: Calcium is stable 10 days at 2-8°C.

**QUALITY CONTROL**

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

**SPIN 800 APPLICATION**

EDIT PARAMETERS			
Test	CA	No.	**
Full name	CA - ARS	Print name	CA
Reac. Type	End Point	Direction	Increase
Pri. Wave.	660	Sec. Wave.	
Unit	mg/dL	Decimal	0.01
Reagent Blank	3 - 4	React. Time	22 - 24
Vol. Sample	3 ul	Increas.	Decrease
			Sample blank
R1	300 ul	R2	R3
			R4
CALIBRATION			
Rule	Linear two points	Water	0
		Calibrator	*

The Calibration is stable until **28 days**. After this period the Calibration must be performed again in order to obtain good results.

**REFERENCE VALUES<sup>1</sup>**

Serum or plasma:

Adults 8.5-10.5 mg /dL  $\cong$  2.1-2.6 mmol/L

Children 10 -12 mg/dL  $\cong$  2.5 - 3 mmol/L

Newborns 8 -13 mg/dL  $\cong$  2 - 3.25 mmol/L

Urine:

Adults 50 - 300 mg/24h  $\cong$  1.25 - 7.5 mmol/24h

Children 80 -160 mg/24h  $\cong$  2 - 4 mmol/24h

These values are for orientation purpose; each laboratory should establish its own reference range.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From detection limit of 0,026 mg/dL to linearity limit of 32 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	8,35	14,28	8,58	14,57
SD	0,08	0,08	0,19	0,34
CV (%)	0,95	0,59	2,24	2,31

**Sensitivity:** 1 mg/dL = 0,0316 A.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r): 0,9506

Regression equation:  $y = 0,8944x + 1,3421$ .

The results of the performance characteristics depend on the analyzer used.

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1. Farrell E C. Calcium. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1051-1255 and 418.
2. Kessler G. et al. Clin Chem 1964; 10 (8); 686-706.
3. Connerty H. V. et al. Am J Clin Path 1996; 45 (3); 200-296.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
6. Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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### BIBLIOGRAPHIE

1. Farell E C. Calcium. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1051-1255 and 418.
2. Kessler G. et al. Clin Chem 1964; 10 (8); 686-706.
3. Connerty H. V. et al. Am J Clin Path 1996; 45 (3); 200-296.
4. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
5. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed. AACC 2001.
6. Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3rd ed. AACC 1999.
7. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed. AACC 1995.

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