
**Quantitative determination of uric acid**

IVD

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

Uric acid is oxidized by uricase to allantoine and hydrogen peroxide ( $2\text{H}_2\text{O}_2$ ), which under the influence of POD, 4-aminophenazone (4-AP) and 2-4 Dichlorophenol sulfonate (DCPS) forms a red quinoneimine compound:

The intensity of the red color formed is proportional to the uric acid concentration in the sample<sup>1,2</sup>.

**CLINICAL SIGNIFICANCE**

Uric acid and its salts are end products of the purine metabolism. With progressive renal insufficiency, there is retention in blood of urea, creatinine and uric acid. Elevate uric acid level may be indicative of renal insufficiency and is commonly associated with gout<sup>1,5,6</sup>. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b> Buffer	Phosphate pH 7.4 2-4 Dichlorophenol sulfonate (DCPS)	50 mmol/L 4 mmol/L
<b>R 2</b> Enzymes	Uricase Peroxidase (POD) Ascorbate oxidase 4 – Aminophenazone (4-AP)	60 U/L 660 U/L 200 U/L 1 mmol/L

**PREPARATION**

All the reagents are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 520 nm ≥ 0.16.

**ADDITIONAL EQUIPMENT**

- SPIN 800 Autoanalyzer.
- General laboratory equipment.

**SAMPLES**

- Serum or plasma<sup>1</sup>: Stability 3-5 days at 2-8°C or 6 months at – 20°C.
- Urine (24 h)<sup>1</sup>: Stability 4 days at 15-25°C, pH >8. Dilute sample 1/50 in distilled water. Mix. Multiply results by 50 (dilution factor); If urine is cloudy; warm the specimen to 60°C for 10 min to dissolve precipitated urates and uric acid. Do not refrigerate.

**QUALITY CONTROL**

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

**SPIN800 APPLICATION**

EDIT PARAMETERS					
Test	<b>URIC</b>	No.	**		
Full name	URIC ACID	Print name	URIC		
Reac. Type	<b>Endpoint</b>	Direction	Increase		
Pri. Wave.	505	Sec. Wave.			
Unit	mg/dL	Decimal	0.01		
Reagent Blank	14 - 16	React. Time	32 - 33		
Vol. Sample	7.5 ul	Increas.	Decrease	Sample blank	
R1	150 ul	R2	150 ul	R3	R4
CALIBRATION					
Rule	Linear two points	Water	0		
		Calibrator	*		

The Calibration is stable until **36 days**. After this period the Calibration must be performed again in order to obtain good results.

**REFERENCE VALUES<sup>1</sup>**

Serum or plasma:

Women 2,5 – 6,8 mg/dL ≈ 149 – 405 μmol/L

Men 3,6 – 7,7 mg/dL ≈ 214 – 458 μmol/L

Urine: 250 - 750 mg/24 h ≈ 1,49 – 4,5 mmol/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From *detection limit* of 0,01647 mg/dL to *linearity limit* of 40 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)	Inter-assay (n=20)
Mean (mg/L)	4,46	4,71
SD	0,02	0,06
CV (%)	0,46	0,44

**Sensitivity:** 1 mg/dL = 0,0323 (A).

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,99734.

Regression equation:  $y=0,816x + 0,319$ .

The results of the performance characteristics depend on the analyzer used.

**NOTES**

1. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
2. Use clean disposable pipette tips for its dispensation.

**BIBLIOGRAPHY**

1. Schultz A. Uric acid. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1261-1266 and 418.
2. Fossati P et al. Clin Chem 1980;26:227-231.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995

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