Quantitative determination of α₁-antitrypsin (α₁-ATRYP) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Anti-α₁-antitrypsin antibodies when mixed with samples containing α₁-antitrypsin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the α₁-antitrypsin concentration of the patient sample, that can be quantified by comparison from a calibrator of known α₁-antitrypsin concentration.

CLINICAL SIGNIFICANCE

α₁-antitrypsin is a glycoprotein synthesized in the hepatic parenchyma cells that circulates in the bloodstream. It is the second highest proteinase inhibitor in plasma after α₂-macroglobulin. α₁-antitrypsin is a strong reactor with elastase, skin collagenase, chymotrypsin, plasmin, and thrombin, and also shows inhibitory activity against fungal and leukocytic proteases.

α₁-antitrypsin deficiency is an inherited disorder, and occurs when both parents pass on an abnormal gene (PiZ) to their child. This deficiency is associated with a very high risk for development of pulmonary emphysema and diseases of the liver - neonatal cholestasis, hepatitis, cirrhosis and hepatocellular carcinoma.

α₁-antitrypsin increases in inflammatory or necrosis process. Serum levels begin to rise after approximately 24 hours and peak at 3 or 4 days if the insult is acute and short-lived.

REAGENTS

R 1 Diluent	Tris buffer 20 mmol/L, PEG 8000, pH 8,3. Preservative.
R 2 Antibody	Goat serum, anti-human α ₁ -antitrypsin, pH 7,5. Preservative.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

The assay has been standardized against the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent.

PREPARATION

Reagents: Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- SPIN 800 autoanalyzer.
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C. Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Newborn: Between 124 - 348 mg/dL.

Adults: 90 - 200 mg/dL.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN 800 APPLICATION

EDIT PARAMETERS

Test	ATRYP	No.	**
Full name	ANTITRIPSIN	Print name	ATRYP
Reac. Type	End Point	Direction	Increase
Pri. Wave.	340	Sec. Wave.	
Unit	mg/dL	Decimal	0.1
Reagent Blank	15 -16	React. Time	24 - 25
Vol. Sample	3 ul	Decrease	Sample blank
R1	240 ul	R2	60 ul
		R3	R4

CALIBRATION

Spline	Water	0	Vol aspirated	Vol net	Vol dil
	Calibrator	*0.1	3	10	90
	Calibrator	*0.25	3	30	90
	Calibrator	*0.495	6	30	90
	Calibrator	*0.735	9	30	90
	Calibrator	(*)Cal value	3	--	--

CAL AUTO DILUTIONS

The Calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

- Measurement range:** Up to 500 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection Limit:** Values less than 16 mg/dL give non-reproducible results.
- Prozone effect:** No prozone effect was detected upon 1000 mg/dL
- Sensitivity:** Δ 3,4 mA / mg/dL.
- Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

EP5	CV (%)		
	34,39 mg/dL	92,7 mg/dL	181,8 mg/dL
Total	4,2%	2,6%	2,8%
Within Run	0,8%	1,1%	1,6%
Between Run	3,8%	2,4%	2,3%
Between Day	1,6%	0%	0%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using the Beckman Array 360 CE. 35 samples ranging from 70 to 250 mg/dL of α₁-antitrypsin were assayed. The correlation coefficient (r) was 0.92 and the regression equation $y = 0,84x + 26,5$.

The results of the performance characteristics depend on the used analyzer.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Sharp HL. Hospital Practice; May 1971: 83-96
- Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clin. laboratory tests, 3th ed. AACC Pres, 1997.

PACKAGING

Ref: MX1102054	Cont.	R1.Diluent:	2 x 40 mL
		R2.Antibody:	1 x 20 mL

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BIBLIOGRAPHIE

- Clinical Guide to Laboratory Tests, Edited by NWTietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- PesceAJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Sharp HL. Hospital Practice; May 1971: 83-96
- Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clin. laboratory tests, 3rd ed. AACC Pres, 1997.

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BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Sharp HL. Hospital Practice; May 1971: 83-96
- Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
- Young DS.Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clin. laboratory tests, 3tn ed. AACC Pres, 1997.

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Ref.: MX1102054	Cont.	R1. Solvente: 2 x 40 mL
		R2. Anticorpo: 1 x 20 mL