Quantitative determination of complement C4 (C4) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

Anti-human C4 antibodies when mixed with samples containing C4, form insoluble complexes. These complexes cause an absorbance change, dependent upon the C4 concentration of the patient sample, that can be quantified by comparison from a calibrator of known C4 concentration.

CLINICAL SIGNIFICANCE¹

C4 is the second component reacting in the classical pathway cascade. Most synthesis occurs in the hepatic parenchymal cells, although some may be synthesized by monocytes or other tissues.

C4 levels in plasma rise modestly after trauma or inflammation and tissue necrosis (acute phase process).

Inherited primary deficiency of C4 is associated with a high prevalence of autoimmune or collagen vascular disease, particularly Systemic Lupus Erythematosus (SLE). Also, levels of C4 are more commonly depressed because of consumption as a consequence of formed immune-complexes.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Preservative.
Antibody (R2)	Goat serum, anti-human C4, pH 7.5. Preservative.
Optional	Cod: 1102003 PROT CAL.

CALIBRATION

The assay is calibrated to the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent.

PREPARATION

Reagents: Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not use.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- SPIN 800 autoanalyzer.
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant.

Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES⁵

Neonates: Between 13 - 38 mg/dL.

Adults: Between 10 - 40 mg/dL.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Cod.:1102004). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

SPIN 800 APPLICATION

EDIT PARAMETERS

Test	C4	No.	**
Full name	C4	Print name	C4
Reac. Type	End Point	Direction	Increase
Pri. Wave.	340	Sec. Wave.	
Unit	mg/dL	Decimal	0
Reagent Blank	15 - 16	React. Time	24 - 25
Vol. Sample	6 ul	Increas.	Decrease
			Sample blank
R1	240 ul	R2	60 ul
		R3	R4

CALIBRATION

Spline	Water	0			
	Calibrator	*0.1	6	13	115
	Calibrator	*0.25	6	35	105
	Calibrator	*0.489	12	35	105
	Calibrator	*0.735	18	35	105
	Calibrator	(*)Cal value	6	--	--

Vol aspirated	Vol net	Vol dil
6	13	115
6	35	105
12	35	105
18	35	105
6	--	--

CAL AUTO DILUTIONS

The Calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

1. Measurement range: Up to 75 mg/dL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and re-tested again. The linearity limit and measurement range depends on the sample to reagent / ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
2. Detection Limit: Values less than 1 mg/dL give non-reproducible results.
3. Prozone effect: No prozone effect was detected upon 500 mg/dL.
4. Sensitivity: Δ 23,6 mA. mg/dL (5 mg/dL), Δ 12,9 mA. mg/dL (37mg/dL).
5. Precision: The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

EP5	CV (%)		
	8,57 mg/dl	22,46 mg/dl	42,98 mg/dl
Total	3,9%	2,4%	1,9%
Within Run	1,6%	1%	1%
Between Run	2,2%	1,6%	1,1%
Between Day	2,8%	1,4%	1,2%

6. Accuracy: Results obtained using this reagent (y) were compared to those obtained using an immunoturbidimetric method from Bayer. 46 samples ranging from 9 to 60 mg/dL of C4 were assayed. The correlation coefficient (r) was 0,97 and the regression equation $y = 1,16x - 1,9$.

The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Yang Y et al. Curr Dir Autoimmun 2004; 7: 98-132.
3. Borque L et al. Clin Biochem 1983; 16: 330-333.
4. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
5. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
6. Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
7. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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EDIT PARAMETERS

Test	C4	No.	**
Full name	C4	Print name	C4
Reac. Type	End Point	Direction	Increase
Pri. Wave.	340	Sec. Wave.	
Unit	mg/dL	Decimal	0
Reagent Blank	15 -16	React. Time	24 - 25
Vol. Sample	6 ul	Decrease	Sample blank
R1	240 ul	R2	60 ul
R3		R4	

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BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Yang Y et al. Curr Dir Autoimmun 2004; 7: 98-132.
- Borque L et al. Clin Biochem 1983; 16: 330-333.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34: 517-520.
- Young DS. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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