APO A-I

Apolipoprotein A-I

Turbidimetry

Quantitative determination of apolipoprotein A-I (APO A-I) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

Turbidimetric test for the measurement of apolipoprotein A-I in human serum or plasma.

Anti- Apo A-I antibodies when mixed with samples containing Apo A-I, form insoluble complexes. These complexes cause an absorbance change, dependent upon the Apo A-I concentration of the patient sample, that can be quantified by comparison from a calibrator of known Apo A-I concentration.

CLINICAL SIGNIFICANCE¹

Apo A-I is the major structural apolipoprotein in HDL and constitutes about 70% of the total protein. Apo A-I is a cofactor for lecithin-cholesterol-acyl-transferase (LCAT), the enzyme responsible for forming cholesteryl esters in plasma and plays an important role in the transport of cholesterol from peripheral tissues to the liver, to be finally excreted. Measurements of Apo A-I concentration is specially important in detecting coronary heart disease risk (CHD) as well as in the diagnostic of hyperlipoproteinemia. Concentrations < 120 mg/L are associated to an increased CHD risk, while concentrations ≥ 160 mg/L may even protect from the same risk. Patients with deficiencies in Apo A-I synthesis may highly increase the CHD risk.

Tanger disease, a consequence of an Apo A-I catabolism defect, is characterized by several reduced plasma HDL cholesterol (HDL-c) concentration, abnormal HDL composition and accumulation of cholesteryl esters in many body tissues. Plasma HDL-c and Apo A-I concentrations in homozygotes are very low, while Apo A-II concentration is less than 10% of its normal concentration. Heterozygotes are characterized by half-normal concentration of HDL-c, Apo AI and Apo -II. Current evidence suggests that these patients have increased incidence of CHD.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, PEG, pH 8,3. Preservative.
Antibody (R2)	Goat serum, anti-human Apo B, tris 50 mmol/L, pH 7,5. Preservative.
Optional	APO CAL ref: 93005

CALIBRATION

The assay and the value of the calibrator concentration have been standardized against the Certified Reference Material WHO/IFCC SP1-01 (CDC, USA). It is recommended the use of the APO CAL Calibrator for calibration.

PREPARATION

Reagents: Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- SPIN 800 autoanalyzer.
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 2 weeks at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES⁵

Between 122 – 161 mg/dL.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact Apolipoprotein Control (Ref.:93006) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN 800 APPLICATION

EDIT PARAMETERS

Test	APO A1	No.	**
Full name	APO A1	Print name	APO A1
Reac. Type	End Point	Direction	Increase
Pri. Wave.	340	Sec. Wave.	
Unit	mg/L	Decimal	0.1
Reagent Blank	15 -16	React. Time	25 - 26
Vol. Sample	2,1 ul	Increas.	Decrease
			Sample blank
R1	240 ul	R2	60 ul
		R3	R4

CALIBRATION	Spline Water	0			
	Calibrator	*0.1	2,1	10	90
	Calibrator	*0.25	2,1	30	90
	Calibrator	*0.497	4,2	30	90
	Calibrator	*0.725	6,3	30	90
	Calibrator	*	2,1	--	--

	<i>Vol aspirated</i>	<i>Vol net</i>	<i>Vol dil</i>
	2,1	10	90
	2,1	30	90
	4,2	30	90
	6,3	30	90
	2,1	--	--
CAL AUTO DILUTIONS			

The Calibration is stable until **27 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

- Measurement range:** Up to 250 mg/dL, under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection Limit:** Values less than 0,1 mg/dL give non-reproducible results.
- Precision:** The reagent has been tested for 20 days, using three levels of serum in a EP5-based study (NCCLS).

EP5	CV (%)		
	27,22mg/dL	65,74 mg/dL	131,07 mg/dL
Total	4%	3,7%	4,8%
Within Run	2,2%	0,8%	1,1%
Between Run	2,3%	1,3%	1,4%
Between Day	2,4%	3,3%	4,5%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained with a Bayer immunoturbidimetric method. 39 samples ranging from 50 to 200 mg/dL of Apo A-I were assayed. The correlation coefficient (r) was 0,92 and the regression equation $y = 1,18x - 37,8$. The results of the performance characteristics depend on the used analyzer.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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