Quantitative determination of Rheumatoid Factors (RF)
IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

The RF-Turbilatex is a quantitative turbidimetric test for the measurement of RF in human serum or plasma.

Latex particles coated with human gammaglobulin are agglutinated when mixed with samples containing RF. The agglutination causes an absorbance change, dependent upon the RF contents of sample that can be quantified by comparison from a calibrator of known RF concentration.

CLINICAL SIGNIFICANCE

Rheumatoid factors are a group of antibodies directed to determinants in the Fc portion of the immunoglobulin G molecule. Although rheumatoid factors are found in a number of rheumatoid disorders, such as systemic lupus erythematosus (SLE) and Sjögren's syndrome, as well as in nonrheumatic conditions, its central role in clinic lies its utility as an aid in the diagnosis of rheumatoid arthritis (RA).

A study of the "American College of Rheumatology" shows that the 80,4% of RA patients were RF positive.

REAGENTS

Diluent (R1)	Tris buffer 20 mmol/L, pH 8,2. Preservative.
Latex (R2)	Latex particles coated with human gammaglobulin, pH 7,4. Preservative.
RF-CAL	Calibrator. Recombinant RF protein. The RF concentration is stated on the vial label.

PRECAUTIONS

R1 and R2: H317 - May cause an allergic skin reaction. Contains 2-Methylisothiazol-3(2H)-one (Proclin 950). Follow the precautionary statements given in MSDS and label of the product.

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION

Use RF Calibrator Reference 1107007.

The sensitivity of the assay and the target value of the calibrator have been standardized against the International Reference NIBSC 64/2 (Rheumatoid Arthritis Serum) WHO.

Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

Ready for use

RF Calibrator: Reconstitute (→) with 2,0 mL of distilled water. Mix gently and bring to room temperature for about 10 minutes before use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Do not freeze; frozen latex and diluent could change the functionality of the test.

Reagent deterioration: Presence of particles (R1, R2) and turbidity (R1).

Reconstituted calibrator: Stable for 1 month at 2-8°C or 3 months at -20°C.

ADDITIONAL EQUIPMENT

- SPIN 800 autoanalyzer.
- Laboratory equipment.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolized or lipemic samples.

REFERENCE VALUES

Normal values up to 8 IU/mL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control Sera are recommended before and after testing samples to monitor the performance of manual and automated assay procedures. It should be used SPINREACT Control ASO/CRP/RF Level L (Ref.: 1102114) and Level H (Ref.: 1102115).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN 800 APPLICATION
EDIT PARAMETERS

Test	RF	No.	**
Full name	RF	Print name	RF
Reac. Type	End Point	Direction	Increase
Pri. Wave.	605	Sec. Wave.	
Unit	IU/mL	Decimal	0.1
Reagent Blank	18 -19	React. Time	24 - 25
Vol. Sample	3 ul	Decrease	Sample blank
R1	240 ul	R2	60 ul
R3		R4	

CALIBRATION

Spline	Water	0	Vol aspirated	Vol net	Vol dil
	Calibrator	*0.0625	3	10	150
	Calibrator	*0.125	3	25	175
	Calibrator	*0.25	3	30	90
	Calibrator	*0.495	6	30	90
	Calibrator	(*)Cal value	3	--	--

Vol aspirated	Vol net	Vol dil
3	10	150
3	25	175
3	30	90
6	30	90
3	--	--

CAL AUTO DILUTIONS

The calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

1. Limit detection: Values less than 6 IU/mL give non-reproducible results.

2. Measurement range: 6-160 IU/mL, under the described assay conditions.

Samples with higher concentrations should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit and measurement range depends on the sample to reagent/ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

3. Prozone effect: No prozone effect was detected upon 800 IU/mL.

4. Sensitivity: Δ 3,34 mA. IU/mL.

5. Precision: The reagent has been tested for 20 days, using three different FR concentrations in a EP5-based study.

EP5	CV (%)		
	35,8 IU/mL	78,05 IU/mL	123,26 IU/mL
Total	4,5%	4,1%	5,9%
Within Run	3,3%	2,6%	3,2%
Between Run	1,7%	2,3%	3,4%
Between Day	2,5%	2,1%	3,6%

6. Accuracy: Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics.

41 samples of different concentrations of FR were assayed. The correlation coefficient (r)² was 0,91 and the regression equation y = 1,2042x + 3,1344.

The results of the performance characteristics depend on the analyzer used.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Frederick Wolfe et al. Arthritis and Rheumatism 1991; 34: 951- 960.
2. Robert W Dorner et al. Clinica Chimica Acta 1987; 167: 1-21.
3. Robert H Shmerling et al. The American Journal of Medicine 1991; 91: 528 - 534.
4. Vladimir Muié et al. Scand J Rheumatology 1972; 1: 181 - 187.
5. Paul R et al. Clin Chem 1979; 25/11: 1909 - 1914.
6. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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