### Quantitative determination of Ferritin

#### IVD

Store 2 - 8°C.

#### PRINCIPLE OF THE METHOD

Ferritin-turbilatex is a quantitative turbidimetric test for the measurement of ferritin in human serum or plasma.

Latex particles coated with specific anti-human ferritin are agglutinated when mixed with samples containing ferritin. The agglutination causes an absorbance change, dependent upon the ferritin contents of the sample that can be quantified by comparison from a calibrator of known ferritin concentration.

#### CLINICAL SIGNIFICANCE

Serum ferritin concentration usually reflects body iron stores and is considered one of the most reliable indicators of iron status of patients. Whereas low serum concentrations of ferritin are always indicative of an iron deficiency, elevated concentrations can occur for variety of reasons. Thus, although elevated concentrations often indicate an excessive iron intake, they are also caused by liver disease, chronic inflammation and malignancies. Pregnant women, blood donors, hemodialysis patients, adolescents and children are groups particularly at risk.

#### REAGENTS

<b>Diluent (R1)</b>	Tris buffer 20 mmol/L, pH 8,2. Preservative.
<b>Latex (R2)</b>	Latex particles coated with rabbit IgG anti-human ferritin, pH 8,2. Preservative.
<b>FERR-CAL</b>	Calibrator. Ferritin concentration is stated on the vial label.

#### PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

#### CALIBRATION

Use Ferritin Calibrator Reference 1107042.

The sensitivity of the assay and the target value of the calibrator have been standardized against the 3<sup>rd</sup> International Standard of Ferritin (94/572, 2008 WHO). Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

#### PREPARATION

**Ferritin Calibrator:** Reconstitute (→) with 3,0 mL of distilled water. Mix gently and incubate at room temperature for about 10 minutes before testing.

**Calibration curve:** Prepare the following FERR calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the FERR calibrator by the corresponding factor stated in table below to obtain the concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator FERR(µL)	--	25	50	100	200	400
NaCl 9 g/L (µL)	400	375	350	300	200	--
Dilution Factor	0	0,0625	0,125	0,25	0,5	1,0

#### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

**Reagent deterioration:** Presence of particles (R1, R2) and turbidity (R1).

#### ADDITIONAL EQUIPMENT

- SPIN 800 autoanalyzer.
- Laboratory equipment.

#### SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolyzed or lipemic samples.

#### REFERENCE VALUES

Men: 30 - 220 µg/L.

Women: 20 - 110 µg/L.

Each laboratory should establish its own reference range.

#### QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT Ferritin Control (Ref.: 1107044).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

#### SPIN 800 APPLICATION

EDIT PARAMETERS			
Test	<b>FERR</b>	No.	**
Full name	FERRITIN	Print name	FERR
Reac. Type	End Point	Direction	Increase
Pri. Wave.	546	Sec. Wave.	
Unit	ug/L	Decimal	0.1
Reagent Blank	18 - 19	React. Time	32 - 33
Vol. Sample	18 ul	Decrease	Sample blank
R1	160 ul	R2	40 ul
		R3	
		R4	
<b>CALIBRATION</b>	Spline	Water	0
		Calibrator	*0.0625
		Calibrator	*0.125
		Calibrator	*0.25
		Calibrator	*0.5
		Calibrator (*)	Cal Value

The Calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

#### PERFORMANCE CHARACTERISTICS

**Measuring range:** Up to 600 µg/L. Samples with higher values should be diluted 1/5 in NaCl 9 g/L and retested. The upper linearity limit increases as the sample volume and the sensitivity decrease.

**Detection limit:** 5,04 µg/L.

**Quantification limit:** Values under 6,6 µg/L may give non-reproducible results.

**Prozone effect:** No prozone effect was detected at least up to 9000 µg/L.

**Precision:** According to the EP5-A2 standards (CLSI), the reagent has been tested for 20 days, measuring each level per duplicate twice a day (n=80):

Mean (µg/L)	Intra-assay (n= 80)			Total (n= 80)		
	33,4	114,5	289,8	33,4	114,5	289,8
SD	1,7	1,4	2,4	2,1	3,4	7,5
CV (%)	5,1	1,2	0,8	6,3	2,9	2,6

**Method comparison:** The reagent was compared to another commercially available Ferritin reagent by testing 144 samples (male and female), with concentrations between 6,97 and 730 µg/L. The coefficient of correlation (r) was 0,988, and the equation  $y = 0,96x + 1,15$

Performance characteristics depend on the analyzer used.

#### NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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