Quantitative determination of microalbumin (μALB) IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

Microalbumin-turbilatex is a quantitative turbidimetric test for the measurement of microalbumin (μALB) in human urine.

Latex particles coated with specific antibodies anti-human albumin are agglutinated when mixed with samples containing μALB. The agglutination causes an absorbance change, dependent upon the μALB contents of the patient sample that can be quantified by comparison from a calibrator of known μALB concentration.

CLINICAL SIGNIFICANCE

Microalbuminuria is at present defined as an excretion rate for albumin between 20 and 200 mg/L, which is already above normal values but still below the values seen in patients with "conventional" proteinuria.

Microalbuminuria is a marker of an increased risk of diabetic nephropathy as well as cardiovascular disease in patients with insulin-dependent diabetes mellitus as well as with non-insulin-dependent diabetes mellitus. More recently, microalbuminuria has been found to be associated with cardiovascular disease also in the non-diabetic population. In fact, microalbuminuria may show to be a risk factor of cardiovascular disease among otherwise apparently healthy people.

REAGENTS

Diluent (R1)	Glycine buffer 100 mmol/L, pH 10,0. Preservative.
Latex (R2)	Particles coated goat IgG with anti -human albumin, pH 8,2. Preservative.
μALB-CAL	Liquid Calibrator. Microalbumin concentration is stated on the vial label.

PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

Use Microalbumin Calibrator Reference 1107072.

The sensitivity of the assay and the target value of the calibrator have been standardized against the International Reference Material ERM-DA 470K/IFCC. Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

Ready for use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Reagent deterioration: Presence of particles (R1, R2) and turbidity (R1).

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- SPIN 800 autoanalyzer.
- Laboratory equipment.

SAMPLES

24 hours or random/ first morning urine specimen. It is recommended to adjust the pH at 7.0 with NaOH/HCl 1 mol/L. Stable 7 days at 2-8°C when sodium azide 1 g/L is added to prevent contamination.

Urine should be centrifuged before testing.

REFERENCE VALUES

Normal values up to 30 mg/24 hrs urine specimen and 20 mg/L in a first morning urine specimen.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT Microalbumin Control Ref: 1107073.

BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.

SPIN 800 APPLICATION

EDIT PARAMETERS

Test	mALB	No.	**
Full name	MICROALBUMIN	Print name	mALB
Reac. Type	Fixed Time	Direction	Increase
Pri. Wave.	546	Sec. Wave.	
Unit	mg/L	Decimal	0.1
Reagent Blank	-	React. Time	18 - 25
Vol. Sample	3 ul	Increas.	Decrease
			Sample blank
R1	240 ul	R2	60 ul
		R3	R4

CALIBRATION

Linear two points	Water	0
	Calibrator	*

The Calibration is stable until 25 days. After this period the Calibration must be performed again in order to obtain good results.

PERFORMANCE CHARACTERISTICS

- Linearity limit:** Up to 150 mg/L, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample reagent ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection limit:** Values less than 2 mg/L give non-reproducible results.
- Prozone effect:** No prozone effect was detected up to 1000 mg/L.
- Sensitivity:** Δ 3,8 mA. mg/L.
- Precision:** The reagent has been tested for 20 days, using three different microalbumin concentrations in a EP5-based study.

EP5	CV (%)		
	+/- 10,36 mg/L	+/- 16,95 mg/L	+/- 57,33 mg/L
Total	4,5%	3,1%	2,5%
Within Run	1,9%	1,4%	1,1%
Between Run	4,1%	2,7%	2,3%
Between Day	0,0%	0,0%	0,0%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 49 samples of different concentrations of microalbumin were assayed. The correlation coefficient (r)² was 0.99 and the regression equation y = 0,424x +10,55.

The results of the performance characteristics depend on the analyzer used.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

- Feldt-Rasmussen B et al. J Diab Comp 1994; 8: 137-145.
- Panuyiotou B N. Journal International Medical Research 1994; 22: 181-201.
- Bar J et al. Diabetic Medicine 1995; 12: 649-656.
- Gilbert R E et al. Diabetic Medicine 1994; 11: 636-645.
- Medcalf E A et al. Clin Chem 1990; 36/3: 446-449.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACCPress, 1995.

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