**Quantitative determination of low levels of C-Reactive Protein IVD**

Store 2 - 8°C.

**PRINCIPLE OF THE METHOD**

The CRP-ultrasensitive is a quantitative turbidimetric test for the measurement of low levels of C- reactive protein (CRP) in human serum or plasma. Latex particles coated with specific anti- human CRP are agglutinated when mixed with samples containing CRP. The agglutination causes an absorbance change, dependent upon the CRP contents of the patient sample that can be quantified by comparison from a calibrator of known CRP concentration.

**CLINICAL SIGNIFICANCE**

CRP is an acute-phase protein present in normal serum, which increases significantly after most forms of tissue injuries, bacterial and virus infections, inflammation and malignant neoplasia. CRP may be also useful in detecting atherosclerotic process and providing important prognostic information about patients with asymptomatic heart disease, unstable angina, and myocardial infarction. Recent studies in apparently healthy people show that CRP concentration in serum rise long before traditional symptoms of heart and vascular diseases are noticed.

**REAGENTS**

<b>Diluent-ultra (R1)</b>	Tris buffer 20 mmol/L, pH 8,2. Preservative.
<b>Latex-ultra (R2)</b>	Latex particles coated with goat IgG anti- human CRP, pH 7,3. Preservative.
<b>U-CRP CAL</b>	Liquid Calibrator. C-Reactive protein concentration is stated on the vial label.
<b>Optional</b>	Ref: 43036 CRP Ultra Control.

**PRECAUTIONS**

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

**CALIBRATION**

Use CRP Ultra Calibrator Reference 43035.

The sensitivity of the assay and the target value of the calibrator have been standardized against the Reference Material ERM-DA 474/IFCC.

Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

**PREPARATION**

Ready for use.

**CRP Calibrator:** Ready for use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

**Reagent deterioration:** Presence of particles (R1, R2) and turbidity (R1).

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

**ADDITIONAL EQUIPMENT**

- SPIN 800 autoanalyzer.
- Laboratory equipment..

**SAMPLES**

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

**REFERENCE VALUES**

Below 3 mg/L is considered normal.

Each laboratory should establish its own reference range.

**QUALITY CONTROL**

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT CRP Ultra Control (Ref.:43036). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**BARCODED REAGENTS LOAD MUST BE PRECEDED OF A SPINREACT "DATABASE" COPY INTO THE ANALYZER SOFTWARE. IT IS AVAILABLE UNDER REQUEST TO SPINREACT.**

**SPIN 800 APPLICATION**

EDIT PARAMETERS					
Test	<b>CRPu</b>	No.	**		
Full name	CRPu	Print name	CRPu		
Reac. Type	End Point	Direction	Increase		
Pri. Wave.	546	Sec. Wave.			
Unit	mg/L	Decimal	0.1		
Reagent Blank	-	React. Time	18 - 31		
Vol. Sample	3 ul	Increas.	Decrease	Sample blank	
R1	240 ul	R2	60 ul	R3	R4

CALIBRATION					
Spline	Water	0	Vol aspirated	Vol net	Vol dil
Calibrator	*0.0625		3	10	150
Calibrator	*0.125		3	25	175
Calibrator	*0.25		3	30	90
Calibrator	*0.4951		6	30	90
Calibrator	*		3	--	--

**CAL AUTO DILUTIONS**

The Calibration is stable until **31 days**. After this period the Calibration must be performed again in order to obtain good results.

**PERFORMANCE CHARACTERISTICS**

- Linearity limit:** Up to 10 mg/L, under the described assay conditions. Samples with higher concentrations should be diluted 1/3 in NaCl 9 g/L and re-tested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection limit:** Values less than 0,05 mg/L give non-reproducible results.
- Prozone effect:** No prozone effect was detected upon 800 mg/L.
- Sensitivity:** Δ 44 mA.mg/L.
- Precision:** The reagent has been tested for 20 days, using three different CRP concentrations in a EP5-based study.

EP5	CV (%)		
	+/- 0,28 mg/L	+/- 3,09 mg/L	+/- 5,95 mg/L
Total	7,7%	2,7%	3,0%
Within Run	4,5%	1,7%	1,4%
Between Run	4,7%	1,9%	2,7%
Between Day	4,1%	0,7%	0,0%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 23 samples of different concentrations of CRP were assayed. The correlation coefficient (r) was 0,99 and the regression equation  $y=1,0028x-0,0625$ . The results of the performance characteristics depend on the analyzer used.

**NOTES**

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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