**One Step Immunochromatographic assay to detect total prostate specific antigen (PSA)**

**For professional use only.**

**Read these instructions carefully before using the test.**

**IVD TEST: Store at 2-30°C**

**INTENDED USE**

SPIN-PSA has been developed for measuring an abnormal PSA levels in human serum. The sensitivity of SPIN-PSA test is 4.0 ng/ml. This value was calculated using a serial dilution of PS-Act standard from Scripps in a buffered solution. It is useful for diagnosis because it is one of the few tumor markers whose detection above a specific concentration presents clinical utility in diagnosing prostate cancer, the second most important cancer in men and the most significant in those of advanced age. This assay provides only a preliminary result. Clinical consideration and professional judgment must be applied.

**CLINICAL SIGNIFICANCE**

The prostate specific antigen (PSA) is a glycoprotein of approximately 34 kDa (free conformation). It has a serin-protease enzymatic activity and plays an important role in liquefaction-gelification of the seminal fluid.

**PRINCIPLE OF THE METHOD**

PSA present in the specimen reacts with the colloidal latex particles coated with specific monoclonal antibodies against PSA. This colloidal complex particles-antibodies-PSA migrates by a chromatographic process to the reaction zone. In this zone there are antibodies to the PSA that will react with the complex colloidal latex particles-antibodies-PSA. This reaction originates the formation of a red deposit.

**REAGENTS AND MATERIAL SUPPLIED**

- Ref. 1504020: 5 test devices, each sealed in a pouch with a dropper pipette.
- Ref. 1504021: 20 test devices, each sealed in a pouch with a dropper pipette.

**MATERIALS INCLUDED IN THE KIT**

- Reaction devices (cassette or strip)
  - Pipettes in the case of the Simple PSA format
- BIOSIGMA S.r.l.**  
Via Valletta, 6 – Zona Ind. Cantarana  
30010 – Cona (VE) – ITALIA

**Indication:**

In the simple PSA, the number of pipettes supplied must equal the number of cassettes. (Example: kit of 20 units of containing 20 pipettes). No pipettes are supplied in the PSA stick.

**MATERIAL REQUIRED BUT NOT PROVIDED**

- Timer

**PRECAUTION**

1. All reagents are for *in vitro* use only.
2. Patient specimens (sera) may contain infectious agents and should be handled and disposed of as potential biohazards.
3. Do not interchange kit components from different kit lot numbers.
4. Allow kit components and specimens to reach the room temperature before use, as cold reagents and/or specimens may decrease assay performances. 20-30 minutes are recommended.
5. Do not use kit components beyond labelled expiration date.
6. In case of packaging breakage, the product can be used if no component has been damaged; however, if the aluminium pouch of the device would be damaged, then the test performance can be affected and therefore results can not be reliable. Then we recommend discard this test.
7. It is very important to add the correct quantity of serum. If it is lower than the suggest one, the sample may not arrive to the reaction area and the test would not run. If it is higher, the reactive could be diluted and the line very slight.
8. All product used should be rejected according to the legislation in force.
9. Do not use the test if any coloured line can be seen in the test zone before performing the test.
10. For a single packaged device, once it is opened it must be run immediately.
11. The specificity of the method is very good, but a false positive reaction might occasionally appear. Therefore, as an established general standard, all samples that give a positive result should be checked with another different method having equal or greater sensitivity. As is true in any test using mouse antibodies, there is a possibility of interference by human anti-mouse antibodies (HAMA) or high rheumatoid factor levels.
12. The final diagnosis should not be based on only the result from a test; the diagnosis should be established by the correlation of test results with other appropriate data and with clinical symptoms.
13. Do not discard the external box of the kit until its content has been totally used. The external box contains essential information regarding the EC marking and lotification.

**STORAGE AND STABILITY**

Store the kit at room temperature 2-30°C (59-86°F) (Not recommended storage in the fridge). Each device may be used until the expiration date printed on the label if it remains sealed in its foil pouch.

Do not freeze and /or expose the kit to temperatures over 30°C (86°F). The expiration date is printed on the package.

**SPECIMEN COLLECTION**

Undiluted human serum, fresh and free from turbidity samples must be used. The samples may be stored in the refrigerator between 2-8°C for 1 or 2 days. For longer storage they must be kept frozen at -20°C. Frozen specimen should be totally thawed, brought to room temperature and homogenised by shaking gently before testing.

**ASSAY PROCEDURE**

1. Remove the SPIN-PSA device from its sealed bag after getting room temperature and just before using.
2. Slowly place 4 drops or 125 µl of clear serum onto the round window marked with an arrow (sample window).
3. Read after 5 minutes.

**INTERPRETATION OF RESULTS**


**NEGATIVE:** Only one **BLUE** band appears across the result window close to the letter "C" control.

**POSITIVE:** In addition to the **BLUE** control band, a distinguishable **PINK-RED** band also appears across the result window (close to the letter "T" test). The intensity of the line depends on the concentration of antigen.


**INVALID:** If no blue line appears, the test has not been correctly run, the reactivities have been damaged or the sample added was incorrect. Run another test.

In order to obtain the correct sensitivity with SPIN-PSA and in order to assure negative results the final reading time should be 5 minutes.

**Any line or colour that should appear after 5 minutes has no diagnostic value. The test has proved to be functional at running temperatures between 20-30°C.**

**NEGATIVE**

**POSITIVE**

**QUALITY CONTROL**

This test contains a built-in control feature, the C line. The presence of the C line indicates that an adequate sample volume was used and that the reagents migrated properly. If a C line does not form, the test is considered invalid. In this case, review the whole procedure and repeat the testing with a new device.

**LIMITATIONS OF THE PROCEDURE**

1. The test is qualitative and no quantitative interpretation should be made with respect to the intensity of the positive line, when reporting the result.
2. Test results should be used in conjunction with information available from the patient clinical evaluation and other diagnostic procedures.
3. It is very important to add the correct quantity of serum. If it is lower than the suggest one, the sample may not arrive to the reaction area and the test would not run. If it is higher, the reactive could be diluted and the line would be very slight.
4. It is very important to control the reaction time. If the reaction time is lower than the suggested, samples in the detection limit would not appear. If the reaction time is higher than the suggested, sensitivity would be affected and interpretations could be wrong.
5. Although specificity is high (higher than 98%), it may show some occasional false positive reaction; consequently, as a general practice, all samples giving positive results should be retested using another different method with the same or better sensitivity. As with any assay employing mouse antibodies, the possibility exists for interference by human anti-mouse antibodies (HAMA) or high levels of RF in the sample.

**PERFORMANCE CHARACTERISTICS**

**Analytical Sensitivity**

SPIN-PSA will give a positive result with specimens containing more than 4.0 ng/ml of PSA, as evaluated with the standards of PSA and PSA-ACT of SCRIPPS Laboratories, USA.

The use of human serum as diluent of PSA antigens for control purposes is not recommended because the reached sensitivity varies between different sources.

ng/ml	PSA	PSA-ACT
64	+	+
32	+	+
16	+	+
8	+	+
4	±	±
2	—	—
1	—	—
0	—	—

Four people carried out a test in order to determine the detection limit. PSA standard was diluted in the PBS-BSA buffer and tested in accordance with the kit instructions. There were no significant variations in the assignment of the detection limit; all persons detect 4 ng/ml without any problem. PSA standard was diluted in the PBS-BSA buffer and tested in accordance with the kit instructions with 3 different batches. Variations between batches are very slight and do not affect significantly the detection limit.

**Diagnostic Sensitivity and Especificity**

The Monoclonal Antibodies used to elaborate the SPIN-PSA recognise epitopes present in different PSA commercial preparations, they do not react with any other human serum proteins and show identical specificity for free PSA and PSA-ACT complex (alpha-1-anti-chymotripsin complexed PSA). To determine test diagnostic sensitivities and specificities, various studies have compared the results of this test with those provided by other tests.

The first study focused on specificity. Tests were performed in duplicate on 198 fresh serum samples – the majority negative – from a hospital in Zaragoza (Spain), using an ELISA Seratec PSA test and the Spin PSA serum from lot 1 (test time: 5 minutes). The results obtained were as follows:

	ELISA positive	negative
SPIN-PSA positive	2	3
SPIN-PSA negative	0	193

Analysis of these results showed a specificity for the Simple PSA or Spin PSA > 98%. In the second study, 109 fresh serum samples from an analysis center in Zaragoza (Spain) were tested. The reference test was the Abbott Architect; there were 28 values <sup>3</sup> 4.00 ng/ml (considered positive) and 81 < 4.00 ng/ml (considered negative). The results obtained from comparing the Architect Total PSA (Abbott) and the Spin PSA test, based on cut-off levels, were as follows:

Abbott Cut-off	Spin- PSA			
	Sensitivity	Specificity	PPV	NPV
6	13/13 = 100%	75/96 = 78%	13/34 = 38%	75/75 = 100%
4	25/28 = 89%	72/81 = 89%	25/34 = 73%	72/75 = 96%
2	34/75 = 45%	34/34 = 100%	34/34 = 100%	34/75 = 45%

**PACKAGING**

Ref. 1504020 Ref. 1504021	Cont	5 cards 20 cards
------------------------------	------	---------------------

We can see that the concordance obtained (4 ng/ml) was 89% sensitivity and 89% specificity. We can also say that all the samples with PSA values above 6 ng/ml (according to the reference test, Architect) as positive. Likewise, nearly all the samples that presented values of less than 2 ng/ml were detected as negative.

**INTERFERING SUBSTANCES**

The substances indicated in the table, at the concentration specified, did not produce any interference in the result for either negative or positivized samples.

Substance	Concentration	Result at [PSA]	
		0 ng/ml	8 ng/ml
Acetaminophen	200 mg/L	-	+
Acetylsalicylic acid	200 mg/L	-	+
Ampicillin	200 mg/L	-	+
Ascórbic acid	200 mg/L	-	+
Atropine	200 mg/L	-	+
Caffeine	200 mg/L	-	+
Gentisic acid	200 mg/L	-	+
Phenilpropanolamine	200 mg/L	-	+
Salicylic acid	200 mg/L	-	+
Glucose	20 mg/ml	-	+
Urea	40 mg/ml	-	+
Uric acid	100 ug/ml	-	+
Protein (BSA)	200 mg/ml	-	+
Bilirubin	0.02 mg/ml	-	+
Estriol	0.02 mg/ml	-	+
Pregnanediol	0.02 mg/ml	-	+

**INTRA-ASSAY PRECISION - REPETEABILIY**

Results within experimental error (differences less than a 1/2 dilution) were obtained on assaying, in quintuplicate, a sensitivity curve with the same lot.

**REPRODUCIBILITY**

**INTER-DAY PRECISION**

Using 1 product lot, 10 replicas of the sensitivity curve were performed over 10 consecutive days. Only a difference of a 1/2 dilution was found.

**INTER-LABORATORY PRECISION**

Six different operators assayed those same samples, maintaining elevated precisions and concordances. Only a difference of a 1/2 dilution was found.

**INTER-LOT PRECISION**

A sensitivity curve was performed in parallel with 3 product lots. The analyses were performed by the same person on the same day. Only a difference of a 1/2 dilution was found. *Observed differences are acceptable and tolerable as it is a qualitative immunochromatographic technique, with an inherent variability.*

**HOOK EFFECT**

The highest PSA amounts detected have been 19,000 ng/ml (strongly positive) and 125,000 ng/ml (clearly positive, although with less intensity). Comparing with a typical cut-off (or sensitivity) of 4 ng/ml, it can be deduced that 2500-25000 times more PSA can be detected without any problems.

**RELATION WITH THE WHO STANDARD**

An ELISA assay was performed, which showed that the relation between the evaluations of the internal Spin PSA standard (Scripps PSA-ACT) and WHO standards (free PSA or 90:10) is close to unity (with differences less than 20%); 4 ng/ml of Scripps PSA-ACT was equivalent to 3.7 ng/ml of WHO PSA 90:10 (96/670) and to 3.4 ng/ml of WHO free PSA (96/668).

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For in vitro diagnostic use

Distributed by

This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices. Notified Body 0197

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OPERON, S.A. - Camino del Plano, 19 - E-50410 Cuarte de Huerva - (Zaragoza) - ESPAÑA