One Step for the qualitative detection of Rotavirus-Adenovirus in faeces
IVD
Store at 2-30°C

INTENDED USE

The Spin-Rotavirus/Adenovirus Test is a rapid coloured chromatographic, for health care professional use only.

It is intended for use at point of care facilities for the qualitative detection of Rotavirus and/or Adenovirus in faeces.

This assay provides only a preliminary result. Clinical consideration and professional judgment must be applied, particularly when preliminary positive results are evaluated.

CLINICAL SIGNIFICANCE

Rotavirus and Adenovirus are the major causes of infectious gastroenteritis in infants and young children, also observed in adults. It is transmitted by fecal-oral contact. The main symptoms of viral gastroenteritis are watery diarrhoea and vomiting. The affected person may also have headache, fever, and abdominal cramps ("stomach ache"). In general, the symptoms begin 1 to 2 days following infection with a virus that causes gastroenteritis and may last for 1 to 10 days, depending on which virus causes the illness (Rotavirus 3 days and Adenovirus 5-8 days).

PRINCIPLE OF THE METHOD

This assay is a chromatographic immunoassay. The membrane is pre-coated with mouse monoclonal antibodies, on the test band region, against viral antigens.

During testing, the sample is allowed to react with the coloured conjugate (anti-rotavirus mouse monoclonal antibodies-red microspheres and anti-adenovirus mouse monoclonal antibodies-blue microspheres) which was pre-dried on the test. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured conjugates. One or two coloured lines will be visible, depending upon the virus content of the sample. These lines are used to interpret the result.

The mixture continues to move across the membrane to the immobilized antibody placed in the control band region, a GREEN coloured band always appears. The presence of this GREEN band serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

MATERIALS SUPPLIED

- 20 test devices, each sealed in a pouch with a desiccant.
- 20 Stool collection tubes, each filled with sample extraction diluent.

MATERIAL REQUIRED BUT NOT PROVIDED

- Specimen collection containers
- Timer
- Disposable gloves

STORAGE AND STABILITY


Store as packaged in the sealed pouch at 2-30°C. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

SPECIMEN COLLECTION AND STORAGE

Collect sufficient quantity of faeces (1-2 g or mL for liquid sample). Stool samples should be collected in clean and dry containers (no preservative or transport media). The samples can be stored in the refrigerator (2-4 °C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed, and brought to room temperature before testing.

Specimen preparation (see illustration):


- Take out the top and add 1 mL (30 drops) of the sample diluent in the stool collection tube.
- Use the stick to pick up a little sample. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.
- For liquid stool samples, aspirate the faecal specimen with a dropper and add 100 µL into the testing tube or vial with buffer.

**PRECAUTIONS**

1. The instructions must be followed to obtain accurate results.
2. Appropriate precautions are necessary in the collection, handling of the specimens and used assay materials as potentially biohazardous.
3. Do not use kit beyond the expiration date, which appears on the package label. Do not mix reagents or components from different lots of test kits.
4. The tests should be discarded in a proper biohazard container after testing.

ASSAY PROCEDURE

1. Refrigerated specimens and other test materials, including devices, must be equilibrated to room temperature (**15-30°C**) before testing to avoid invalid results.
2. Use a separate stool collection tube and device for each sample or control. Label them with specimen identification.
3. Proceed to shake the stool collection tube in order to assure good sample dispersion.
4. Cut the end of the top (4).
5. Remove the device from its sealed bag just before using. **Do not open pouches until ready to perform the assay.**
6. Dispense exactly 4 drops or 150 µL into the circular window marked with an arrow, avoiding to add solid particles with the liquid (5).
7. Read the result at **10 minutes** (the coloured bands appear).

(4)

INTERPRETATION OF RESULTS

Adenovirus	Rotavirus	Adenovirus-Rotavirus	NEGATIVE	INVALID

NEGATIVE: only one GREEN band appears across the central window in the site marked with the letter C (control line).

ROTOVIRUS POSITIVE: in addition to the GREEN control band, a RED band (Rotavirus test line) also appears in the site marked with the letter T2.

ADENOVIRUS POSITIVE: in addition to the GREEN control band, a RED band (Adenovirus test line) also appears in the site marked with the letter T1.

ROTOVIRUS-ADENOVIRUS POSITIVE: All the lines above described (a GREEN control band in the control region, two REDS bands in the result region) could appear at the same time during the test performance due to a simultaneous infection of Rotavirus and Adenovirus.

INVALID: A total absence of the control coloured band (GREEN) regardless the appearance or not of the results lines (REDS). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit and contact your local distributor.

Notes

The intensity of the red coloured band(s) in the result line region (T) will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

Any other result different from the ones described, must be considered as INVALID.

QUALITY CONTROL PROCEDURE**• Built-in Control Features**

This test contains a built-in quality control feature, the green line appearing in the control region (C). It confirms sufficient specimen volume and correct procedural technique. A clear background is an internal negative background control. If the test is working properly, the background in the result area should be clear and not interfere with the ability to read the result.

• External Quality Control

External controls are recommended, positive and negative, to monitor the performance of the assay.

LIMITATIONS

1. This test is a qualitative assay for professional *in vitro* diagnostic use only.
2. The test must be carried out within 2 hours of opening the sealed bag.
3. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
4. After one week of infection, the number of viruses in feces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
5. This test provides a presumptive diagnosis for Rotavirus infections. A confirmed infection diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.

EXPECTED VALUES

Negative results are expected in healthy infants and young children, also in healthy adults.

PERFORMANCE

The evaluation was conducted comparing the results obtained using the SPIN-ROTA+ADENOVIRUS to another commercial available Rotavirus and Adenovirus membrane assay.

Sensitivity

The detection of Rotavirus showed >99% of concordance in sensitivity.

The detection of Adenovirus showed a 90% of concordance in sensitivity.

Specificity

The detection of Rotavirus showed a 98% of concordance in specificity.

The detection of Adenovirus showed >99% of concordance in specificity.

The use of mouse monoclonal antibodies in the elaboration of SPIN-ROTA+ADENOVIRUS assures high degree of specificity for the detection of this virus.

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