
**Quantitative determination of lipase IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

The pancreatic lipase in presence of colipase, desoxycholate and calcium ions, hydrolyses the substrate 1-2-O-dilauryl-rac-glycero-3-glutaric acid-(6'-methylresorufin)-ester. The sequence of reactions involved in the enzymatic direct lipase determination is the following:

The rate of methylresorufin formation, measured photometrically, is proportional to the catalytic concentration of lipase present in the sample.

**CLINICAL SIGNIFICANCE**

Lipase (LPS) is a pancreatic enzyme necessary for the absorption and digestion of nutrients that catalyzes the hydrolysis of glycerol esters of fatty acids. Determination of LPS is used for diagnosis of diseases of pancreas such as acute and chronic pancreatitis and obstruction of the pancreatic duct<sup>1,7,8</sup>. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b> Buffer	TRIS pH 8.3	40 mmol/L
	Colipase	≥ 1 mg/L
	Desoxycholate	1,8 mmol/L
	Taurodesoxycholate	7,2 mmol/L
<b>R 2</b> Substrate (micro-emulsion)	Tartrate pH 4,0	15 mmol/L
	Lipase Substrate	≥ 0,7 mmol/L
	Calcium chloride (CaCl <sub>2</sub> )	0,1 mmol/L

**PREPARATION**

The reagent is ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 580 nm ≥ 1,00.
- R 2 is a turbid orange-colored micro-emulsion, discard if turning to red.

**ADDITIONAL EQUIPMENT**

- MINDRAY BS-120 / BS-200E Autoanalyzer.
- General laboratory equipment <sup>(Note 2)</sup>.

**SAMPLES**

Serum or plasma with sodium citrate, EDTA or heparin<sup>1</sup>.  
Avoid repeated frozen and unfrozen. Stability: 2 days at 2-8°C.

**REFERENCE VALUES<sup>1</sup>**

≤ 38 U/L (U/L methylresorufin at 37°C).  
These values are for orientation purpose; each laboratory should establish its own reference range.

**QUALITY CONTROL**

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**NOTES**

- In some storage conditions (i.e. storage at a temperature lower than the one indicate) a precipitate may appear in the vial that will not influence that the reagent performance; however, it is recommended to resuspend the product with a slight rotation.
- In order to avoid contamination it is recommended to use disposable material.

**MINDRAY BS-120 / BS-200E APPLICATION**

PARAMETERS			
Test	LIPA / LIPA	R1	300 / 300
Nº	**	R2	60 / 60
Full Name	LIPA / LIPA	Sample volume	3 / 3
Standard Nº		R1 Blank	
Reac. Type	Kinet / Kinet	Mixed Rgt Blank	
Pri. Wavelength	578 / 570	Linearity Range	1.0 U/L 250.0 U/L
Sec. Wavelength		Linearity Limit	*
Direction	Increase / Increase	Substrate Limit	*
Reac. Time	4_11 / 4_11	Factor	*
Incuba. Time		Prozone check	*
Units	U/L / U/L	q1	q2
Precision	0.1 / 0.1	q3	q4
		PC	Abs
CALIBRATION (Cal + Rgt Blk)			
Rule	One-point Linear / Two-point Linear		
Sensitivity	1 / 1		
Replicates	2 / 2		
Interval (days)	0 / 0		
Difference Limit			
SD			
Blank Response			
Error Limit			
Correlation Coefficient			

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **35 days**. After this period the blank parameter must be performed again in order to validate the calibration.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From detection limit of 5 U/L to linearity limit of 250 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (U/L)	40,2	59,35	38,5	58,9
SD	0,410	0,875	1,10	1,25
CV (%)	1,02	1,47	2,86	2,13

**Sensitivity:** 1 U/L= 0,00059792 (A)

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 101 samples were the following:

Correlation coefficient (r)<sup>2</sup>:0,99732.

Regression equation: y= 0,50054x + 3,9443.

The results of the performance characteristics depend on the analyzer used.

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