# Alkaline phosphatase

p-Nitrophenylphosphate. Kinetic. AMP buffer (IFCC)


## Quantitative determination of alkaline phosphatase (ALP) IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Kinetic photometric test, according to the International Federation of Clinical Chemistry and Laboratory Medicines (IFCC)<sup>7</sup>.

Alkaline phosphatase (ALP) catalyses the transfer of the phosphate group from p-nitrophenylphosphate to 2-amino-2-methyl-1-propanol (AMP), liberating p-nitrophenol according to the following reaction:

The rate of p-Nitrophenol formation, measured photometrically, is proportional to the catalytic concentration of alkaline phosphatase present in the sample<sup>1,2</sup>.

### CLINICAL SIGNIFICANCE

Alkaline phosphatase is an enzyme present in almost all weaves of the organism, being particularly high in bone, liver, placenta, intestine and kidney.

Both increases and decreases of plasma ALP are of importance clinically. Causes of increased plasma ALP: Paget's disease of bone, obstructive liver disease, hepatitis, hepatotoxicity caused by drugs or osteomalacia.

Causes of decreased plasma ALP: Cretinism and vitamin C deficiency<sup>1,5,6</sup>. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

|                  |  |            |
|------------------|--|------------|
| R 1<br>Buffer    | 2-Amino-2-methyl-1-propanol                        | 0,35 mol/L |
|                  | Zinc sulfate                                       | 1 mmol/L   |
|                  | Magnesium acetate                                  | 2 mmol/L   |
|                  | N-hydroxyethylethylenediaminetriacetic acid (EDTA) | 2 mmol/L   |
| R 2<br>Substrate | p-Nitrophenylphosphate (pNPP)                      | 10 mmol/L  |

### PREPARATION

The reagent is ready to use.

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not freeze the reagents.

Do not use reagents over the expiration date.

### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm > 1,50.

### ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200E Autoanalyzer.
- General laboratory equipment.

### SAMPLES

Serum or heparinized plasma<sup>1</sup>. Use unhemolyzed serum, separated from the clot as soon as possible. Stability: 3 days at 2-8°C.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### REFERENCE VALUES<sup>1</sup>

|        | 25°C        | 30°C        | 37°C         |
|--------|-------------|-------------|--------------|
| Adults | 17 - 77 U/L | 21 - 94 U/L | 26 - 117 U/L |

Factors affecting ALP activities in a normal population include exercise, periods of repaid growth in children and pregnancy.

These values are for orientation purpose; each laboratory should establish its own reference range.

## MINDRAY BS-120 / BS-200E APPLICATION

### PARAMETERS

|                 |                     |                 |                |
|-----------------|---------------------|-----------------|----------------|
| Test            | ALP / ALP           | R1              | 240 / 240      |
| Nº              | **                  | R2              | 60 / 60        |
| Full Name       | ALP / ALP           | Sample volume   | 6 / 6          |
| Standard Nº     |                     | R1 Blank        |                |
| Reac. Type      | Kinet / Kinet       | Mixed Rgt Blank |                |
| Pri. Wavelength | 405 / 412           | Linearity Range | 4 U/L 1400 U/L |
| Sec. Wavelength |                     | Linearity Limit | *              |
| Direction       | Increase / Increase | Substrate Limit | *              |
| Reac. Time      | 2_10 / 2_10         | Factor          | *              |
| Incuba. Time    |                     | Prozone check   | *              |
| Units           | U/L / U/L           | q1              | q2             |
| Precision       | Interger / Interger | q3              | q4             |
|                 |                     | PC              | Abs            |

### CALIBRATION (Cal + Rgt Blk)

|                         |                                     |
|-------------------------|-------------------------------------|
| Rule                    | One-point Linear / Two-point Linear |
| Sensitivity             | 1 / 1                               |
| Replicates              | 2 / 2                               |
| Interval (days)         | 0 / 0                               |
| Difference Limit        |                                     |
| SD                      |                                     |
| Blank Response          |                                     |
| Error Limit             |                                     |
| Correlation Coefficient |                                     |

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **7 days**. After this period the blank parameter must be performed again in order to validate the calibration.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From detection limit of 1,307 U/L to linearity limit of 1400 U/L. If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

### Precision:

|            | Intra-assay (n=20) |      | Inter-assay (n=20) |      |
|------------|--------------------|------|--------------------|------|
|            | Mean (U/L)         | SD   | CV (%)             |      |
| Mean (U/L) | 73                 | 1,67 | 2,27               | 78   |
| SD         | 194                | 3,03 | 1,58               | 209  |
| CV (%)     | 2,27               | 1,86 | 2,13               | 4,90 |
|            |                    |      |                    | 2,34 |

**Sensitivity:** 1 U/L = 0.0004 ΔA / min.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,98929.

Regression equation: y= 2,214x + 2,131.

The results of the performance characteristics depend on the analyzer used.

### NOTES

**SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

### BIBLIOGRAPHY

1. Wenger C. et al. Alkaline phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1094-1098.
2. Rosalki S et al. Clin Chem 1993; 39/4: 648-652.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. IFCC methods for the measurement of catalytic concentration of enzymes. J.Clin.Chem.Clin.Biochem. 1983; 21: 731-748.

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### BIBLIOGRAPHIE

- Wenger C. et al. Alkaline phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1094-1098.
- Rosalki S et al. Clin Chem 1993; 39/4: 648-652.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995. IFCC methods for the measurement of catalytic concentration of enzymes. J.Clin.Chem.Clin.Biochem. 1983; 21: 731-748.

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### BIBLIOGRAFIA

- Wenger C. et al. Alkaline phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1094-1098.
- Rosalki S et al. Clin Chem 1993; 39/4: 648-652.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
- IFCC methods for the measurement of catalytic concentration of enzymes. J.Clin.Chem.Clin.Biochem. 1983; 21: 731-748.

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Ref: MI41245

Cont.

R1:

R2:

5 x 25 mL

2 x 16 mL