CK-MB-LQ (Creatine kinase – MB)

Anti CK-M. Immunoinhibition. Kinetic UV. Liquid

Quantitative determination of creatine kinase MB (CK-MB)

IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The procedure involves measurement of CK activity in the presence of an antibody to CK-M monomer. This antibody completely inhibits the activity of CK-MM and half of the activity of CK-MB while not affecting the B subunit activity of CK-MB and CK-BB. Then it's used the CK method to quantitatively determine CK-B activity^{1,2}. The CK-MB activity is obtained by multiplying the CK-B activity by two.

CLINICAL SIGNIFICANCE

CK-MB is an enzyme formed by the association of two subunits from muscle (M) and nerve cells (B). CK-MB is usually present in serum at low concentration; it is increased after an acute infarct of myocardium and later descends at normal levels. Also is increased, rarely, in skeletal muscle damage^{5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Imidazol, pH 6.7	125 mmol/L
	D-Glucose	25 mmol/L
	N-Acetyl-L-Cysteine	25 mmol/L
	Magnesium acetate	12.5 mmol/L
	NADP	2.52 mmol/L
	EDTA	2.02 mmol/L
Hexokinase		≥6 800 U/L
Anti-human polyclonal CK-M antibody (sheep) sufficient to inhibit up to 2 000 U/L of CK-MM		
R 2	ADP	15.2 mmol/L
	AMP	25 mmol/L
	di-Adenosine-5- pentaphosphate	103 mmol/L
	Glucose-6-phosphate dehydrogenase	≥8 800 U/L
	Creatine phosphate	250 mmol/L

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm ≥ 1,20.

ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200E Autoanalyzer.
- General laboratory equipment.

SAMPLES

Serum free of hemolysis or heparin plasma¹: Stability 7 days at 2-8°C, protected from light.

CK-MB activity decreases a 10% after 24 hours at 4°C or 1 hour at 25°C.

Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1.00	1.53	2.38
30°C	0.65	1.00	1.56
37°C	0.42	0.64	1.00

REFERENCE VALUES

The suspicion of myocardial damage is based on the three following factors:

CK-MB	25°C	30°C	37°C
	> 10 U/L	> 15 U/L	> 24 U/L
TOTAL CK	25°C	30°C	37°C
Men, up to	80 U/L	130 U/L	195 U/L
Women, up to	70 U/L	110 U/L	170 U/L

$$\frac{\text{CK-MB Activity}}{\text{CK Total Activity}} \times 100: 6-25\% \text{ CK-MB Activity in the sample}$$

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

CK-Nac/CK-MB specific control sera (Ref. 1002260) are recommended to monitor the performance of assay procedures.

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

MINDRAY BS-120 / BS-200E APPLICATION

PARAMETERS				
Test	CKMB / CKMB	R1	240 / 240	
Nº	**	R2	60 / 60	
Full Name	CKMB / CKMB	Sample volume	12 / 12	
Standard Nº		R1 Blank		
Reac. Type	Fixed T / Fixed T	Mixed Rgt Blank		
Pri. Wavelength	340 / 340	Linearity Range	2 U/L	1000 U/L
Sec. Wavelength		Linearity Limit	*	
Direction	Increase / Increase	Substrate Limit	*	
Reac. Time	4_17 / 5_17	Factor	*	
Incuba. Time		Prozone check	*	
Units	U/L / U/L	q1	q2	
Precision	Interger / Interger	q3	q4	
		PC	Abs	
CALIBRATION (Cal + Rgt Blk)				
Rule	One-point Linear / Two-point Linear			
Sensitivity	1 / 1			
Replicates	2 / 2			
Interval (days)	0 / 0			
Difference Limit				
SD				
Blank Response				
Error Limit				
Correlation Coefficient				

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **35 days**. After this period the blank parameter must be performed again in order to validate the calibration.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 1,9 U/L to linearity limit of 1000 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/1 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay		Inter-assay	
	Mean (U/L)	SD	CV (%)	
Mean (U/L)	33,7	166,5	31,3	161,0
SD	1,00	3,76	1,19	3,47
CV (%)	2,96	2,26	3,81	2,15

Sensitivity: 1 U/L = 0,000134 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained were the following:

Correlation coefficient (r)²: 0,999.

Regression equation: y = 0,976 x - 0,269.

The results of the performance characteristics depend on the analyzer used.

LIMITATION OF THE PROCEDURE

- The method will also measure any CK-BB isoenzyme present in serum. The activity of the isoenzyme is usually negligible, however, if a significant amount of CK-BB activity is present the CK-MB activity will be overestimated.
- A macro form of BB (immunoglobulin complexed) has been observed which will be measured as B in the assay. If the measured CK-B activity exceeds 20% of the total CK activity, the presence of macro BB should be suspected.

BIBLIOGRAPHY

- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
- Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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CALIBRATION (Cal + Rgt Blk)

Rule	One-point Linear / Two-point Linear
Sensitivity	1 / 1
Replicates	2 / 2
Interval (days)	0 / 0
Difference Limit	
SD	
Blank Response	
Error Limit	
Correlation Coefficient	

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BIBLIOGRAPHIE

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W. et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: MI41254	Cont.	R1: 5 x 25 mL
		R2: 1 x 32 mL

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BIBLIOGRAFIA

- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
- Gerhardt W. et al. Creatine kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: MI41254

Cont.

R1: 5 x 25 mL

R2: 1 x 32 mL