Quantitative determination of Transferrin (TRF)

IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

Anti-transferrin antibodies when mixed with samples containing TRF, form insoluble complexes. These complexes cause an absorbance change, dependent upon the TRF concentration of the patient sample, that can be quantified by comparison from a calibrator of known TRF concentration.

CLINICAL SIGNIFICANCE

Transferrin is a plasma protein that contains a single polypeptide chain with approximately 6% carbohydrate. It is synthesized in the liver and transfers iron through the serum.

Evaluation of plasma TRF levels is useful for the differential diagnosis of anemia and for monitoring its treatment. In the hypochromic anemia (iron deficiency), the TRF level is increased. On the other hand, if the anemia is due to a failure to incorporate iron into erythrocytes, the TRF level is normal or low but the protein is highly saturated with iron. In iron overload, the TRF concentration is normal but saturation exceeds 55% and may be as great as 90%.

TRF concentration may, in fact, be used for assessing nutritional status. In the congenital defect atransferrinemia, very low level of TRF is accompanied by iron overload and a severe hypochromic anemia that is resistant to iron therapy. High levels of TRF occur in pregnancy and during estrogen administration.

REAGENTS

R 1 Diluent	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Preservative.
R 2 Antibody	Goat serum, anti-human transferrin, pH 7.5. Preservative.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

The assay has been standardized against the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL Calibrator dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the transferrin calibrator by the corresponding factor stated in table below to obtain the transferrin concentration of each dilution:

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	6,25	12,5	25	50	100
NaCl 9 g/L (µL)	100	93,75	87,5	75	50	-
Factor	0	0,0625	0,125	0,25	0,5	1,0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200E autoanalyzer.
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C. The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Between 200 - 360 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

MINDRAY BS-120 / BS-200E APPLICATION

PARAMETERS

Test	TRF / TRF	R1	360 / 240
Nº	**	R2	90 / 60
Full Name	TRF / TRF	Sample volume	3 / 2
Standard Nº	6 / 6	R1 Blank	
Reac. Type	Endpoint / Endpoint	Mixed Rgt Blank	
Pri. Wavelength	340 / 340	Linearity Range	*
Sec. Wavelength		Linearity Limit	750 mg/dL
Direction	Increase / Increase	Substrate Limit	*
Reac. Time	-1_7 / -1_10	Factor	
Incuba. Time		Prozone check	*
Units	mg/dL / mg/dL	q1	q2
Precision	Interger / Interger	q3	q4
		PC	Abs

CALIBRATION (Cal + Rgt Blk)

BS120 (MANUAL) / BS200E (AUTODILUTION)

Rule	Spline / Spline
Sensitivity	1 / 1
Replicates	2 / 2
Interval (days)	0 / 0

BS200E CALIBRATION (DILUTION) → 5 CAL LEVELS + 1 WATER LEVEL

Nº CAL DIL	CONCENTRATION	DIL SAMPLE	DIL VOL	SAMPLE VOL
1	CAL *0.0625	9.5	115	2.0
2	CAL *0.125	19.0	115	2.0
3	CAL *0.25	38.0	115	2.0
4	CAL *0.5	38.0	115	4.0
5	CAL *1			2.0

*Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **35 days**. After this period the blank parameter must be performed again in order to validate the calibration.*

PERFORMANCE CHARACTERISTICS

1. Measurement range: Up to the calibrator value (aprox. 800 mg/dL) under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

2. Detection Limit: 0 mg/dL

3. Prozone effect: No prozone effect was detected upon 2000 mg/dL

4. Sensitivity: Δ 3,0 mA / mg/dL (94 mg/dL).

5. Precision: The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

EP5	CV (%)		
	77,02 mg/dL	206,99 mg/dL	377 mg/dL
Total	5,4%	2,5%	5,4%
Within Run	1%	0,8%	1,2%
Between Run	1,7%	1,3%	2,1%
Between Day	5%	2%	4,9%

6. Accuracy: Results obtained using this reagent (y) were compared to those obtained using the Immage from Beckman. 100 samples ranging from 50 to 700 mg/dL of TRF were assayed. The correlation coefficient (r) was 0,95 and the regression equation y = 1,046x +3,843.

The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

1. McNeely Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Kreutzer HJH. J Clin Chem Clin Biochem 1976; 14: 401-406
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAFIA

- McNeely Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Kreutzer HJH. J Clin Chem Clin Biochem 1976; 14: 401-406
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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BIBLIOGRAPHIE

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Kreutzer HJH. J Clin Chem Clin Biochem 1976; 14: 401-406
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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