Antithrombin III

Turbidimetry

Quantitative determination of Antithrombin III (ATHROM-III) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Anti-antithrombin III antibodies when mixed with samples containing antithrombin-III, form insoluble complexes. These complexes cause an absorbance change, dependent upon the antithrombin-III concentration of the patient sample, that can be quantified by comparison from a calibrator of known antithrombin-III concentration.

CLINICAL SIGNIFICANCE

Antithrombin III is a protein synthesized in the liver, normally present in the human plasma. It is the major inhibitor of the thrombin, and inhibits coagulation and limits the forming of blood clots. Antithrombin-III is also capable of activating other components of the coagulation cascade (eg, factor Xa), as well as plasmin. Antithrombin-III deficiency can cause or lead to thrombosis, a clot forming in a blood vessel. Clots forming in the legs and pulmonary embolism are most commonly reported. Antithrombin-III deficiency is usually inherited and affects males and females equally. All family members should be tested if there is history of the disease.

Acquired antithrombin-III deficiency can occur as a result of other conditions. It has been reported in patients with liver diseases, patients receiving certain kinds of chemotherapy, and patients using oral contraceptives.

REAGENTS

R 1 Diluent	Tris buffer 20 mmol/L, PEG 8000, pH, 8,3. Preservative.
R 2 Antibody	Goat serum, anti-human antithrombin III, pH 7,5. Preservative.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

It must be used the PROT CAL to calibrate the reagent.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the ATHROM-III calibrator by the corresponding factor stated in table below to obtain the antithrombin-III concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (µL)	--	10	25	50	75	100
NaCl 9 g/L (µL)	100	90	75	50	25	-
Factor	0	0,1	0,25	0,5	0,75	1,0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity.

Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200E autoanalyzer.
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. Sodium citrate should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C. Do not use highly hemolyzed or lipemic samples.

REFERENCE VALUES

Between 17 - 30 mg/dL. Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

MINDRAY BS-120 / BS-200E APPLICATION

PARAMETERS

Test	AT III / AT III	R1	240 / 240
Nº	**	R2	60 / 60
Full Name	AT III / AT III	Sample volume	6 / 6
Standard Nº	6 / 6	R1 Blank	
Reac. Type	Endpoint / Endpoint	Mixed Rgt Blank	
Pri. Wavelength	340 / 340	Linearity Range	7mg/dL 70mg/dL
Sec. Wavelength		Linearity Limit	*
Direction	Increase / Increase	Substrate Limit	*
Reac. Time	1_17 / -1_17	Factor	
Incuba. Time		Prozone check	*
Units	mg/dL / mg/dL	q1	q2
Precision	0.01 / 0.01	q3	q4
		PC	Abs

CALIBRATION (Cal + Rgt Blk)

BS120 (MANUAL) / BS200E (AUTODILUTION)

Rule	Spline / Spline
Sensitivity	1 / 1
Replicates	2 / 2
Interval (days)	0 / 0

BS200E CALIBRATION (DILUTION) → 5 CAL LEVELS + 1 WATER LEVEL

Nº CAL DIL	CONCENTRATION	DIL SAMPLE	DIL VOL	SAMPLE VOL
1	CAL *0.1	13.0	115	6.0
2	CAL *0.25	38.0	115	6.0
3	CAL *0.5	38.0	115	12.0
4	CAL *0.75	38.0	115	18.0
5	CAL *1			6.0

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **35 days**. After this period the blank parameter must be performed again in order to validate the calibration.

PERFORMANCE CHARACTERISTICS

- 1. Measurement range:** Up to 70 mg/dL (Nota 1), under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.
- 2. Limit detection:** Values less than 7,4 mg/dL give non-reproducible results.
- 3. Prozone effect:** No prozone effect was detected upon 200 mg/dL.
- 4. Sensitivity:** Δ 7,5 mA / mg/dL.
- 5. Precision:** The reagent has been tested for 20 days, using two levels of serum in a EP5-based study.

EP5	CV (%)	
	17,85 mg/dL	35,93 mg/dL
Total	3,2%	3,5%
Within Run	0,8%	1%
Between Run	2,4%	2,4%
Between Day	2%	2,3%

- 6. Accuracy:** Results obtained using this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

The results of the performance characteristics depend on the used analyzer.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
3. Buller HR et al. Critical care Medicine 1982; 10: 311.
4. Kauffman et al. Am J Med 1978; 65: 607.
5. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
6. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

PACKAGING

Ref: MI1102014	Cont.	R1.Diluent:	2 x 30 mL
		R2.Antibody:	1 x 15 mL

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BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Buller HR et al. Critical care Medicine 1982; 10: 311.
- Kauffman et al. Am J Med 1978; 65: 607.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Pres, 1997.

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