**Quantitative determination of creatinine IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

The assay is based on the reaction of creatinine with sodium picrate as described by Jaffé.

Creatinine reacts with alkaline picrate forming a red complex. The time interval chosen for measurements avoids interferences from other serum constituents.

 The intensity of the color formed is proportional to the creatinine concentration in the sample<sup>1</sup>.

**CLINICAL SIGNIFICANCE**

Creatinine is the result of the degradation of the creatine, component of muscles; it can be transformed into ATP, which is a source of high energy for the cells. The creatinine production depends on the modification of the muscular mass, and it varies little and the levels usually are very stable.

Is excreted by the kidneys. With progressive renal insufficiency there is retention in blood of urea, creatinine and uric acid.

 Elevate creatinine level may be indicative of renal insufficiency<sup>1,4,5</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b> Picric Reagent	Picric acid	17,5 mmol/L
<b>R 2</b> Alkaline Reagent	Sodium hydroxide	0,29 mol/L

**PRECAUTIONS**

R2: H314-Causes severe skin burns and eye damage.

Follow the precautionary statements given in MSDS and label of the product.

**PREPARATION**

All the reagents are ready to use.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm  $\geq$  1,80.

**ADDITIONAL EQUIPMENT**

- MINDRAY BS-120 / BS-200E Autoanalyzer.
- General laboratory equipment.

**SAMPLES**

- Serum or heparinized plasma<sup>1</sup>.

Creatinine stability: 24 hours at 2-8°C.

- Urine (24 h)<sup>1</sup>: Dilute sample 1/50 with distilled water. Mix. Multiply results by 50 (dilution factor);

Creatinine stability: 7 days at 2-8°C.

**REFERENCE VALUES<sup>1</sup>**

Serum or plasma:

Male	0,7 - 1,4 mg/dL	$\cong$ 61,8 – 123,7 $\mu$ mol/L
Female	0,6 - 1,1 mg/dL	$\cong$ 53,0 – 97,2 $\mu$ mol/L

Urine: 15-25 mg/Kg/24 h

Male	10 - 20 mg/Kg/24 h
Female	8 – 18 mg/Kg/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

**QUALITY CONTROL**

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINROL H Calibrator, SPINROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**MINDRAY BS-120 / BS-200E APPLICATION**

PARAMETERS			
Test	CREA / CREA	R1	180 / 180
Nº	**	R2	180 / 180
Full Name	CREA / CREA	Sample volume	36 / 36
Standard Nº		R1 Blank	
Reac. Type	Fixed T / Fixed T	Mixed Rgt Blank	
Pri. Wavelength	510 / 505	Linearity Range	0.20 mg/dL 15.00 mg/dL
Sec. Wavelength		Linearity Limit	*
Direction	Increase / Increase	Substrate Limit	*
Reac. Time	2_6 / 2_6	Factor	*
Incuba. Time		Prozone check	*
Units	mg/dL / mg/dL	q1	q2
Precision	0.01 / 0.01	q3	q4
		PC	Abs
CALIBRATION (Cal + Rgt Blk)			
Rule	One-point Linear / Two-point Linear		
Sensitivity	1 / 1		
Replicates	2 / 2		
Interval (days)	0 / 0		
Difference Limit			
SD			
Blank Response			
Error Limit			
Correlation Coefficient			

Blank parameter must be performed in order to get good results in CALIB screen from main menu. The blank calibration is stable until **5 days**. After this period the blank parameter must be performed again in order to validate the calibration.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** From detection limit of 0,000 mg/dL to linearity limit of 35 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	0,92	3,43	0,96	3,50
SD	0,03	0,07	0,04	0,09
CV (%)	2,76	1,90	3,97	2,51

**Sensitivity:** 1 mg/dL = 0,0407  $\Delta$ Abs/min .

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,99584

Regression equation: y = 0,953x + 0,075

The results of the performance characteristics depend on the analyzer used.

**NOTES**

1. Calibration with the aqueous Standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
2. Use clean disposable pipette tips for its dispensation.

**BIBLIOGRAPHY**

1. Murray R.L. Creatinine. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1261-1266 and 418.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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