Quantitative determination of calprotectin in feces IVD

Store at 2 – 8 °C

PRINCIPLE OF THE METHOD

Calprotectin-Turbilatex is a quantitative turbidimetric test for the measurement of human calprotectin in feces.

Latex particles coated with specific anti-human calprotectin are agglutinated when mixed with samples containing calprotectin. The agglutination causes an absorbance change, dependent upon the calprotectin contents in the sample that can be quantified by comparison from a calibrator of known calprotectin concentration.

CLINICAL SIGNIFICANCE

Inflammatory Bowel Disease (IBD) includes Crohn's disease and ulcerative colitis. There is an overlap in symptoms between IBD and Irritable Bowel Syndrome (IBS). Calprotectin-turbilatex allows to rule out IBD in patients under suspicion of IBS. Moreover, it is useful to monitor treatment in patients with IBD ^{1,2,3,4}.

REAGENTS

Diluent (R1)	Buffered solution. Preservative.
Latex (R2)	Latex particles coated with polyclonal IgY antibody anti-human calprotectin. Preservative.
Related products	Control L and Control H Ref.: 1107093 Calprotectin Calibrator Ref.: 1107092 Sample collection tubes Ref.: 1107195

PRECAUTIONS

R1 and R2: H317 - May cause an allergic skin reaction. Contains 2-Methylisothiazol-3(2H)-one (ProClin 950). Follow the precautionary statements given in MSDS and label of the product.

Contains products of animal origin. Handle cautiously as potentially infectious.

CALIBRATION

Use Calprotectin Calibrator Ref. 1107092. Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

Mix the R2 reagent gently before use. The calibrators and controls are ready to be used. The concentrations of calibrators and controls are indicated on the vial label.

STORAGE AND STABILITY

All the components of the kit are stable until the expiry date on the label when stored tightly closed at 2 - 8 °C, protected from light and contaminations prevented during their use. Reagents should not be left inside the analyser after use, they must be stored refrigerated at 2 – 8 °C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date. Do not freeze. Reagent deterioration: presence of particles (R1, R2) and turbidity (R1).

Calibrators/controls stability: stable until the expiry date indicated on the vial label when the vial is kept tightly closed at 2 - 8 °C.

ADDITIONAL EQUIPMENT

- MINDRAY BS-120 / BS-200E autoanalyser.
- Laboratory equipment.

SAMPLES

Use only Spinreact's sample collection tube (ref. 1107195, 1 x 100 tubes), following the instructions included in the package. Store the collection tube with fecal samples at 2 – 30 °C during no more than 28 days.

REFERENCE VALUES

Each laboratory should determine its own normal reference value. A value <50 µg/g might be considered negative, between 50 and 100 µg/g weakly positive, and >100 µg/g highly positive ³.

QUALITY CONTROL

Control solutions are recommended to monitor the performance of assay procedures (Control L and Control H). If control values are found outside the defined range, check the instrument, reagents and technique for problems. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

INTERFERENCES

No significant interferences were detected with bilirubin (<0.05 mg/dL), hemoglobin (<20 mg/dL), Intralipid (<100 mg/dL), ascorbic acid (<10 mg/dL), BaSO₄ (<5 mg/dL), azathioprine (<0.7 mg/dL), mesalazine (<100 mg/dL), prednisolone (<1.5 mg/dL).

MINDRAY BS-120 / BS-200E APPLICATION

PARAMETERS

Test	CALPRO	R1	180
No.	**	R2	60
Full Name	CALPRO	Sample volume	6
Standard No.	6	R1 Blank	
Reac. Type	Fixed Time	Mixed Rgt Blank	
Pri. Wavelength	510/505 nm	Linearity Range	*
Sec. Wavelength		Linearity Limit	*
Direction	Increase	Substrate Limit	*
Reac. Time	1 13	Factor	
Incuba. Time		Prozone check	*
Unit	µg/g	q1	q2
Precision	0.01	q3	q4
		PC	Abs

CALIBRATION (Cal + Rgt Blk)

Rule	Spline
Sensitivity	1
Replicates	2
Interval (days)	0
Difference Limit	
SD	
Blank Response	
Error Limit	
Correlation Coefficient	

*With this parameter the blank must be performed in order to get good results in CALIB screen from main menu. The calibration is stable until **26 days**. After this period a new calibration must be performed.*

PERFORMANCE CHARACTERISTICS

- Linearity limit:** up to the maximum value of the calibrator (approx. 2 000 µg/g), under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in sample buffer (obtained by emptying new collection tubes) and re-tested. The linearity limit depends on the sample/reagent ratio and the analyzer used.
- Quantification limit:** values < 22.6 µg/g may give non-reproducible results.
- Prozone effect:** no prozone effect was detected up to 3400 µg/g.
- Precision:** two calprotectin levels were tested on 19 consecutive working days, two runs per day and two replicates per run:

	% CV (n= 76)	
	53 µg/g	201 µg/g
Between Run	3.1 %	2.1 %
Repeatability	4 %	1.8 %
Total	5.3 %	2.8 %

- Accuracy:** 150 fecal samples were tested comparing with Liaison-calprotectin. The regression curve is fitted to the lineal equation $y = 0.703x + 18.4$ with a regression coefficient $r = 0.815$.

The results of the performance characteristics depend on the analyzer used.

NOTES

Calprotectin levels in pediatric patients are normally higher. The cutoff of 212 µg/g was proposed ⁵.

The clinical diagnosis should not be based on a single test result; the complete medical history of the patient should be taken into consideration.

BIBLIOGRAPHY

- von Roon AC. et al., Diagnostic precision of fecal calprotectin for IBD and colorectal malignancy. Am. J. Gastroenterol 2007.
- van Rhee PF et al., Fecal calprotectin for screening of patients with suspected IBD: diagnostic metanalysis, BMJ 2010.
- D'Angelo F et al., Calprotectin in daily practice: where do we stand in 2017. Digestion, 2017.
- Sipponen T et al., Fecal Calprotectin in diagnosis and clinical assessment of IBD, Scand. J. Gastroenterol., 2015.
- Degrauwe P.L.J. et al., Faecal calprotectin in suspected paediatric IBD. Gastroenterology, 2015.

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