
Quantitative determination of alanine aminotransferase GPT (ALT) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Alanine aminotransferase (ALT) o Glutamate pyruvate transaminase (GPT) catalyses the reversible transfer of an amino group from alanine to α -ketoglutarate forming glutamate and pyruvate.

The pyruvate produced is reduced to lactate by lactate dehydrogenase (LDH) and NADH:

The rate of decrease in concentration of NADH, measured photometrically, is proportional to the catalytic concentration of ALT present in the sample¹.

CLINICAL SIGNIFICANCE

The ALT is a cellular enzyme, found in highest concentration in liver and kidney. High levels are observed in hepatic disease like hepatitis, diseases of muscles and traumatism, its better application is in the diagnosis of the diseases of the liver.

When they are used in conjunction with AST aid in the diagnosis of infarcts in the myocardium, since the value of the ALT stays within the normal limits in the presence of elevated levels of AST^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	TRIS pH 7,8	100 mmol/L
Buffer	L-Alanine	500 mmol/L
R 2	NADH	0,18 mmol/L
Substrate	Lactate dehydrogenase (LDH)	1200 U/L
	α -Ketoglutarate	15 mmol/L

PRECAUTIONS

R1: H290-May be corrosive to metals.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

DUAL MODE: Ready to use.

MONO MODE: Pour reagent 2 content over reagent 1. Mix thoroughly avoiding foam forming and it will be ready to use (WR).

Working reagent (WR) stability: 21 days at 2-8°C or 72 hours at room temperature (15-25°C).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1,00.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Thermostatic bath at 25°C, 30°C ó 37°C (\pm 0,1°C).
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment.

SAMPLES

Serum or plasma¹: Stability 7 days at 2-8°C.

INTERFERENCES

Anticoagulants currently in use like heparin, EDTA, oxalate and fluoride do not affect the results. Hemolysis interferes with the assay¹.

A list of drugs and other interfering substances with ALT determination has been reported^{2,3}.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINLAB 180

Name	GPT	Ref. male low	0 U/L
Abbr. Name	GPT	Ref. male high	40 U/L
Mode	Kinetic	Ref. female low	0 U/L
Wavelength	340 nm	Ref. female high	32 U/L
Units	U/L	Ref. Ped. Low	*
Decimals	0 U/L	Ref. Ped. High	*
Low Conc.	5 U/L	Panic value low	*
High Conc.	500 U/L	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1,000
		Correlat. offset	0,000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	240 μ L	normal volume	- μ L
rerun volume	240 μ L	rerun volume	- μ L
Sample		Sample	
normal volume	30,0 μ L	normal volume	- μ L
rerun volume	15,0 μ L	rerun volume	- μ L
R2 bottle (mL)	5 mL		
normal volume	60,0 μ L		
rerun volume	60,0 μ L		
Predilución	No		
Slope blank	No		
Delay, min. time	50, 186 sec.	Delay, min. time	- sec.
Linearity limit	10,0 %	Linearity limit	10,0 %
Factor	**	Factor	-
Reagent blank	No	Reagent blank	No
Low Absorbance	-0,100 Abs	Low Absorbance	-0,100 Abs
High Absorbance	3,000 Abs	High Absorbance	3,000 Abs
R. Abs. L. Limit	-0,100 Abs	R. Abs. L. Limit	-0,100 Abs
R. Abs. H. Limit	3,000 Abs	R. Abs. H. Limit	3,000 Abs
R. Abs. Deviation	3,000 Abs	R. Abs. Deviation	3,000 Abs

REFERENCE VALUES^{4,5}

	25°C	30°C	37°C
Men	up to 22 U/L	29 U/L	40 U/L
Women	up to 18 U/L	22 U/L	32 U/L

Normal newborns have been reported to show a reference range of up to double the adult, attributed to the neonate's hepatocytes. These values decline to adult levels by approximately 3 months of age.

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0 U/L to *linearity limit* of 400 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	Mean	SD
Mean (U/L)	42,0	116	41,1	115
SD	0,47	0,42	0,76	1,61
CV (%)	1,11	0,36	1,85	1,40

Sensitivity: 1 U/L = 0,00052 Δ A / min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r^2): 0,99597.

Regression equation: $y=1,1209x + 1,390$.

The results of the performance characteristics depend on the analyzer used.

NOTES

SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- Murray R. Alanine aminotransferase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1088-1090.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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APPLICATION AU SPINLAB 180

Name	GPT	Ref. male low	0 U/L
Abbr. Name	GPT	Ref. male high	40 U/L
Mode	Kinetic	Ref. female low	0 U/L
Wavelength	340 nm	Ref. female high	32 U/L
Units	U/L	Ref. Ped. Low	*
Decimals	0 U/L	Ref. Ped. High	*
Low Conc.	5 U/L	Panic value low	*
High Conc.	500 U/L	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1,000
		Correlat. offset	0,000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	240 µL	normal volume	- µL
rerun volume	240 µL	rerun volume	- µL
Sample		Sample	
normal volume	30,0 µL	normal volume	- µL
rerun volume	15,0 µL	rerun volume	- µL
R2 bottle (mL)	5 mL		
normal volume	60,0 µL		
rerun volume	60,0 µL		
Predilución	No		
Slope blank	No		
Delay, min. time	50, 186 sec.	Delay, min. time	- sec.
Linearity limit	10,0 %	Linearity limit	10,0 %
Factor	**	Factor	
Reagent blank	No	Reagent blank	No
Low Absorbance	-0,100 Abs	Low Absorbance	-0,100 Abs
High Absorbance	3,000 Abs	High Absorbance	3,000 Abs
R. Abs. L. Limit	-0,100 Abs	R. Abs. L. Limit	-0,100 Abs
R. Abs. H. Limit	3,000 Abs	R. Abs. H. Limit	3,000 Abs
R. Abs. Deviation	3,000 Abs	R. Abs. Deviation	3,000 Abs

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