
Quantitative determination of aspartate aminotransferase GOT (AST) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Aspartate aminotransferase (AST) formerly called glutamate oxaloacetate (GOT) catalyses the reversible transfer of an amino group from aspartate to α -ketoglutarate forming glutamate and oxalacetate. The oxalacetate produced is reduced to malate by malate dehydrogenase (MDH) and NADH:

The rate of decrease in concentration of NADH, measured photometrically, is proportional to the catalytic concentration of AST present in the sample¹.

CLINICAL SIGNIFICANCE

The AST is a cellular enzyme, is found in highest concentration in heart muscle, the cells of the liver, the cells of the skeletal muscle and in smaller amounts in other weaves.

Although an elevated level of AST in the serum is not specific of the hepatic disease, is used mainly to diagnostic and to verify the course of this disease with other enzymes like ALT and ALP.

Also, it is used to control the patients after myocardial infarction, in skeletal muscle disease and other^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1 Buffer	TRIS pH 7.8	80 mmol/L
	Lactate dehydrogenase (LDH)	800 U/L
	Malate dehydrogenase (MDH)	600 U/L
	L-Aspartate	200 mmol/L
R 2 Substrate	NADH	0.18 mmol/L
	α -Ketoglutarate	12 mmol/L

PRECAUTIONS

R1: H290-May be corrosive to metals.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

The reagent is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1.00.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Thermostatic bath at 25°C, 30°C o 37°C ($\pm 0.1^\circ\text{C}$)
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

Serum or plasma¹: Stability 7 days at 2-8°C.

INTERFERENCES

Anticoagulants currently in use like heparin, EDTA, oxalate and fluoride do not affect the results. Haemolysis interferes with the assay¹

A list of drugs and other interfering substances with AST determination has been reported^{2,3}.

APPLICATION SPINLAB 180

Name	GOT / AST	Ref. male low	0
Abbr. Name	GOT	Ref. male high	38
Mode	Kinetic	Ref. female low	0
Wavelength	340 nm	Ref. female high	31
Units	U/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	5 U/L	Control 1	*
High Conc.	500 U/L	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	-
R1 bottle (mL)	25 mL	R1 bottle (mL)	- mL
normal volume	240 μL	normal volume	- μL
rerun volume	240 μL	rerun volume	- μL
Sample		Sample	
normal volume	30,0 μL	normal volume	- μL
rerun volume	15,0 μL	rerun volume	- μL
R2 bottle (mL)	5 mL		
normal volume	60,0 μL		
rerun volume	60,0 μL		
Predilución	No		
Slope blank	No		
Delay, min. time	50, 186 sec.	Delay, min. time	- sec.
Linearity limit	10,0 %	Linearity limit	-
Factor	**	Factor	-
Reagent blank	No	Reagent blank	No
Low Absorbance	-0,100 Abs	Low Absorbance	-0,100 Abs
High Absorbance	3,000 Abs	High Absorbance	3,000 Abs
R. Abs. L. Limit	-0,100 Abs	R. Abs. L. Limit	-0,100 Abs
R. Abs. H. Limit	3,000 Abs	R. Abs. H. Limit	3,000 Abs
R. Abs. Deviation	3,000 Abs	R. Abs. Deviation	3,000 Abs

REFERENCE VALUES¹

	25°C	30°C	37°C
Men	up to 19 U/L	26 U/L	38 U/L
Women	up to 16 U/L	22 U/L	31 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0 U/L to *linearity limit* of 467 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	Mean	SD
Mean (U/L)	48,1	159	47,4	156
SD	0,56	0,57	1,42	4,35
CV (%)	1,16	0,36	3,00	2,79

Sensitivity: 1 U/L = 0,00053 $\Delta\text{A}/\text{min}$.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r^2): 0,99956.

Regression equation: $y = 1,042x - 0,342$

The results of the performance characteristics depend on the analyzer used.

NOTES

SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- Murray R. Aspartate aminotransferase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1112-1116.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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APPLICATION AU SPINLAB 180

Name	GOT / AST	Ref. male low	0
Abbr. Name	GOT	Ref. male high	38
Mode	Kinetic	Ref. female low	0
Wavelength	340 nm	Ref. female high	31
Units	U/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	5 U/L	Control 1	*
High Conc.	500 U/L	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1,000
		Correlat. offset	0,000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	-
R1 bottle (mL)	25 mL	R1 bottle (mL)	- mL
normal volume	240 µL	normal volume	- µL
rerun volume	240 µL	rerun volume	- µL
Sample		Sample	
normal volume	30,0 µL	normal volume	- µL
rerun volume	15,0 µL	rerun volume	- µL
R2 bottle (mL)	5 mL		
normal volume	60,0 µL		
rerun volume	60,0 µL		
Predilución	No		
Slope blank	No		
Delay, min. time	50, 186 sec.	Delay, min. time	- sec.
Linearity limit	10,0 %	Linearity limit	-
Factor	**	Factor	-
Reagent blank	No	Reagent blank	No
Low Absorbance	-0,100 Abs	Low Absorbance	-0,100 Abs
High Absorbance	3,000 Abs	High Absorbance	3,000 Abs
R. Abs. L. Limit	-0,100 Abs	R. Abs. L. Limit	-0,100 Abs
R. Abs. H. Limit	3,000 Abs	R. Abs. H. Limit	3,000 Abs
R. Abs. Deviation	3,000 Abs	R. Abs. Deviation	3,000 Abs

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