
## Quantitative determination of gamma-glutamyl transferase (γ-GT) IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Gamma-glutamyl transferase (γ-GT) catalyses the transfer of γ-glutamyl group from γ-glutamyl-p-nitroanilide to acceptor glycyglycine, according to the following reaction:

The rate of 2-nitro-5-aminobenzoic acid formation, measured photometrically, is proportional to the catalytic concentration of γ-GT present in the sample<sup>1,2</sup>.

### CLINICAL SIGNIFICANCE

Gamma-glutamyl transferase (γ-GT) is a cellular enzyme with wide tissue distribution in the body, primarily in the kidney, pancreas, liver and prostate.

Measurements of gamma-glutamyl transferase (γ-GT) activity are used in the diagnosis and treatment of hepatobiliary diseases such as biliary obstruction, cirrhosis or liver tumours<sup>1,2,5,6</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

<b>R 1</b>	TRIS	100 mmol/L
Buffer	Glycylglycine	100 mmol/L
<b>R 2</b>	L-γ-glutamyl-3-carboxy-4-nitroanilide	3 mmol/L
Substrate		

### PREPARATION

DUAL MODE: Ready to use.

MONO MODE: Pour reagent 2 content over reagent 1. Mix thoroughly avoiding foam forming and it will be ready to use (WR).

Working reagent (WR) stability: 21 days at 2-8°C or 5 days at room temperature (15-25°C).

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm  $\geq 1,20$ .

### ADDITIONAL EQUIPMENT

Spectrophotometer or colorimeter measuring at 405 nm.

- Thermostatic bath at 25°C, 30°C or 37°C ( $\pm 0,1^\circ\text{C}$ )
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment.

### SAMPLES

Serum<sup>1</sup>. γ-GT is stable for at least 3 days at 2-8°C, 8 hours at 15-25°C and 1 month at -20°C.

### INTERFERENCES

Plasma should not be used, anticoagulants inhibit the enzyme. Gross haemolysis interferes in the assay<sup>1</sup>.

A list of drugs and other interfering substances with γ-GT determination has been reported<sup>3,4</sup>.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## APPLICATION SPINLAB 180

Name	GGT	Ref. male low	11 U/L
Abbr. Name	GGT	Ref. male high	50 U/L
Mode	Kinetic	Ref. female low	7 U/L
Wavelength	405 nm	Ref. female high	32 U/L
Units	U/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	2 U/L	Panic value low	*
High Conc.	250 U/L	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1,000
		Correlat. offset	0,000
<b>DUAL MODE</b>		<b>MONO MODE</b>	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
Normal volume	240 μL	Normal volume	- μL
Rerun volume	240 μL	Rerun volume	- μL
Sample		Sample	
Normal volume	30,0 μL	Normal volume	- μL
Rerun volume	15,0 μL	Rerun volume	- μL
R2 bottle (mL)	5 mL		
Normal volume	60,0 μL		
Rerun volume	60,0 μL		
Predilución	No		
Slope blank	No		
Delay, min. time	50, 186 sec.	Delay, min. time	- sec.
Linearity limit	10,0 %	Linearity limit	10,0 %
Factor	**	Factor	-
Reagent blank	No	Reagent blank	No
Low Absorbance	-0,100 Abs	Low Absorbance	-0,100 Abs
High Absorbance	3,000 Abs	High Absorbance	3,000 Abs
R. Abs. L. Limit	-0,100 Abs	R. Abs. L. Limit	-0,100 Abs
R. Abs. H. Limit	3,000 Abs	R. Abs. H. Limit	3,000 Abs
R. Abs. Deviation	3,000 Abs	R. Abs. Deviation	3,000 Abs

### REFERENCE VALUES<sup>1</sup>

	25°C	30°C	37°C
Men	4-18 U/L	5-25 U/L	7-32 U/L
Women	6-28 U/L	8-38 U/L	11-50 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From detection limit of 2,0 U/L to linearity limit of 300U/L. If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

### Precision:

Mean (U/L)	Intra-assay (n=20)		Inter-assay (n=20)	
	SD	CV (%)	SD	CV (%)
38,3	0,39	1,03	40,1	2,05
190	0,53	0,28	198	1,16

**Sensitivity:** 1 U/L = 0,0008 ΔA/min.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,99990.

Regression equation: y = 1,334x - 1,493.

The results of the performance characteristics depend on the analyzer used.

### NOTES

**SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

### BIBLIOGRAPHY

- Gendler S. γ-GT. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1120-1123.
- Persijn J P et al. J Clin Chem Clin Biochem 1976; (14) 9: 421-427.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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**APPLICATION AU SPINLAB 180**

Name	GGT	Ref. male low	11 U/L
Abbr. Name	GGT	Ref. male high	50 U/L
Mode	Kinetic	Ref. female low	7 U/L
Wavelength	405 nm	Ref. female high	32 U/L
Units	U/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	2 U/L	Panic value low	*
High Conc.	250 U/L	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1,000
		Correlat. offset	0,000
<b>DUAL MODE</b>		<b>MONO MODE</b>	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
Normal volume	240 µL	Normal volume	- µL
Rerun volume	240 µL	Rerun volume	- µL
Sample		Sample	
Normal volume	30,0 µL	Normal volume	- µL
Rerun volume	15,0 µL	Rerun volume	- µL
R2 bottle (mL)	5 mL		
Normal volume	60,0 µL		
Rerun volume	60,0 µL		
Predilución	No		
Slope blank	No		
Delay, min. time	50, 186 sec.	Delay, min. time	- sec.
Linearity limit	10,0 %	Linearity limit	10,0 %
Factor	**	Factor	-
Reagent blank	No	Reagent blank	No
Low Absorbance	-0,100 Abs	Low Absorbance	-0,100 Abs
High Absorbance	3,000 Abs	High Absorbance	3,000 Abs
R. Abs. L. Limit	-0,100 Abs	R. Abs. L. Limit	-0,100 Abs
R. Abs. H. Limit	3,000 Abs	R. Abs. H. Limit	3,000 Abs
R. Abs. Deviation	3,000 Abs	R. Abs. Deviation	3,000 Abs

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	25°C	30°C	37°C
Femmes	4-18 U/L	5-25 U/L	7-32 U/L
Hommes	6-28 U/L	8-38 U/L	11-50 U/L

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