# CK-MB-LQ (Creatine kinase-MB)

Anti CK-M. Immunoinhibition. Kinetic UV. Liquid

## Quantitative determination of creatine kinase MB (CK-MB) IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

The procedure involves measurement of CK activity in the presence of an antibody to CK-M monomer. This antibody completely inhibits the activity of CK-MM and half of the activity of CK-MB while not affecting the B subunit activity of CK-MB and CK-BB. Then it's used the CK method to quantitatively determine CK-B activity<sup>1,2</sup>. The CK-MB activity is obtained by multiplying the CK-B activity by two.

### CLINICAL SIGNIFICANCE

CK-MB is an enzyme formed by the association of two subunits from muscle (M) and nerve cells (B). CK-MB is usually present in serum at low concentration; it is increased after an acute infarct of myocardium and later descends at normal levels. Also, is increased, rarely, in skeletal muscle damage<sup>5,6</sup>. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

|  |                                   |             |
|--|-----------------------------------|-------------|
| R 1  | Imidazol, pH 6.7                  | 125 mmol/L  |
|  | D-Glucose                         | 25 mmol/L   |
|  | N-Acetyl-L-Cysteine               | 25 mmol/L   |
|  | Magnesium acetate                 | 12,5 mmol/L |
|  | NADP                              | 2,52 mmol/L |
|  | EDTA                              | 2,02 mmol/L |
|  | Hexokinase                        | ≥6 800 U/L  |
| Anti-human polyclonal CK-M antibody (sheep) sufficient to inhibit up to 2 000 U/L of CK-MM |                                   |             |
| R 2  | ADP                               | 15,2 mmol/L |
|  | AMP                               | 25 mmol/L   |
|  | di-Adenosine-5- pentaphosphate    | 103 mmol/L  |
|  | Glucose-6-phosphate dehydrogenase | ≥8 800 U/L  |
|  | Creatine phosphate                | 250 mmol/L  |
|  |                                   |             |

### Optional

|                        |                         |              |
|------------------------|-------------------------|--------------|
| CK-Nac / CK-MB CONTROL | Lyophilized human serum | Ref: 1002260 |
|------------------------|-------------------------|--------------|

### PREPARATION

MONO MODE: Pour reagent 2 content over reagent 1. Mix thoroughly avoiding foam forming and it will be ready to use (WR).

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented.

Do not use reagents over the expiration date.

#### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm ≥ 1.60.

### ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Thermostatic bath at 25°C, 30°C ó 37° C (± 0.1°C).
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

### LIMITATION OF THE PROCEDURE

- 1- The method will also measure any CK-BB isoenzyme present in serum. The activity of the isoenzyme is usually negligible, however, if a significant amount of CK-BB activity is present the CK-MB activity will be overestimated.
- 2- A macro form of BB (immunoglobulin complexed) has been observed wich will be measured as B in the assay. If the measured CK-B activity exceeds 20% of the total CK activity, the presence of macro BB should be suspected.

### SAMPLES

Serum free of hemolysis or heparin plasma<sup>1</sup>: Stability 7 days at 2-8°C, protected from light.

CK-MB activity decreases a 10% after 24 hours at 4°C or 1 hour at 25°C.

### REFERENCE VALUES

The suspicion of myocardial damage is based on the three following factors:

|          |              |          |          |         |
|----------|--------------|----------|----------|---------|
| CK-MB    | 25°C         | 30°C     | 37°C     |         |
|          | > 10 U/L     | > 15 U/L | > 24 U/L |         |
|          | TOTAL CK     |          |          |         |
| TOTAL CK | 25°C         | 30°C     | 37°C     |         |
|          | Men, up to   | 80 U/L   | 130 U/L  | 195 U/L |
|          | Women, up to | 70 U/L   | 110 U/L  | 170 U/L |

$\frac{\text{CK - MB Activity}}{\text{CK Total Activity}} \times 100$ : 6 - 25 % CK - MB Activity in the sample

These values are for orientation purpose; each laboratory should establish its own reference range.

### APPLICATION TO SPINLAB 180

|                   |               |                   |               |
|-------------------|---------------|-------------------|---------------|
| Name              | CK-MB.        | Ref. male low     | 0.0           |
| Abbr. Name        | CKMB          | Ref. male high    | 24.0          |
| Mode              | Kinetic       | Ref. female low   | 0.0           |
| Wavelength        | 340 nm        | Ref. female high  | 24.0          |
| Units             | U/L           | Ref. Ped. Low     | *             |
| Decimals          | 1             | Ref. Ped. High    | *             |
| Low Conc.         | 2.0 U/L       | Control 1         | *             |
| High Conc.        | 1000.0 U/L    | Control 2         | *             |
| Calibrator name   | CAL           | Control 3         | *             |
| Prozone check     | No            | Correlat. factor  | 1.000         |
|                   |               | Correlat. offset  | 0.000         |
| <b>DUAL MODE</b>  |               | <b>MONO MODE</b>  |               |
| Sample blank      | No            | Sample blank      | No            |
| R1 bottle (mL)    | 25 mL         | R1 bottle (mL)    | 25 mL         |
| normal volume     | 300 µL        | normal volume     | 300 µL        |
| rerun volume      | 300 µL        | rerun volume      | 300 µL        |
| Sample            |               | Sample            |               |
| normal volume     | 12.0 µL       | normal volume     | 12.0 µL       |
| rerun volume      | 6.0 µL        | rerun volume      | 6.0 µL        |
| R2 bottle (mL)    | 5 mL          |                   |               |
| normal volume     | 0 µL          |                   |               |
| rerun volume      | 0 µL          |                   |               |
| Predilución       | No            |                   |               |
| Slope blank       | No            |                   |               |
| Delay, min. time  | 360, 159 sec. | Delay, min. time  | 304, 118 sec. |
| Linearity limit   | 10.0 %        | Linearity limit   | 10.0 %        |
| Factor            |               | Factor            |               |
| Reagent blank     | No            | Reagent blank     | No            |
| Low Absorbance    | -0.100 Abs    | Low Absorbance    | -0.100 Abs    |
| High Absorbance   | 3.000 Abs     | High Absorbance   | 3.000 Abs     |
| R. Abs. L. Limit  | -0.100 Abs    | R. Abs. L. Limit  | -0.100 Abs    |
| R. Abs. H. Limit  | 3.000 Abs     | R. Abs. H. Limit  | 3.000 Abs     |
| R. Abs. Deviation | 3.000 Abs     | R. Abs. Deviation | 3.000 Abs     |

### QUALITY CONTROL

CK-Nac/CK-MB specific control sera (Ref. 1002260) are recommended to monitor the performance of assay procedures.

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From *detection limit* of 1,9 U/L to *linearity limit* of 318 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/1 with NaCl 9 g/L and multiply the result by 2.

#### Precision:

| Mean (U/L) | Intra-assay |       | Inter-assay |       |
|------------|-------------|-------|-------------|-------|
|            | 33,7        | 166,5 | 31,3        | 161,0 |
| SD         | 1,00        | 3,76  | 1,19        | 3,47  |
| CV (%)     | 2,96        | 2,26  | 3,81        | 2,15  |

**Sensitivity:** 1 U/L = 0,000134 (A).

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained were the following:

Correlation coefficient (r)<sup>2</sup>: 0,999.

Regression equation: y = 0,976 x - 0,269.

The results of the performance characteristics depend on the analyzer used.

### BIBLIOGRAPHY

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommendation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.
8. Neumeier, D., Prellwitz, W., Würzburg, U. et coll. Determination of creatine kinase isoenzyme MB activity in serum using immunological inhibition of creatine kinase M subunit activity. Activity kinetics and diagnostic significance in myocardial infarction, Clin. Chim. Acta, 73, (1976), 445.

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- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
- Gerhardt W. et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
- Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la créatinine kinase dans le sérum humain. Ann. Biol. Clin.,40, (1482), 87.
- Neumeier, D., Prellwitz, W., Würzburg, U. et coll. Determination of creatine kinase isoenzyme MB activity in serum using immunological inhibition of creatine kinase M subunit activity. Activity kinetics and diagnostic significance in myocardial infarction. Clin. Chim. Acta, 73, (1976), 445.

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### BIBLIOGRAPHIE

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la créatine kinase dans le sérum humain. Ann. Biol. Clin.,40, (1482), 87.
8. Neumeier, D., Prellwitz, W., Würzburg, U. et coll. Determination of creatine kinase isoenzyme MB activity in serum using immunological inhibition of creatine kinase M subunit activity. Activity kinetics and diagnostic significance in myocardial infarction, Clin. Chim. Acta, 73, (1976), 445.

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|              |       |                |
|--------------|-------|----------------|
| Ref: SP41254 | Cont. | R1: 10 x 20 mL |
|              |       | R2: 10 x 5 mL  |

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### BIBLIOGRAFIA

- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
- Gerhardt W. et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
- Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la créatinine kinase dans le sérum humain. Ann. Biol. Clin.,40, (1482), 87.
- Neumeier, D., Prellwitz, W., Würzburg, U. et coll. Determination of creatine kinase isoenzyme MB activity in serum using immunological inhibition of creatine kinase M subunit activity. Activity kinetics and diagnostic significance in myocardial infarction, Clin. Chim. Acta, 73, (1976), 445.

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Ref: SP41254

Cont.

R1: 10 x 20 mL

R2: 10 x 5 mL