
Quantitative determination of lipase
IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The pancreatic lipase in presence of colipase, desoxycholate and calcium ions, hydrolyses the substrate 1-2-O-dilauryl-rac-glycero-3-glutaric acid-(6'-methylresorufin)-ester. The sequence of reactions involved in the enzymatic direct lipase determination is the following:

The rate of methylresorufin formation, measured photometrically, is proportional to the catalytic concentration of lipase present in the sample.

CLINICAL SIGNIFICANCE

Lipase (LPS) is a pancreatic enzyme necessary for the absorption and digestion of nutrients that catalyzes the hydrolysis of glycerol esters of fatty acids. Determination of LPS is used for diagnosis of diseases of pancreas such as acute and chronic pancreatitis and obstruction of the pancreatic duct^{1,7,8}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1 Buffer	TRIS pH 8.3	40 mmol/L
	Colipase	> 1 mg/L
	Desoxycholate	1,8 mmol/L
R 2 Substrate (micro-emulsion)	Tartrate pH 4,0	15 mmol/L
	Lipase	> 0,7 mmol/L
	Calcium chloride (CaCl ₂)	0,1 mmol/L
OPTIONAL	SPINTROL H CAL	

PREPARATION

- R1 – R2 Ready to use. Stability after opening 90 days at 2-8°C. R2 Mix gently before use (Nota 1).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 580 nm \geq 1.00.
- R 2 is a turbid orange-colored micro-emulsion, discard if turning to red.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 580 nm.
- Thermostatic bath at 37° C (\pm 0.1°C)
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment (Nota 2).

SAMPLES

Serum or plasma with sodium citrate, EDTA or heparin¹.
Avoid repeated frozen and unfrozen.
Stability: 2 days at 2-8°C.

REFERENCE VALUES¹

\leq 38 U/L (U/L methylresorufin at 37°C).

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINLAB 180

Name	LIPASE	Ref. male low	0
Abbr. Name	LIPA	Ref. male high	30
Mode	Kinetic	Ref. female low	0
Wavelength	578 nm	Ref. female high	30
Units	U/L	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	1.0 U/L	Panic value low	*
High Conc.	250.0 U/L	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
Normal volume	300 μ L		
Rerun volume	300 μ L		
Sample			
Normal volume	3.0 μ L		
Rerun volume	2.0 μ L		
R2 bottle (mL)	5 mL		
Normal volume	60.0 μ L		
Rerun volume	60.0 μ L		
Predilución	No		
Slope blank	No		
Delay, min. time	50, 186 sec.		
Linearity limit	10.0 %		
Factor			
Reagent blank	Yes (0.000)		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
R. Abs. Deviation	3.000 Abs		

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 5 U/L to *linearity limit* of 250 U/L.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	CV (%)	
Mean (U/L)	40,2	59,35	38,5	58,9
SD	0,410	0,875	1,10	1,25
CV (%)	1,02	1,47	2,86	2,13

Sensitivity: 1 U/L = 0,00059792 (A)

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 101 samples were the following:

Correlation coefficient (r)²: 0,99732.

Regression equation: y = 0,50054x + 3,9443.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. In some storage conditions (i.e. storage at a temperature lower than the one indicate) a precipitate may appear in the vial that will not influence that the reagent performance; however, it is recommended to resuspend the product with a slight rotation.
2. In order to avoid contamination, it is recommended to use disposable material.
3. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

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APPLICATION AU SPINLAB 180

Name	LIPASE	Ref. male low	0
Abbr. Name	LIPA	Ref. male high	30
Mode	Kinetic	Ref. female low	0
Wavelength	578 nm	Ref. female high	30
Units	U/L	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	1.0 U/L	Panic value low	*
High Conc.	250.0 U/L	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
Normal volume	300 µL		
Rerun volume	300 µL		
Sample			
Normal volume	3.0 µL		
Rerun volume	2.0 µL		
R2 bottle (mL)	5 mL		
Normal volume	60.0 µL		
Rerun volume	60.0 µL		
Predilución	No		
Slope blank	No		
Delay, min. time	50, 186 sec.		
Linearity limit	10.0 %		
Factor			
Reagent blank	Yes (0.000)		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
R. Abs. Deviation	3.000 Abs		

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