
**Quantitative determination of urea IVD**

Store at 2-8°C

**PRINCIPLE OF THE METHOD**

 Urea in the sample is hydrolyzed enzymatically into ammonia (NH<sub>4</sub><sup>+</sup>) and carbon dioxide (CO<sub>2</sub>).

 Ammonia ions formed reacts with α-ketoglutarate in a reaction catalysed by glutamate dehydrogenase (GLDH) with simultaneous oxidation of NADH to NAD<sup>+</sup>:

 The decrease in concentration of NADH, is proportional to urea concentration in the sample<sup>1</sup>.

**CLINICAL SIGNIFICANCE**

Urea is the final result of the metabolism of proteins; It is formed in the liver from their destruction.

 It can appear the urea elevated in blood (uremia) in: diets with excess of proteins, renal diseases, heart failure, gastrointestinal hemorrhage, dehydration or renal obstruction<sup>1,4,5</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**REAGENTS**

<b>R 1</b>	TRIS pH 7.8	80 mmol/L
<b>Buffer</b>	α-Ketoglutarate	6 mmol/L
	Urease	75000 U/L
<b>R 2</b>	GLDH	60000 U/L
<b>Enzymes</b>	NADH	0.32 mmol/L
<b>Optional</b>	SPINTROL H CAL	

**PREPARATION**

DUAL MODE: Ready to use.

MONO MODE: Pour reagent 2 content over reagent 1. Mix thoroughly avoiding foam forming and it will be ready to use (WR). The Working reagent (WR) is stable for 1 month at 2-8°C.

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

**Signs of reagent deterioration:**

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm < 1.00.

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or colorimeter measuring at 340 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment<sup>(Note 1)</sup>.

**SAMPLES**

 - Serum or heparinized plasma<sup>1</sup>: Do not use ammonium salts or fluoride as anticoagulants.

 - Urine<sup>1</sup>: Dilute sample 1/50 in distilled water. Mix. Multiply the results by 50 (dilution factor). Preserve urine samples at pH < 4.

Urea is stable at 2-8°C for 5 days.

**INTERFERENCES**

 It is recommended to use heparin as anticoagulant. Do not use ammonium salts or fluoride<sup>1</sup>.

 A list of drugs and other interfering substances with urea determination has been reported by Young et. al<sup>2,3</sup>.

**QUALITY CONTROL**

Control Sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagent and calibration for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**APPLICATION SPINLAB 180**

Name	UREA	Ref. male low	15 mg/dL
Abbr. Name	UREA	Ref. male high	45 mg/dL
Mode	Twopoint	Ref. female low	15 mg/dL
Wavelength	340 nm	Ref. female high	45 mg/dL
Units	mg/dL	Ref. Ped. Low	15 mg/dL
Decimals	0	Ref. Ped. High	45 mg/dL
Low Conc.	15	Panic value low	*
High Conc.	250	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>		<b>MONO MODE</b>	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	30 mL
Normal volume	240 µL	Normal volume	300 µL
Rerun volume	240 µL	Rerun volume	300 µL
Sample		Sample	
Normal volume	3.0 µL	Normal volume	3.0 µL
Rerun volume	2.0 µL	Rerun volume	2.0 µL
R2 bottle (mL)	5 mL		
Normal volume	60.0 µL		
Rerun volume	60.0 µL		
Predilución	No		
Slope blank	No		
1 <sup>st</sup> , 2 <sup>nd</sup> time	24, 103 sec.	1 <sup>st</sup> , 2 <sup>nd</sup> time	24, 103 sec.
Factor		Factor	
Reagent blank	Yes	Reagent blank	Yes
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

**REFERENCE VALUES<sup>4,5</sup>**

Serum or plasma:

15-45 mg/dL ≅ 2.5-7.5 mmol/L

Urine:

26 – 43 g/24 h ≅ 428-714 mmol/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

**PERFORMANCE CHARACTERISTICS**
**Measuring range:** From *detection limit* 0,743 mg/dL to *linearity limit* 400 mg/dL.

If the concentration is greater than linearity limit dilute 1/2 the sample with NaCl 9 g/L and multiply the result by 2.

**Precision:**

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mg/dL)	SD	Mean (mg/dL)	SD
Mean (mg/dL)	37,5	120	40,0	126
SD	1,05	0,92	1,06	2,07
CV (%)	2,79	0,77	2,65	1,65

**Sensitivity:** 1 mg/dL = 0,00180 A.

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagent (x).

The results obtained using 50 samples was the following:

 Correlation coefficient (r)<sup>2</sup>: 0,98209.

Regression equation y= 1,0343x – 1,2105.

The results of the performance characteristics depend on the analyzer used.

**NOTES**

1. Glassware and distilled water must be free of ammonia and ammonium salts<sup>1</sup>.
2. Use clean disposable pipette tips for its dispensation.
3. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

**BIBLIOGRAPHY**

1. Kaplan A. Urea. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1257-1260 and 437 and 418.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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