
Quantitative determination of uric acid IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Uric acid is oxidized by uricase to allantoin and hydrogen peroxide (2H₂O₂), which under the influence of POD, 4-aminophenazone (4-AP) and 2-4 Dichlorophenol sulfonate (DCPS) forms a red quinoneimine compound:

The intensity of the red color formed is proportional to the uric acid concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Uric acid and its salts are end products of the purine metabolism. With progressive renal insufficiency, there is retention in blood of urea, creatinine and uric acid. Elevated uric acid level may be indicative of renal insufficiency and is commonly associated with gout^{1,5,6}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1 Buffer	Phosphate pH 7.4	50 mmol/L
	2-4 Dichlorophenol sulfonate (DCPS)	4 mmol/L
R 2 Enzymes	Uricase	60 U/L
	Peroxidase (POD)	660 U/L
	Ascorbate oxidase	200 U/L
	4 - Aminophenazone (4-AP)	1 mmol/L
Optional	SPINTROL H CAL	

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 520 nm ≥ 0.16.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 520 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

- Serum or plasma¹: Stability 3-5 days at 2-8°C or 6 months at -20°C.
- Urine (24 h)¹: Stability 4 days at 15-25°C, pH >8. Dilute sample 1/50 in distilled water. Mix. Multiply results by 50 (dilution factor); If urine is cloudy; warm the specimen to 60°C for 10 min to dissolve precipitated urates and uric acid. Do not refrigerate.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINLAB 180

Name	Uric Acid	Ref. male low	3.40
Abbr. Name	URI	Ref. male high	7.00
Mode	Two point	Ref. female low	2.50
Wavelength	505 nm	Ref. female high	6.00
Units	mg/dL	Ref. Ped. Low	2.00
Decimals	2	Ref. Ped. High	5.50
Low Conc.	0.03 mg/dL	Panic value low	*
High Conc.	25.00 mg/dL	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	160 µL		
rerun volume	160 µL		
Sample			
normal volume	8.0 µL		
rerun volume	4.0 µL		
R2 bottle (mL)	5 mL		
normal volume	160 µL		
rerun volume	160 µL		
Predilución	No		
Slope blank	No		
Point one, two	-3, 236 sec.		
Factor			
Reagent blank	Yes (0.000)		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		

REFERENCE VALUES¹

Serum or plasma:

Women	2.5 - 6.8 mg/dL	≅	149 - 405 µmol/L
Men	3.6 - 7.7 mg/dL	≅	214 - 458 µmol/L

Urine: 250 - 750 mg/24 h ≅ 1.49 - 4.5 mmol/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0,01647 mg/dL to *linearity limit* of 40 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

Mean (mg/L)	Intra-assay (n=20)		Inter-assay (n=20)	
	4,46	10,37	4,71	11,02
SD	0,02	0,05	0,06	0,15
CV (%)	0,46	0,44	1,20	1,37

Sensitivity: 1 mg/dL = 0,0323 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99734.

Regression equation: y=0,816x + 0,319.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. Use clean disposable pipette tips for its dispensation.
2. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

1. Schultz A. Uric acid. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1261-1266 and 418.
2. Fossati P et al. Clin Chem 1980;26:227-231.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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APPLICATION AU SPINLAB 180

Name	Uric Acid	Ref. male low	3.40
Abbr. Name	URI	Ref. male high	7.00
Mode	Two point	Ref. female low	2.50
Wavelength	505 nm	Ref. female high	6.00
Units	mg/dL	Ref. Ped. Low	2.00
Decimals	2	Ref. Ped. High	5.50
Low Conc.	0.03 mg/dL	Panic value low	*
High Conc.	25.00 mg/dL	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	160 μL		
rerun volume	160 μL		
Sample			
normal volume	8.0 μL		
rerun volume	4.0 μL		
R2 bottle (mL)	5 mL		
normal volume	160 μL		
rerun volume	160 μL		
Predilución	No		
Slope blank	No		
Point one, two	-3, 236 sec.		
Factor			
Reagent blank	Yes (0.000)		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		

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