
Quantitative determination of glucose
IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Glucose oxidase (GOD) catalyses the oxidation of glucose to gluconic acid. The formed hydrogen peroxide (H₂O₂), is detected by a chromogenic oxygen acceptor, phenol, 4 – aminophenazone (4-AP) in the presence of peroxidase (POD):

The intensity of the color formed is proportional to the glucose concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Glucose is a major source of energy for most cells of the body; insulin facilitates glucose entry into the cells. Diabetes is a disease manifested by hyperglycemia; patients with diabetes demonstrate an inability to produce insulin^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	TRIS pH 7.4	92 mmol/L
	Phenol	0,3 mmol/L
	Glucose oxidase (GOD)	15000 U/L
	Peroxidase (POD)	1000 U/L
	4 – Aminophenazone (4-AP)	2,6 mmol/L
Optional	SPINTROL H CAL	

PREPARATION

Reagent and calibrator provided are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 505 nm \geq 0.32.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 505 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

Serum or plasma, free of hemolysis¹:

Serum should be removed from the clot as quickly as possible.

Stability of the sample: Glucose in serum or plasma is stable at 2-8° for 3 days.

REFERENCE VALUES¹

Serum or plasma:

60 – 110 mg/dL \cong 3.33 – 6.10 mmol/L

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagent and calibration for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION TO SPINLAB-180

Name	Glucose	Ref. male low	60 mg/dL
Abbr. Name	GLU	Ref. male high	110 mg/dL
Mode	Endpoint	Ref. female low	60 mg/dL
Wavelength	505 nm	Ref. female high	110 mg/dL
Units	mg/dL	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	20 mg/dL	Control 1	*
High Conc.	500 mg/dL	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 μ L	normal volume	300 μ L
rerun volume	300 μ L	rerun volume	300 μ L
Sample		Sample	
normal volume	3.0 μ L	normal volume	3.0 μ L
rerun volume	2.0 μ L	rerun volume	2.0 μ L
R2 bottle (mL)	5 mL		
normal volume	0 μ L		
rerun volume	0 μ L		
Predilución	No		
Incubation time	4.5 min.	Incubation time	4.5 min.
Factor		Factor	
Reagent blank	Yes (0.000)	Reagent blank	Yes (0.000)
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* 0,3709 mg/dL to *linearity limit* 500 mg/dL.

If the concentration is greater than linearity limit dilute 1/2 the sample with NaCl 9 g/L and multiply the result by 2.

Precision:

Mean (mg/dL)	Intra-assay (n=20)		Inter-assay (n=20)	
		98,5	265	92,5
SD	0,58	1,27	2,76	6,44
CV (%)	0,59	0,48	2,98	2,57

Sensitivity: 1 mg/dL = 0,0039 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagent (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99492.

Regression equation: y=1,104x - 1,249.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. Use clean disposable pipette tips for its dispensation.
2. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

1. Kaplan L.A. Glucose. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1032-1036.
2. Trinder P. Ann Clin Biochem 1969; 6 24-33.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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APPLICATION AU SPINLAB 180

Name	Glucose	Ref. male low	60 mg/dL
Abbr. Name	GLU	Ref. male high	110 mg/dL
Mode	Endpoint	Ref. female low	60 mg/dL
Wavelength	505 nm	Ref. female high	110 mg/dL
Units	mg/dL	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	20 mg/dL	Control 1	*
High Conc.	500 mg/dL	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 µL	normal volume	300 µL
rerun volume	300 µL	rerun volume	300 µL
Sample		Sample	
normal volume	3.0 µL	normal volume	3.0 µL
rerun volume	2.0 µL	rerun volume	2.0 µL
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volume	0 µL		
Predilución	No		
Incubation time	4.5 min.	Incubation time	4.5 min.
Factor			
Reagent blank	Yes (0.000)	Reagent blank	Yes (0.000)
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs

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