Quantitative determination of total protein IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Proteins give an intensive violet-blue complex with copper salts in an alkaline medium. Iodide is included as an antioxidant.

 The intensity of the color formed is proportional to the total protein concentration in the sample^{1,4}.

CLINICAL SIGNIFICANCE

The proteins are macromolecular organic compounds, widely distributed in the organism. They act like structural and transport elements. The proteins of the serum are divided in two fractions, albumin and globulins

The determination of total proteins is useful in the detection of:

- High protein levels caused by hemoconcentration like in the dehydrations or increase in the concentration of specific proteins.
- Low protein level caused by hemodilution by an impaired synthesis or loss (as by hemorrhage) or excessive protein catabolism^{4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R	Sodium potassium tartrate	15 mmol/L	
	Sodium iodide	100 mmol/L	
	Biuret	Potassium iodide	5 mmol/L
	Copper (II) sulphate	5 mmol/L	
Optional	SPINTROL H CAL		

PRECAUTION

R: H314-Causes severe skin burns and eye damage. H412-Harmful to aquatic life with long lasting effects. Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

The reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 540 nm \geq 0.22.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 540 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

SAMPLES

 Serum or heparinized plasma¹:

Stability of the sample: 1 month at refrigerator (2-8°C).

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINLAB 180

Name	TOTAL PROT.	Ref. male low	6.7 g/dL
Abbr. Name	TP	Ref. male high	8.7 g/dL
Mode	Endpoint	Ref. female low	6.7 g/dL
Wavelength	546 nm	Ref. female high	8.7 g/dL
Units	g/dL	Ref. Ped. Low	5.2 g/dL
Decimals	2	Ref. Ped. High	9.1 g/dL
Low Conc.	1.0 g/dL	Panic value low	*
High Conc.	15.0 g/dL	Panic value high	*
Calibrator	CAL	Control 1	*
Cutt-off	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
I			
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	280 μ L	normal volume	280 μ L
rerun volume	280 μ L	rerun volume	280 μ L
Sample		Sample	
normal volume	7.0 μ L	normal volume	7.0 μ L
rerun volume	5.0 μ L	Vol. normal	5.0 μ L
R2 bottle (mL)	5 mL		
normal volume	0 μ L		
rerun volume	0 μ L		
Predilución	No		
Incubation time	11.5 min.	Incubation time	11.5 min.
Factor	*	Factor	*
Reagent blank	Yes	Reagent blank	Yes
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs

REFERENCE VALUES¹

Adults: 6.6 – 8.3 g/dL

Newborn: 5.2 – 9.1 g/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS
Measuring range: From *detection limit* of 0,007 g/dL to *linearity limit* of 14 g/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (g/dL)	6,53	4,89	6,77	5,08
SD	0,01	0,01	0,07	0,05
CV (%)	0,21	0,24	1,05	0,94

Sensitivity: 1 g/dL = 0,0825 A.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

 Correlation coefficient (r)²: 0,97002

Regression equation y= 0,954x +0,511.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. Use clean disposable pipette tips for its dispensation.
2. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

1. Koller A. Total serum protein. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1316-1324 and 418.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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APPLICATION AU SPINLAB 180

Name	TOTAL PROT.	Ref. male low	6.7 g/dL
Abbr. Name	TP	Ref. male high	8.7 g/dL
Mode	Endpoint	Ref. female low	6.7 g/dL
Wavelength	546 nm	Ref. female high	8.7 g/dL
Units	g/dL	Ref. Ped. Low	5.2 g/dL
Decimals	2	Ref. Ped. High	9.1 g/dL
Low Conc.	1.0 g/dL	Panic value low	*
High Conc.	15.0 g/dL	Panic value high	*
Calibrator	CAL	Control 1	*
Cutt-off	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
I			
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	280 µL	normal volume	280 µL
rerun volume	280 µL	rerun volume	280 µL
Sample		Sample	
normal volume	7.0 µL	normal volume	7.0 µL
rerun volume	5.0 µL	Vol. normal	5.0 µL
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volume	0 µL		
Predilución	No		
Incubation time	11.5 min.	Incubation time	11.5 min.
Factor	*	Factor	*
Reagent blank	Yes	Reagent blank	Yes
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs

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