
Quantitative determination of chloride ion IVD.

Store at 2-8°C

PRINCIPLE OF THE METHOD

The quantitative displacement of thiocyanate by chloride from mercuric thiocyanate and subsequent formation of a red ferric thiocyanate complex is measured colorimetrically:

The intensity of the color formed is proportional to the chloride ion concentration in the sample^{1,2,3,4}.

CLINICAL SIGNIFICANCE

It is important clinically the determination of chloride due regulation of osmotic pressure of extra cellular fluid and to its significant role in acid-base balance. Increases in chloride ion concentration may be found in severe dehydration, excessive intake of chloride, severe renal tubular damage and in patients with cystic fibrosis.

Decrease in chloride ion concentration may be found in metabolic acidosis, loss from prolonged vomiting and chronic pyelonephritis^{2,7,8}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R Thiocyanate-Hg	Mercuric thiocyanate	2 mmol/L
	Ferric nitrate	40 mmol/L
	Mercuric nitrate	0,15 mmol/L
	Nitric acid	45 mmol/L
Opcional	SPINTROL H CAL	

PRECAUTIONS

R: H314-Causes severe skin burns and eye damage.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 480 nm $\geq 0,15$.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 480 nm.
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment (Note 1, 2, 3).

SAMPLES

- Serum or plasma free of haemolysis and separated from cells as rapidly as possible. Anticoagulants such as oxalate or EDTA are not acceptable they will interfere with results.

- Urine¹: Collect 24-hour urine specimen in chloride free containers. Dilute a sample 1/2 in distilled water. Mix. Multiply results by 2 (dilution factor).

Stability of the sample: Ion chloride is stable 1 week at room temperature (15-25°C), in refrigerator (2-8°C) or frozen (-20°C) temperatures.

REFERENCE VALUES¹

Serum or plasma:	95 - 115 mmol/L
Urine:	110 - 250 mmol/24h

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINLAB-180

Name	Chloride	Ref. male low	98
Abbr. Name	CL	Ref. male high	107
Mode	Endpoint	Ref. female low	98
Wavelength	505 nm	Ref. female high	107
Units	mmol/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	0 mmol/L	Control 1	*
High Conc.	130 mmol/L	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 µL	normal volume	300 µL
rerun volume	300 µL	rerun volume	300 µL
Sample		Sample	
normal volume	3,0 µL	normal volume	3,0
rerun volume	2,0 µL	rerun volume	2,0
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volume	0		
Predilución	No		
Incubation time	4,5 min.	Incubation time	11,5 min.
Factor		Factor	
Vol. repet.		Vol. repet.	30,0 µL
Reagent blank	Yes (0,000)	Reagent blank	Yes (0,000)
Low Absorbance	-0,100 Abs	Low Absorbance	-0,100 Abs
High Absorbance	3,000 Abs	High Absorbance	3,000 Abs
R. Abs. L. Limit	-0,100 Abs	R. Abs. L. Limit	-0,100 Abs
R. Abs. H. Limit	3,000 Abs	R. Abs. H. Limit	3,000 Abs

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0,454 mmol/L to *linearity limit* of 190 mmol/L.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with distilled water and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (mmol/L)	SD	Mean (mmol/L)	SD
Mean (mmol/L)	84,2	114	82,5	111
SD	0,81	0,62	1,07	1,87
CV (%)	0,96	0,55	1,30	1,68

Sensitivity: 1 mmol/L = 0,00471 A.

Accuracy: Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagents. The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,96731.

Regression equation: y=0,990x + 0,100.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. It is recommended to use disposable material. If glassware is used the material should be scrupulously cleaned with H₂SO₄ - K₂Cr₂O₇ Solution and then thoroughly rinsed it with distilled water.
2. Most of the detergents and water softening products used in the laboratories contains chelating agents. A defective rinsing will invalidate the procedure.
3. Avoid the contact with metal materials.
4. Use clean disposable pipette tips for its dispensation.
5. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

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