
Quantitative determination of iron IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The iron is dissociated from transferring-iron complex in weakly acid medium. Liberated iron is reduced into the bivalent form by means of ascorbic acid. Ferrous ions give with FerroZine a coloured complex:

 The intensity of the color formed is proportional to the iron concentration in the sample^{1,2}.

CLINICAL SIGNIFICANCE

The iron is the component of a great number of enzymes. The myoglobin, muscular protein, contains iron, as well as the liver.

Iron is necessary for the hemoglobin production, molecule that transports oxygen inside red globules. Their deficit in the last causes the ferropenic anemia. High levels of iron are found in hemochromatosis, cirrhosis, hepatitis and in increased transferrin levels.

 The variation day to day is quite marked in healthy people^{1,5,6}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Buffer	Acetate pH 4.9	100 mmol/L
R 2	Reductant	Ascorbic acid	99.7%
R 3	Color	FerroZine	40 mmol/L
Optional		SPINTROL H CAL	

PREPARATION

Dissolve (→) the contents of one tube of reductant (R2) in one bottle of R 1 Buffer. Cap and mix gently to dissolve contents. Stability: 3 months at 2-8°C or 1 month at 15-25°C.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 562 nm \geq 0.020.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 562 nm.
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment (Note 1).

SAMPLES

 Serum or heparinized plasma. Free of hemolysis and separated from cells as rapidly as possible. Stability of the sample: 2-8°C for 7 days¹.

REFERENCE VALUES⁵

Male	65 - 175 µg/dL	\cong	11.6 - 31.3 µmol/L (Note 4)
Female	40 - 150 µg/dL	\cong	7.16 - 26.85 µmol/L (Note 4)

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

SPINLAB 180 APPLICATION

Name	IRON	Ref. male low	65 µg/dL
Abbr. Name	FE	Ref. male high	175 µg/dL
Mode	Endpoint	Ref. female low	40 µg/dL
Wavelength	578 nm	Ref. female high	150 µg/dL
Units	µg/dL	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	10 µg/dL	Panic value low	*
High Conc.	1000 µg/dL	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	250 µL		
rerun volume	250 µL		
Sample			
normal volume	30 µL		
rerun volume	15 µL		
R3 bottle (mL)	5 mL		
normal volume	12.0 µL		
rerun volume	12.0 µL		
Predilución	No		
Incubation time	6.5 min.		
Factor			
Reagent blank	Yes		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		

INTERFERENCES

 Hemolyzed samples are rejected, since erythrocytes contain iron and therefore falsely elevate the serum results^{1,2}.

 A list of drugs and other interfering substances with iron determination has been reported by Young et. al^{3,4}.

PERFORMANCE CHARACTERISTICS
Measuring range: From *detection limit* of 0,850 µg/dL to *linearity limit* of 1000 µg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (µg/dL)	SD	CV (%)	
Mean (µg/dL)	113	250	111	249
SD	0,89	0,72	3,51	6,29
CV (%)	0,79	0,29	3,17	2,52

Sensitivity: 1 µg/dL = 0,00104 A.

Accuracy: Results obtained using SPINREACT reagents did not show systematic differences when compared with other commercial reagents.

The results obtained using 50 samples were the following:

 Correlation coefficient (r)²: 0,9934.

Regression equation: y= 1,0243x - 3,877.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. It is recommended to use disposable material. If glassware is used the material should be soaking for 6 h in diluted HCl (20% v/v) and then thoroughly rinsed with distilled water and dried before use.
2. Use clean disposable pipette tips for its dispensation.
3. Strongly method dependent.
4. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

1. Perrotta G. Iron and iron-binding capacity. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1063-1065.
2. Itano M M D. Cap Serum Iron Survey 1978 (70): 516-522.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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BIBLIOGRAPHIE

- Perrotta G. Iron and iron-binding capacity. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1063-1065.
- Itano M M D. Cap Serum Iron Survey 1978 (70): 516-522.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: SP1001247	Cont.	R1: 10 x 25 mL, R2: 10 x 500 mg, R3: 3 x 5 mL
----------------	-------	---

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