Quantitative determination of HDL cholesterol




IVD

Store at 2-8°C


PRINCIPLE OF THE METHOD

Directly determination of serum HDLc (high-density lipoprotein cholesterol) levels without the need for any pre-treatment or centrifugation of the sample^{3,5}
The assay takes place in two steps.

- 1° Elimination of lipoprotein no-HDL

- 2° Measurement of HDLc

The intensity of the color formed is proportional to the HDLc concentration in the sample.

CLINICAL SIGNIFICANCE

HDL particles are high-density lipoproteins that transport cholesterol from the body tissues to the liver. Since HDL can remove cholesterol from the arteries and carry it back to the liver for their excretion, HDL is known as "good cholesterol" because high levels are thought to lower the risk of heart disease and coronary artery disease.

A low HDL cholesterol levels, is considered a greater heart disease risk^{1,2,4}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 6,6	100 mM
	N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline (HDAOS)	0,7 mM
	Cholesterol Esterase	≥ 800 U/L
	Cholesterol oxidase	≥ 500 U/L
	Catalase	≥ 300 U/L
	Ascorbic oxidase	≥ 3000 U/L
R 2	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 7,0	1,1 mmol/L
	4 - Aminoantipyrine (4-AA)	100 mM
	Peroxidase	≥ 3500 U/L
	HDLc/ LDLc CAL	Calibrator. Lyophilized human serum.

PRECAUTIONS

HDLc/ LDLc CAL

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

PREPARATION

DUAL MODE: Ready to use.

- **HDLc/ LDLc CAL:** Dissolve the contents with 1 mL. of distilled water. Cap vial and mix gently to dissolve contents.

STORAGE AND STABILITY (Note 1)

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not freeze the reagents.

- **HDLc/ LDLc CAL:** Once reconstitute 30 hours at 20-25°C, 2 weeks at 2-8°C or 3 months at -20°C.

- Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 600 nm.

- Matched cuvettes 1.0 cm light path.

- General laboratory equipment.

SAMPLES

Serum, heparinized plasma or EDTA plasma. If any samples show precipitates, centrifuge before using.

Stability of the sample: 6 days at 2-8°C and 1 year when stored at -70°C.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINLAB 180

Name	HDL Cholesterol	Ref. male low	35.0
Abbr. Name	HDL	Ref. male high	55.0
Mode	Twopoint	Ref. female low	45.0
Wavelength	578 nm	Ref. female high	65.0
Units	mg/dL	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	3.0 mg/dL	Panic value low	*
High Conc.	120.0 mg/dL	Panic value high	*
Calibrator name	*	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	225 µL		
rerun volume	225 µL		
Sample			
normal volume	3 µL		
rerun volume	2 µL		
R2 bottle (mL)	5 mL		
normal volume	75.0 µL		
rerun volume	75.0 µL		
Predilution	No		
Slope blank	No		
Point one, two	24, 236 sec.		
Factor	**		
Reagent blank	Yes (0.000)		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Sustrate depletion	3.000 Abs		

REFERENCE VALUES

	Men	Women
Low risk	> 50 mg/dL	> 60 mg/dL
Normal risk	35 - 50 mg/dL	45 - 60 mg/dL
High risk	< 35 mg/dL	< 45 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 5,0 mg/dL to *linearity limit* of 151 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

Mean (mg/dL)	Intra-assay (n=20)		Inter-assay (n=20)	
	28,0	76,1	27,5	75,3
SD	0,25	0,81	1,26	2,04
CV (%)	0,89	1,06	4,60	2,71

Sensitivity: 1mg/L = 0,001399 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,938.

Regression equation: y = 0,9825x - 1,41606.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. The reagent 2 presents yellowish coloration due to the peroxidase, but it does not affect its functionality.

2. **SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

BIBLIOGRAPHY

- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 - 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 – 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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Ref: SP1001096	Cont.	R1: 10 x 24 mL, R2:10 x 8 mL, CAL: 1 x 1 mL
Ref: SP1001097		R1: 2 x 24 mL, R2: 2 x 8 mL, CAL: 1 x 1 mL

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BIBLIOGRAPHIE

1. National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
2. Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
3. Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 - 364
4. Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
5. Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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Ref: SP1001096	Cont.	R1: 10 x 24 mL, R2:10 x 8 mL, CAL: 1 x 1 mL
Ref: SP1001097		R1: 2 x 24 mL, R2: 2 x 8 mL, CAL: 1 x 1 mL

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