LDL Cholesterol D

Enzymatic colorimetric. Liquid

Quantitative determination of LDL cholesterol IVD




Store at 2-8°C

PRINCIPLE OF THE METHOD




Direct determination of serum LDLc (low-density lipoprotein cholesterol) levels without the need for any pre-treatment or centrifugation steps.

The assay takes place in two steps.

- 1° Elimination of lipoprotein no-LDL

- 2° Measurement of LDLc

The intensity of the color formed is proportional to the LDLc concentration in the sample.

CLINICAL SIGNIFICANCE

The LDLc particle is lipoproteins that transport cholesterol to the cells. Often called "bad cholesterol" because high levels are risk factor for coronary heart disease and are associated with obesity, diabetes and nephrosis^{1,2,9}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1 Enzymes	PIPES Buffer	50 mmol/L
	Cholesterol esterase (CHE)	≥600 U/L
	Cholesterol oxidase (CHOD)	≥500 U/L
	Catalase	≥600 KU/L
R 2 Enzymes	TOOS	2 mmol/L
	PIPES Buffer	50 mmol/L
	4 - Aminoantipyrine (4-AA)	4 mmol/L
	Peroxidase (POD)	≥ 4 KU/L
HDLc/LDLc CAL	Calibrator. Lyophilized human serum	

PRECAUTIONS

HDLc/LDLc CAL: Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

PREPARATION

DUAL MODE: Ready to use.

- **HDLc/LDLc CAL:** Dissolve the contents with 1 mL of distilled water. Cap vial and mix gently to dissolve contents.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use.

- **R 1 and R 2:** Once opened is stable 4 weeks at 2-8°C.

- **HDLc/LDLc CAL:** Once reconstitute 30 hours at 20-25°C, 2 weeks at 2-8°C or 3 months -20°C.

- Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 600 nm.

- Matched cuvettes 1,0 cm light path.

- General laboratory equipment.

SAMPLES

Serum, heparinized plasma or EDTA plasma. If any sample show precipitates, centrifuge before using⁵.

Serum stable 6 days at 2-8°C. Do not freeze the samples.

REFERENCE VALUES^{6,7,8}

Optimal	< 100 mg/dL
Near or above optimal	100-129 mg/dL
Borderline high	130-160 mg/dL
High	> 160 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

APPLICATION SPINLAB 180

Name	LDL Colesterol	Ref. male low	49.0
Abbr. Name	LDL	Ref. male high	172.0
Mode	Twopoint	Ref. female low	63.0
Wavelength	578 nm	Ref. female high	167.0
Units	mg/dL	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	0.0 mg/dL	Panic value low	*
High Conc.	250.0 mg/dL	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	225 µL		
rerun volume	225 µL		
Sample			
normal volume	3.0 µL		
rerun volume	2.0 µL		
R2 bottle (mL)	5 mL		
normal volume	75.0 µL		
rerun volume	75.0 µL		
Predilution	No		
Slope blank	No		
Point one, two	24, 236 sec.		
Factor			
Reagent blank	Yes (0.000)		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substrate depletion	3.000 Abs		

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINCONTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 10 mg/dL to *linearity limit* of 976 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

Media (mg/dL)	Intraserie (n= 20)		Interserie (n= 20)	
		31,4	67,8	32,1
SD	0,42	1,11	0,92	2,02
CV (%)	1,35	1,64	2,87	2,97

Sensibility: 1 mg/dL = 0,001784 (A).

Accuracy^{10,11}: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,99123.

Regression equation: y = 0,914x + 1,58283.

The results of the performance characteristics depend on the analyzer used.

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APPLICATION AU SPINLAB 180

Name	LDL Colesterol	Ref. male low	49.0
Abbr. Name	LDL	Ref. male high	172.0
Mode	Twopoint	Ref. female low	63.0
Wavelength	578 nm	Ref. female high	167.0
Units	mg/dL	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	0.0 mg/dL	Panic value low	*
High Conc.	250.0 mg/dL	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	225 µL		
rerun volume	225 µL		
Sample			
normal volume	3.0 µL		
rerun volume	2.0 µL		
R2 bottle (mL)	5 mL		
normal volume	75.0 µL		
rerun volume	75.0 µL		
Predilution	No		
Slope blank	No		
Point one, two	24, 236 sec.		
Factor			
Reagent blank	Yes (0.000)		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Sustrate depletion	3.000 Abs		

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