Quantitative determination of creatinine IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

In the first reaction, creatinase and sarcosine oxidase were used in the enzymatic hydrolysis of endogenous creatine to produce hydrogen peroxide, that is eliminated by catalase. In the second reaction, the catalase is inhibited by sodium azide, and creatinase and 4- aminoantipyrine (4-AA) were added, and only the creatine generated from creatinine by creatininase was hydrolyzed sequentially by creatinase and sarcosine oxidase to produce hydrogen peroxide. This newly-formed hydrogen peroxide was measured in a coupled reaction catalyzed by peroxidase, with N-ethyl-n-sulphopropyl-mtoluidine (TOPS)/4-AA as a chromogen.

CLINICAL SIGNIFICANCE

Creatinine measurements are used in the diagnosis and treatment of renal diseases, in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes. Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

REAGENTS

R 1	MOPS 25 mmol/L, TOPS 0.5 mmol/L, Creatinase 10 KU/L, Sarcosine Oxidase 5 KU/L, Catalase 3 KU/L, EDTA 1mmol/L, pH 7,5.
R 2	MOPS 90 mmol/L, Creatininase 30 KU/L, peroxidase 10 KU/L, pH 7,5. Azida sódica 0,5 g/L.

PREPARATION

R1 and R2 are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.
R1 and R2 are stable 8 weeks after opening bottle.

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer measuring at 545±20 nm
- Cell holder thermostatable at 37°C
- General laboratory equipment.

SAMPLES

- Serum or plasma¹.
- Urine (24 h)¹: Dilute fresh urine 1/50 with distilled water. Multiply the result by 50 (sample dilution factor). Creatinine is stable 1 day at 2-8°C.

REFERENCE VALUES¹

Serum or plasma:

Men 0,9 - 1,3 mg/dL

Women 0,6 - 1,1 mg/dL

Urine:

Men 14- 26 mg/Kg/24 h

Women 11-20 mg/Kg/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

INTERFERENCES

No interferences were observed with hemoglobin until 5 g/dL, bilirubin 40 mg/dL.
Other drugs and substances may interfere^{3,4}.

APPLICATION SPINLAB 180

Name	Creatinine	Ref. male low	0.70
Abbr. Name	CREA	Ref. male high	1.40
Mode	Twopoints	Ref. female low	0.60
Wavelength	546 nm	Ref. female high	1.10
Units	mg/dL	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	0.03 mg/dL	Control 1	*
High Conc.	150.0 mg/dL	Control 2	*
Calibrator name	CAL	Control 3	*
Cut-off	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	270 µL		
rerun volume	270 µL		
Sample			
normal volume	6.0 µL		
rerun volume	6.0 µL		
R2 bottle (mL)	25 mL		
normal volume	90.0 µL		
rerun volume	90.0 µL		
Predilución	No		
Point one, two	24, 236 sec.		
Factor	**		
Reagent blank	Yes (0.000)		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
R. Abs. Dev.	3.000 Abs		

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 y 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 0,00 mg/dL to *linearity limit* of 180 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
Mean (mg/dL)	0,87	3,82	0,87	3,75
SD	0,01	0,06	0,02	0,06
CV (%)	1,63	1,44	2,31	1,72

Sensitivity: 1 mg/dL = 0,0226 (ΔA)

Accuracy: Results obtained using SPINREACT these reagents did not show systematic differences when compared with other commercial reagents or with HPLC method.

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,9730

Regression equation: y= 1,066x - 0,020.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

1. Fossati et al. Clin Chem 1983;29:1494-1496.
2. Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co.,1999.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.

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APPLICATION AU SPINLAB 180

Name	Creatinine	Ref. male low	0.70
Abbr. Name	CREA	Ref. male high	1.40
Mode	Twopoints	Ref. female low	0.60
Wavelength	546 nm	Ref. female high	1.10
Units	mg/dL	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	0.03 mg/dL	Control 1	*
High Conc.	150.0 mg/dL	Control 2	*
Calibrator name	CAL	Control 3	*
Cut-off	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	270 µL		
rerun volume	270 µL		
Sample			
normal volume	6.0 µL		
rerun volume	6.0 µL		
Vol. repet.	3.0 µL		
R2 bottle (mL)	25 mL		
normal volume	90.0 µL		
rerun volume	90.0 µL		
Predilución	No		
Point one, two	24, 236 sec.		
Factor	**		
Reagent blank	Yes (0.000)		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		

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