# ALP-LQ (Alkaline phosphatase)


p-Nitrophenylphosphate. kinetic. Liquid. DGKC

## Quantitative determination of alkaline phosphatase (ALP) IVD

Store at 2-8°C

### PRINCIPLE OF THE METHOD

Alkaline phosphatase (ALP) catalyses the hydrolysis of p-nitrophenyl phosphate at pH 10.4, liberating p-nitrophenol and phosphate, according to the following reaction:

The rate of p-nitrophenol formation, measured photometrically, is proportional to the catalytic concentration of alkaline phosphatase present in the sample<sup>1,2</sup>.

### CLINICAL SIGNIFICANCE

Alkaline phosphatase is an enzyme present in almost all tissues of the organism, being particularly high in bone, liver, placenta, intestine and kidney. Both increases and decreases of plasma ALP are of importance clinically.

Causes of increased plasma ALP: Paget's disease of bone, obstructive liver disease, hepatitis, hepatotoxicity caused by drugs or osteomalacia.

Causes of decreased plasma ALP: Cretinism and vitamin C deficiency<sup>1,5,6</sup>.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

### REAGENTS

<b>R 1</b>	Diethanolamine (DEA) pH 10,4	1 mmol/L
Buffer	Magnesium chloride	0,5 mmol/L
<b>R 2</b>	p-Nitrophenylphosphate (pNPP)	10 mmol/L
Substrate		
<b>Optional</b>	SPINTROL H CAL	

### PRECAUTIONS

R1: H315- Contains Diethanolamine (HN(CH<sub>2</sub>CH<sub>2</sub>OH)<sub>2</sub>). Causes skin irritation. H318- Causes serious eye damage. H373- May cause damage to organs through prolonged or repeated exposure. Follow the precautionary statements given in MSDS and label of the product.

### PREPARATION

DUAL MODE: Ready to use.

MONO MODE: Pour reagent 2 content over reagent 1. Mix thoroughly avoiding foam forming and it will be ready to use (WR).

Working reagent (WR) stability: 1 month at 2-8°C or 10 days at room temperature (15-25°C).

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

### Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 405 nm  $\geq 1.50$ .

### ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 405 nm.
- Thermostatic bath at 25°C, 30°C or 37°C ( $\pm 0.1^\circ\text{C}$ )
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

### SAMPLES

Serum or heparinized plasma<sup>1</sup>. Use unhemolyzed serum, separated from the clot as soon as possible. Stability: 3 days at 2-8°C.

### INTERFERENCES

Fluoride, oxalate, citrate and EDTA inhibit alkaline phosphate activity and should therefore not be used as anticoagulants. Haemolyses interferes due to the high concentration of alkaline phosphatase in red cells<sup>1,2</sup>.

A list of drugs and other interfering substances with acid phosphatase determination has been reported<sup>3,4</sup>.

### QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

## APPLICATION SPINLAB 180

Name	ALP	Ref. male low	98
Abbr. Name	ALP	Ref. male high	279
Mode	Kinetic	Ref. female low	98
Wavelength	405 nm	Ref. female high	279
Units	U/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	4 U/L	Panic value low	*
High Conc.	825 U/L	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>		<b>MONO MODE</b>	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	312 $\mu\text{L}$	normal volume	390 $\mu\text{L}$
rerun volume	312 $\mu\text{L}$	rerun volume	390 $\mu\text{L}$
Sample		Sample	
normal volume	6.5 $\mu\text{L}$	normal volume	6.5 $\mu\text{L}$
rerun volume	6 $\mu\text{L}$	rerun volume	6 $\mu\text{L}$
R2 bottle (mL)	5 mL		
normal volume	78 $\mu\text{L}$		
rerun volume	78 $\mu\text{L}$		
Predilución	No		
Slope blank	No		
Delay, min. time	24.212 sec.	Delay, min. time	24.212 sec.
Linearity limit	10.0 %	Linearity limit	10.0 %
Factor		Factor	
Reagent blank	No	Reagent blank	No
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

### REFERENCE VALUES<sup>1</sup>

	25°C	30°C	37°C
Children (1-14 years)	< 400 U/L	< 480 U/L	< 645 U/L
Adults	60 - 170 U/L	73 - 207 U/L	98 - 279 U/L

Factors affecting ALP activities in a normal population include exercise, periods of repaid growth in children and pregnancy.

These values are for orientation purpose; each laboratory should establish its own reference range.

### PERFORMANCE CHARACTERISTICS

**Measuring range:** From *detection limit* of 0,6845 U/L to *linearity limit* of 1200 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

### Precision:

Mean (U/L)	Intra-assay (n=20)		Inter-assay (n=20)	
	174	443	175	434
SD	0,72	1,56	6,88	11,93
CV (%)	0,41	0,35	3,93	2,75

**Sensitivity:** 1 U/L = 0,0003  $\Delta\text{A}/\text{min}$ .

**Accuracy:** Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x). The results obtained using 50 samples were the following:

Correlation coefficient (r)<sup>2</sup>: 0,99938.

Regression equation:  $y = 1,025x - 1,105$ .

The results of the performance characteristics depend on the analyzer used.

### NOTES

**SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.**

### BIBLIOGRAPHY

- Wenger C. et al. Alkaline phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1094-1098.
- Rosalki S et al. Clin Chem 1993; 39/4: 648-652.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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## APLICACION AL SPINLAB 180

Name	ALP	Ref. male low	98
Abbr. Name	ALP	Ref. male high	279
Mode	Kinetic	Ref. female low	98
Wavelength	405 nm	Ref. female high	279
Units	U/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	4 U/L	Panic value low	*
High Conc.	825 U/L	Panic value high	*
Calibrator name	CAL	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>		<b>MONO MODE</b>	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	312 $\mu$ L	normal volume	390 $\mu$ L
rerun volume	312 $\mu$ L	rerun volume	390 $\mu$ L
Sample		Sample	
normal volume	6,5 $\mu$ L	normal volume	6,5 $\mu$ L
rerun volume	6 $\mu$ L	rerun volume	6 $\mu$ L
R2 bottle (mL)	5 mL		
normal volume	78 $\mu$ L		
rerun volume	78 $\mu$ L		
Predilución	No		
Slope blank	No		
Delay, min. time	24.212 sec.	Delay, min. time	24.212 sec.
Linearity limit	10.0 %	Linearity limit	10.0 %
Factor		Factor	
Reagent blank	No	Reagent blank	No
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

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	25°C	30°C	37°C
Enfants (1-14 ans)	< 400 U/L	< 480 U/L	< 645 U/L
Adultes	60 - 170 U/L	73 - 207 U/L	98 - 279 U/L

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### BIBLIOGRAPHIE

1. Wenger C. et al. Alkaline phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1094-1098.
2. Rosalki S et al. Clin Chem 1993; 39/4: 648-652.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Ref: SP41233	Cont.	R 1: 10 x 20 mL
		R 2: 10 x 5 mL

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### BIBLIOGRAFIA

- Wenger C. et al. Alkaline phosphatase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1094-1098.
- Rosalki S et al. Clin Chem 1993; 39/4: 648-652.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995

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Ref: SP41233

Con

R1: 10 x 20 mL

R2: 10 x 5 mL