


Quantitative determination of creatine kinase liquid (CK) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Kinetic determination of the creatine kinase based upon IFCC and DGKC recommendations.

Creatine kinase (CK) catalyses the reversible transfer of a phosphate group from phosphocreatine to ADP. This reaction is coupled to those catalysed by hexokinase (HK) and glucose-6-phosphate dehydrogenase (G6P-DH):

 The rate of NADPH formation, measured photometrically, is proportional to the catalytic concentration of CK present in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Creatine kinase is a cellular enzyme with wide tissue distribution in the body. Its physiological role is associated with adenosine triphosphate (ATP) generation for contractile or transport systems.

 Elevated CK values are observed in diseases of skeletal muscle and after myocardial infarction^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Imidazol, pH 6.7	125 mmol/L
	D-Glucose	25 mmol/L
	N-Acetyl-L-Cysteine	25 mmol/L
	Magnesium acetate	12.5 mmol/L
	NADP	2.52 mmol/L
	EDTA	2.02 mmol/L
R 2	ADP	15.2 mmol/L
	AMP	15.2 mmol/L
	di-Adenosine-5- pentaphosphate	25 mmol/L
		103 mmol/L
	Glucose-6-phosphate dehydrogenase	≥8 800 U/L
	Creatine phosphate	250 mmol/L

Optional

CK-Nac / CK-MB CONTROL	Lyophilized human serum	Ref: 1002260
-------------------------------	-------------------------	--------------

PRECAUTIONS

R1: H360-May damage fertility or the unborn child. Contains: Imidazole (C3H4N2)

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

R1 and R2 ready to be used.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use the tablets if appears broken.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm ≥ 1.60.

ADDITIONAL EQUIPMENT

- SPINLAB 180 autoanalyzer
- Laboratory equipment.

SAMPLES

Serum free of hemolysis or heparin plasma.

Stability 7 days at 2-8°C, protected from light.

The creatin kinase activity decreases 10% after 1 day at 2-5°C or after 1 hour at 15-25°C.

REFERENCE VALUES¹

	25°C	30°C	37°C
Men, up to	80 U/L	130 U/L	195 U/L
Women, up to	70 U/L	110 U/L	170 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINLAB 180

Name	CK-NAC	Ref. male low	24
Abbr. Name	CK	Ref. male high	195
Mode	Kinetic	Ref. female low	24
Wavelength	340 nm	Ref. female high	170
Units	U/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	10 U/L	Control 1	*
High Conc.	2000 U/L	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	240 µL	normal volume	300 µL
rerun volume	240 µL	rerun volume	300 µL
Sample		Sample	
normal volume	5 µL	normal volume	6 µL
rerun volume	5 µL	rerun volume	6 µL
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60		
Predilution	No		
Slope blank	No		
Delay, min. time	103.133 sec.	Delay, min. time	183, 177 sec.
Linearity limit	10.0 %	Linearity limit	10.0 %
Factor	-	Factor	-
Vol. repet.			30.0 µL
Reagent blank	No	Reagent blank	No
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

PERFORMANCE CHARACTERISTICS
Measuring range: From detection limit of 2,12 U/L to linearity limit of 2000 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay		Inter-assay	
	Mean (U/L)	SD	CV (%)	
Mean (U/L)	147	494	145	485
SD	1,23	3,60	2,91	8,97
CV (%)	0,84	0,73	2,01	1,85

Sensitivity: 1 U/L = 0,00012 ΔA/min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained were the following:

 Correlation coefficient (r)²: 0,9995

Regression equation: y = 1,0846x - 0,3512.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

 No interferences were observed with glucose until 7 g/L, hemoglobin until 5 g/L and triglycerides 7 mmol/L. A list of drugs and other interfering substances with CK determination has been reported by Young et. al^{3,4}.

BIBLIOGRAPHY

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommendation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommendation pour la mesure de la concentration catalytique de la créatinine kinase dans le sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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APPLICATION AU SPINLAB 180

Name	CK-NAC	Ref. male low	24
Abbr. Name	CK	Ref. male high	195
Mode	Kinetic	Ref. female low	24
Wavelength	340 nm	Ref. female high	170
Units	U/L	Ref. Ped. Low	*
Decimals	0	Ref. Ped. High	*
Low Conc.	10 U/L	Control 1	*
High Conc.	2000 U/L	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	240 µL	normal volume	300 µL
rerun volume	240 µL	rerun volume	300 µL
Sample		Sample	
normal volume	5 µL	normal volume	6µL
rerun volume	5 µL	rerun volume	6µL
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60 µL		
Predilution	No		
Slope blank	No		
Delay, min. time	103.133 sec.	Delay, min. time	183, 177 sec.
Linearity limit	10.0 %	Linearity limit	10.0 %
Factor	-	Factor	-
Reagent blank	No	Reagent blank	No
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

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BIBLIOGRAPHIE

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la créatinine kinase dans le sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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BIBLIOGRAFIA

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W. et al. Creatine kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommendation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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Ref: SP41250

Cont.

R1: 10 x 20 mL

R2: 10 x 5 mL