CK-MB-LQ (Creatine kinase-MB)

Anti CK-M. Immunoinhibition. Kinetic UV. Liquid

Quantitative determination of creatine kinase MB (CK-MB)

IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The procedure involves measurement of CK activity in the presence of an antibody to CK-M monomer. This antibody completely inhibits the activity of CK-MM and half of the activity of CK-MB while not affecting the B subunit activity of CK-MB and CK-BB. Then it's used the CK method to quantitatively determine CK-B activity^{1,2}. The CK-MB activity is obtained by multiplying the CK-B activity by two.

CLINICAL SIGNIFICANCE

CK-MB is an enzyme formed by the association of two subunits from muscle (M) and nerve cells (B). CK-MB is usually present in serum at low concentration; it is increased after an acute infarct of myocardium and later descends at normal levels. Also, is increased, rarely, in skeletal muscle damage^{5,6}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Imidazol, pH 6.7	125 mmol/L
	D-Glucose	25 mmol/L
	N-Acetyl-L-Cysteine	25 mmol/L
	Magnesium acetate	12,5 mmol/L
	NADP	2,52 mmol/L
	EDTA	2,02 mmol/L
	Hexokinase	≥6 800 U/L
Anti-human polyclonal CK-M antibody (sheep) sufficient to inhibit up to 2 000 U/L of CK-MM		
R 2	ADP	15,2 mmol/L
	AMP	25 mmol/L
	di-Adenosine-5- pentaphosphate	103 mmol/L
	Glucose-6-phosphate dehydrogenase	≥8 800 U/L
	Creatine phosphate	250 mmol/L

Optional

CK-Nac / CK-MB CONTROL	Lyophilized human serum	Ref: 1002260
------------------------	-------------------------	--------------

PREPARATION

MONO MODE: Pour reagent 2 content over reagent 1. Mix thoroughly avoiding foam forming and it will be ready to use (WR).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm ≥ 1.60.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 340 nm.
- Thermostatic bath at 25°C, 30°C ó 37° C (± 0.1°C).
- Matched cuvettes 1.0 cm light path.
- General laboratory equipment.

LIMITATION OF THE PROCEDURE

- 1- The method will also measure any CK-BB isoenzyme present in serum. The activity of the isoenzyme is usually negligible, however, if a significant amount of CK-BB activity is present the CK-MB activity will be overestimated.
- 2- A macro form of BB (immunoglobulin complexed) has been observed which will be measured as B in the assay. If the measured CK-B activity exceeds 20% of the total CK activity, the presence of macro BB should be suspected.

SAMPLES

Serum free of hemolysis or heparin plasma¹: Stability 7 days at 2-8°C, protected from light.

CK-MB activity decreases a 10% after 24 hours at 4°C or 1 hour at 25°C.

REFERENCE VALUES

The suspicion of myocardial damage is based on the three following factors:

CK-MB	25°C	30°C	37°C	
	> 10 U/L	> 15 U/L	> 24 U/L	
TOTAL CK	25°C	30°C	37°C	
	Men, up to	80 U/L	130 U/L	195 U/L
	Women, up to	70 U/L	110 U/L	170 U/L

CK - MB Activity
 CK Total Activity × 100: 6 - 25 % CK - MB Activity in the sample

These values are for orientation purpose; each laboratory should establish its own reference range.

APPLICATION TO SPINLAB 180

Name	CK-MB.	Ref. male low	0.0
Abbr. Name	CKMB	Ref. male high	24.0
Mode	Kinetic	Ref. female low	0.0
Wavelength	340 nm	Ref. female high	24.0
Units	U/L	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	2.0 U/L	Control 1	*
High Conc.	1000.0 U/L	Control 2	*
Calibrator name	CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE		MONO MODE	
Sample blank	No	Sample blank	No
R1 bottle (mL)	25 mL	R1 bottle (mL)	25 mL
normal volume	300 µL	normal volume	300 µL
rerun volume	300 µL	rerun volume	300 µL
Sample		Sample	
normal volume	12.0 µL	normal volume	12.0 µL
rerun volume	6.0 µL	rerun volume	6.0 µL
R2 bottle (mL)	5 mL		
normal volume	0 µL		
rerun volume	0 µL		
Predilución	No		
Slope blank	No		
Delay, min. time	360, 159 sec.	Delay, min. time	304, 118 sec.
Linearity limit	10.0 %	Linearity limit	10.0 %
Factor		Factor	
Reagent blank	No	Reagent blank	No
Low Absorbance	-0.100 Abs	Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs	High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs	R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs	R. Abs. H. Limit	3.000 Abs
R. Abs. Deviation	3.000 Abs	R. Abs. Deviation	3.000 Abs

QUALITY CONTROL

CK-Nac/CK-MB specific control sera (Ref. 1002260) are recommended to monitor the performance of assay procedures.

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 1,9 U/L to *linearity limit* of 1000 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/1 with NaCl 9 g/L and multiply the result by 2.

Precision:

Mean (U/L)	Intra-assay		Inter-assay	
	SD	33,7	166,5	31,3
CV (%)	1,00	3,76	1,19	3,47
	2,96	2,26	3,81	2,15

Sensitivity: 1 U/L = 0,000134 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained were the following:

Correlation coefficient (r)²: 0,999.

Regression equation: y = 0,976 x - 0,269.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.
8. Neumeier, D., Prellwitz, W., Würzburg, U. et coll. Determination of creatine kinase isoenzyme MB activity in serum using immunological inhibition of creatine kinase M subunit activity. Activity kinetics and diagnostic significance in myocardial infarction. Clin. Chim. Acta, 73, (1976), 445.

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- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
- Gerhardt W. et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
- Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.
- Neumeier, D., Prellwitz, W., Würzburg, U. et coll. Determination of creatine kinase isoenzyme MB activity in serum using immunological inhibition of creatine kinase M subunit activity. Activity kinetics and diagnostic significance in myocardial infarction. Clin. Chim. Acta, 73, (1976), 445.

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BIBLIOGRAPHIE

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la créatine kinase dans le sérum humain. Ann. Biol. Clin.,40, (1482), 87.
8. Neumeier, D., Prellwitz, W., Würzburg, U. et coll. Determination of creatine kinase isoenzyme MB activity in serum using immunological inhibition of creatine kinase M subunit activity. Activity kinetics and diagnostic significance in myocardial infarction, Clin. Chim. Acta, 73, (1976), 445.

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Ref: SP41254	Cont.	R1: 10 x 20 mL
		R2: 10 x 5 mL

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BIBLIOGRAFIA

- Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
- Gerhardt W. et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
- Mathieu M. et coll. Recommendation pour la mesure de la concentration catalytique de la creatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.
- Neumeier, D., Prellwitz, W., Würzburg, U. et coll. Determination of creatine kinase isoenzyme MB activity in serum using immunological inhibition of creatine kinase M subunit activity. Activity kinetics and diagnostic significance in myocardial infarction, Clin. Chim. Acta, 73, (1976), 445.

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Ref: SP41254

Cont.

R1: 10 x 20 mL

R2: 10 x 5 mL