Quantitative determination of α₁-antitrypsin (α₁-ATRYP)

IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Anti- α₁-antitrypsin antibodies when mixed with samples containing α₁-antitrypsin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the α₁-antitrypsin concentration of the patient sample, that can be quantified by comparison from a calibrator of known α₁-antitrypsin concentration.

CLINICAL SIGNIFICANCE

α₁-antitrypsin is a glycoprotein synthesized in the hepatic parenchyma cells that circulates in the bloodstream. It is the second highest proteinase inhibitor in plasma after α₂-macroglobulin. α₁-antitrypsin is a strong reactor with elastase, skin collagenase, chymotrypsin, plasmin, and thrombin, and also shows inhibitory activity against fungal and leukocytic proteases.

α₁-antitrypsin deficiency is an inherited disorder, and occurs when both parents pass on an abnormal gene (PiZ) to their child. This deficiency is associated with a very high risk for development of pulmonary emphysema and diseases of the liver - neonatal cholestasis, hepatitis, cirrhosis and hepatocellular carcinoma.

α₁-antitrypsin increases in inflammatory or necrosis process. Serum levels begin to rise after approximately 24 hours and peak at 3 or 4 days if the insult is acute and short-lived.

REAGENTS

R 1 Diluent	Tris buffer 20 mmol/L, PEG 8000, pH 8.3. Preservative.
R 2 Antibody	Goat serum, anti-human α ₁ -antitrypsin, pH 7.5. Preservative.
Optional	Ref: 1102003 PROT CAL.

CALIBRATION

The assay has been standardized against the Reference Material ERM-DA470k/IFCC. It must be used the PROT CAL to calibrate the reagent. The reagent (both monoreagent and bireagent) should be recalibrated every month, when the controls are out of specifications, and when changing the reagent lot or the instrument settings.

PREPARATION

Reagents: Ready to use.

Calibration Curve: Prepare the following PROT CAL dilutions in NaCl 9 g/L as diluent. Multiply the concentration of the α₁-antitrypsin calibrator by the corresponding factor stated in table below to obtain the α₁-antitrypsin concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator (μL)	--	10	25	50	75	100
NaCl 9 g/L (μL)	100	90	75	50	25	-
Factor	0	0.1	0.25	0.5	0.75	1.0

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

Reagent deterioration: The presence of particles and turbidity. Do not freeze; frozen Antibody or Diluent could change the functionality of the test.

ADDITIONAL EQUIPMENT

- Spinlab 180 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum or plasma. EDTA or heparin should be used as anticoagulant. Stable 7 days at 2-8°C or 3 months at -20°C.

REFERENCE VALUES

Newborn: Between 124 - 348 mg/dL.

Adults: 90 - 200 mg/dL.

Each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of manual and automated assay procedures. Spinreact PROT CONTROL (Ref.:1102004) is available. Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

SPINLAB 180 APPLICATION

Name	Antitrypsin	Ref. male low	90 mg/dL
Abbr. Name	ATRYP	Ref. male high	200 mg/dL
Mode	Twopoints	Ref. female low	90 mg/dL
Wavelength	340 nm	Ref. female high	200 mg/dL
Units	mg/dL	Ref. Ped. Low	124 mg/dL
Decimals	0	Ref. Ped. High	348 mg/dL
Low Conc.	16 mg/dL	Control 1	*
High Conc.	400 mg/dL	Control 2	*
Calibrator name	SP CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000

DUAL MODE	
Sample blank	No
R1 bottle (mL)	25 mL
normal volume	318 μL
rerun volume	318 μL
Sample	
normal volume	2.0 μL
rerun volume	2.0 μL
R2 bottle (mL)	5 mL
normal volume	80 μL
rerun volume	80 μL
Predilución	No
Slope blank	No
Point one,two	-3.130 sec.

Reagent blank	No
Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs
Substr.Depletion	3.000 Abs

PERFORMANCE CHARACTERISTICS

Measurement range: Up to 500 mg/dL under the described assay conditions. Samples with higher concentrations, should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit depends on the sample / reagent ratio. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

Detection Limit: Values less than 16 mg/dL give non-reproducible results.

Prozone effect: No prozone effect was detected upon 750 mg/dL

Precision: The reagent has been tested for 20 days, using three levels of serum in a EP5-based study.

EP5	CV (%)		
	34.39 mg/dl	92.7 mg/dl	181.8 mg/dl
Total	4.2%	2.6%	2.8%
Within Run	0.8%	1.1%	1.6%
Between Run	3.8%	2.4%	2.3%
Between Day	1.6%	0%	0%

Sensitivity: Δ 3.4 mA / mg/dL.

Accuracy: Results obtained using this reagent (y) were compared to those obtained using the Beckman Array 360 CE. 35 samples ranging from 70 to 250 mg/dL of α₁antitrypsin were assayed. The correlation coefficient (r) was 0.92 and the regression equation y = 0.8567x +26.5.

The results of the performance characteristics depend on the used analyzer.

NOTES

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Sharp HL. Hospital Practice; May 1971: 83-96
5. Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
7. Friedman and Young. Effects of disease on clin. laboratory tests, 3th ed. AACC Pres, 1997.

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APPLICATION AU SPINLAB 180

Name	Antitrypsin	Ref. male low	90 mg/dL
Abbr. Name	ATRYP	Ref. male high	200 mg/dL
Mode	Twopoints	Ref. female low	90 mg/dL
Wavelength	340 nm	Ref. female high	200 mg/dL
Units	mg/dL	Ref. Ped. Low	124 mg/dL
Decimals	0	Ref. Ped. High	348 mg/dL
Low Conc.	16 mg/dL	Control 1	*
High Conc.	400 mg/dL	Control 2	*
Calibrator name	SP CAL	Control 3	*
Prozone check	No	Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	318 μL		
rerun volume	318 μL		
Sample			
normal volume	2.0 μL		
rerun volume	2.0 μL		
R2 bottle (mL)	5 mL		
normal volume	80 μL		
rerun volume	80 μL		
Predilución	No		
Slope blank	No		
Point one,two	-3.130 sec.		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

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BIBLIOGRAPHIE

1. Clinical Guide to Laboratory Tests, Edited by NWTietz W B Saunders Co., Philadelphia, 483, 1983.
2. Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
3. PesceAJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
4. Sharp HL. Hospital Practice; May 1971: 83-96
5. Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
7. Friedman and Young. Effects of disease on clin. laboratory tests, 3th ed. AACC Pres, 1997.

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Ref: SP1102054	Cont.	R1: 2 x 20 mL
		R2: 2 x 5 mL

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BIBLIOGRAFIA

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W B Saunders Co., Philadelphia, 483, 1983.
- Dati F et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520.
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
- Sharp HL. Hospital Practice; May 1971: 83-96
- Carrel RW et al. Assays Med Biochem 1978; 4: 83-119
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Pres, 1995.
- Friedman and Young. Effects of disease on clin. laboratory tests, 3th ed. AACC Pres, 1997.

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Ref: SP1102054	Cont.	R1: 2 x 20 mL
		R2: 2 x 5 mL