**Quantitative determination of Rheumatoid Factors (RF) IVD**

Store 2 - 8°C.

**PRINCIPLE OF THE METHOD**

The RF-Turbilatex is a quantitative turbidimetric test for the measurement of RF in human serum or plasma.

Latex particles coated with human gammaglobulin are agglutinated when mixed with samples containing RF. The agglutination causes an absorbance change, dependent upon the RF contents of sample that can be quantified by comparison from a calibrator of known RF concentration.

**CLINICAL SIGNIFICANCE**

Rheumatoid factors are a group of antibodies directed to determinants in the Fc portion of the immunoglobulin G molecule. Although rheumatoid factors are found in a number of rheumatoid disorders, such as systemic lupus erythematosus (SLE) and Sjögren's syndrome, as well as in nonrheumatic conditions, its central role in clinic lies its utility as an aid in the diagnosis of rheumatoid arthritis (RA).

A study of the "American College of Rheumatology" shows that the 80.4% of RA patients were RF positive.

**REAGENTS**

<b>Diluent (R1)</b>	Tris buffer 20 mmol/L, pH 8.2. Preservative.
<b>Latex (R2)</b>	Latex particles coated with human gammaglobulin, pH 7.4. Preservative.
<b>RF-CAL</b>	Calibrator. Recombinant RF protein. The RF concentration is stated on the vial label.
<b>Optional</b>	Ref.:1102114 Control serum ASO/CRP/RF Level L Ref.:1102115 Control serum ASO/CRP/RF Level H.

**PRECAUTIONS**

R1 and R2: H317 - May cause an allergic skin reaction. Contains 2-Methylisothiazol-3(2H)-one (Proclin 950). Follow the precautionary statements given in MSDS and label of the product.

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

**CALIBRATION**

Use RF Calibrator Reference 1107007.

The sensitivity of the assay and the target value of the calibrator have been standardized against the International Reference NIBSC 64/2 (Rheumatoid Arthritis Serum) WHO.

Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

**PREPARATION**

Ready for use

**RF Calibrator:** Reconstitute (→) with 2.0 mL of distilled water. Mix gently and bring to room temperature for about 10 minutes before use.

**Calibration Curve:** Prepare the following RF calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the RF calibrator by the corresponding factor stated in table below to obtain the RF concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator RF (µL)	--	25	50	100	200	400
NaCl 9 g/L (µL)	400	375	350	300	200	-
Factor	0	0,0625	0,125	0,25	0,5	1,0

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Do not freeze; frozen latex and diluent could change the functionality of the test.

**Reagent deterioration:** Presence of particles (R1, R2) and turbidity (R1).

**Reconstituted calibrator:** Stable for 1 month at 2-8°C or 3 months at -20°C.

**ADDITIONAL EQUIPMENT**

- SPINLAB 180 autoanalyzer
- Laboratory equipment.

**SAMPLES**

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

**REFERENCE VALUES**

Normal values up to 8 IU/mL. Each laboratory should establish its own reference range.

**QUALITY CONTROL**

Control Sera are recommended before and after testing samples to monitor the performance of manual and automated assay procedures. It should be used SPINREACT Control ASO/CRP/RF Level L (Ref.: 1102114) and Level H (Ref.: 1102115).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**SPINLAB 180 APPLICATION**

Name	RF	Ref. male low	*
Abbr. Name	RF	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	620 nm	Ref. female high	*
Units	IU/mL	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	6 IU/mL	Panic value low	*
High Conc.	160 IU/mL	Panic value high	*
Calibrator name**	CAL RF**	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000

**DUAL MODE**

Sample blank	No
R1 bottle (mL)	25 mL
normal volume	240 µL
rerun volume	240 µL
Sample	
normal volume	3.0 µL
rerun volume	3.0 µL
R2 bottle (mL)	5 mL
normal volume	60 µL
rerun volume	60 µL
Predilución	No
Slope blank	No
Point one,two	-3, 130 sec.

Reagent blank	Yes
Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs
Substr.Depletion	3.000 Abs

\*\* Adjust maximum inaccuracy calibration % until the curve cover all the Abs found.

**PERFORMANCE CHARACTERISTICS**

**1. Limit detection:** Values less than 6 IU/mL give non-reproducible results.

**2. Measurement range:** 6-160 IU/mL, under the described assay conditions. Samples with higher concentrations should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity limit and measurement range depends on the sample to reagent/ratio, as well as the analyzer used. It will be higher by decreasing the sample volume, although the sensitivity of the test will be proportionally decreased.

**3. Prozone effect:** No prozone effect was detected upon 800 IU/mL.

**4. Sensitivity:** Δ 3.34 mA. IU/mL.

**5. Precision:** The reagent has been tested for 20 days, using three different FR concentrations in a EP5-based study.

EP5	CV (%)		
	35,8 IU/mL	78,05 IU/mL	123,26 IU/mL
Total	4,5%	4,1%	5,9%
Within Run	3,3%	2,6%	3,2%
Between Run	1,7%	2,3%	3,4%
Between Day	2,5%	2,1%	3,6%

**6. Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 41 samples of different concentrations of FR were assayed. The correlation coefficient (r<sup>2</sup>) was 0.91 and the regression equation  $y = 1,2042x + 3,1344$ . The results of the performance characteristics depend on the analyzer used.

**INTERFERENCES**

Hemoglobin (10 g/L), bilirubin (20 mg/dL) and lipemia (10 g/L), do not interfere. Other substances may interfere<sup>6</sup>.

**NOTES**

1. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

**BIBLIOGRAPHY**

1. Frederick Wolfe et al. Arthritis and Rheumatism 1991; 34: 951 - 960.
2. Robert W Dorner et al. Clinica Chimica Acta 1987; 167: 1-21.
3. Robert H Shmerling et al. The American Journal of Medicine 1991; 91: 528 - 534.
4. Vladimir Muié et al. Scand J Rheumatology 1972; 1: 181 - 187.
5. Paul R et al. Clin Chem 1979; 25/11: 1909 - 1914.
6. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AAC Press, 1995.

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Name	RF	Ref. male low	*
Abbr. Name	RF	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	620 nm	Ref. female high	*
Units	IU/mL	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	6 IU/mL	Panic value low	*
High Conc.	160 IU/mL	Panic value high	*
Calibrator name**	CAL RF**	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 µL		
rerun volume	240 µL		
Sample			
normal volume	3.0 µL		
rerun volume	3.0 µL		
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60 µL		
Predilución	No		
Slope blank	No		
Point one,two	-3, 130 sec.		
Reagent blank	Yes		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

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**BIBLIOGRAFIA**

- Frederick Wolfe et al. Arthritis and Rheumatism 1991; 34: 951 - 960.
- Robert W Dornier et al. Clinica Chimica Acta 1987; 167: 1-21.
- Robert H Shmerling et al. The American J. of Medicine 1991; 91: 528 - 534.
- Vladimir Muié et al. Scand J Rheumatology 1972; 1: 181 - 187.
- Paul R et al. Clin Chem 1979; 25/11: 1909 - 1914.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.

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Ref.: SP1107005	Cont.	R1. Diluíente: 2 x 20mL
		R2. Latex: 2 x 5 mL
		RF-CAL: 1 x 2 mL