β<sub>2</sub>-m TURBI

# β<sub>2</sub>-m turbilatex

Latex turbidimetry

## Quantitative determination of β<sub>2</sub>-microglobulin (β<sub>2</sub>-m)

### IVD

Store 2 - 8°C.

### PRINCIPLE OF THE METHOD

The β<sub>2</sub>-m Turbilatex is a quantitative turbidimetric test for the measurement of β<sub>2</sub>-microglobulin (β<sub>2</sub>-m) in human serum, plasma or urine. Latex particles coated with anti-human β<sub>2</sub>-m are agglutinated when mixed with samples containing β<sub>2</sub>-m. The agglutination causes an absorbance change, dependent upon the β<sub>2</sub>-m contents of the patient sample that can be quantified by comparison from a calibrator of known concentration.

### CLINICAL SIGNIFICANCE

β<sub>2</sub>-m is a protein located on the surface of human lymphocytes and other nucleated cells. Free β<sub>2</sub>-m is filtered by the glomerulus and subsequently reabsorbed in the proximal tubular cells. Increased urinary excretion of β<sub>2</sub>-m is a sensitive indicator of renal insufficiency. Also, the β<sub>2</sub>-m level in serum is a useful marker of other diseases including carcinomas, lymphoid tumors, rheumatoid arthritis and AIDS.

### REAGENTS

<b>β<sub>2</sub>-m Diluent (R1)</b>	Tris buffer 20 mmol/L, pH 8.2. Preservative.
<b>β<sub>2</sub>-m Latex (R2)</b>	Particles coated with goat IgG anti-human β <sub>2</sub> -m, pH 8.0. Preservative.
<b>β<sub>2</sub>-m CAL</b>	Calibrator. β <sub>2</sub> -m concentration is stated on the vial label.
<b>Optional:</b>	Ref: 1107034 β <sub>2</sub> -m Control.

### PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

### CALIBRATION

Use β<sub>2</sub>- Microglobulin Calibrator Reference 1107032. The sensitivity of the assay and the target value of the calibrator have been standardized against the 1<sup>st</sup> International ERM-DA470k/IFCC. Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

### PREPARATION

Ready for use.

#### β<sub>2</sub>-m Calibrator:

Serum method: Reconstitute (→) with 1.0 mL of distilled water. Mix gently and bring to room temperature for about 10 minutes before use.

Urine method: Dilute reconstituted calibrator 1/6 with NaCl 9 g/L (50 μL calibrator + 250 μL NaCl 9 g/L).

### STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

**Reagent deterioration:** Presence of particles (R1, R2) and turbidity (R1).

**β<sub>2</sub>-m Calibrator:** Stable for 1 month at 2-8°C or 3 months at -20°C.

### ADDITIONAL EQUIPMENT

- SPINLAB 180 autoanalyzer
- Laboratory equipment.

### SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

Fresh urine. Adjust samples to pH 7-8 by the addition of K<sub>2</sub>HPO<sub>4</sub>. Stable 2 days at 2-8°C or 2 months at -20°C.

The samples with particles or fibrin should be centrifuged before testing. Do not use hemolyzed or lipemic samples.

### QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT β<sub>2</sub>-m Control (Ref: 1107034). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

### REFERENCE VALUES

**Serum:** from 1.0 to 3.0 mg/L.

**Urine:** from 0.1 to 0.3 mg/L.

Each laboratory should establish its own reference range

## SPINLAB 180 APPLICATION

Name	b2M	Ref. male low	*
Abbr. Name	b2M	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	mg/L	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	0.22 mg/L	Panic value low	*
High Conc.	18.00 mg/L	Panic value high	*
Calibrator name	CAL b2M	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 μL		
rerun volume	240 μL		
Sample			
normal volume	3.0 μL		
rerun volume	3.0 μL		
R2 bottle (mL)	5 mL		
normal volume	60 μL		
rerun volume	60 μL		
Predilución	No		
Slope blank	No		
Point one,two	6, 183 sec.		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

### INTERFERENCES

**Serum method:** bilirubin (20 mg/L), hemoglobin (10 g/L) and lipids (10 g/L), do not interfere. Rheumatoid factors (150 IU/mL), interfere.

**Urine method:** urea (urine) (50 g/L), uric ac. (20 g/L) and glucose (100 g/L), do not interfere.

Other substances may interfere<sup>7</sup>.

### PERFORMANCE CHARACTERISTICS

- Linearity limit:** Up to 18 mg/L (serum) and 3 mg/L (urine), under the described assay conditions. Samples higher results should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity depends on the sample-reagent ratio, as well as the analyzer used. It will be higher by decreasing sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection limit:** Values less than 0.09 mg/L (serum) and 0.01 mg/L (urine) give non-reproducible results.
- Prozone effect:** No prozone effect was detected up to 100 mg/L (serum) and 20 mg/L (urine).
- Sensitivity:** Δ 0.048 A. mg/L (serum) and Δ 0.228 A. mg/L (urine).
- Precision:** The reagent has been tested for 20 days, using three different β<sub>2</sub>-m concentrations in a EP5-based study.

EP5	CV (%)		
	+/- 1 mg/L	+/- 3.2 mg/L	+/- 8.5 mg/L
Total	4.0%	3.4%	1.7%
Within Run	2.8%	2.0%	1.2%
Between Run	1.7%	1.5%	1.2%
Between Day	2.2%	2.4%	0.0%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 36 samples of different concentrations of β<sub>2</sub>-m were assayed. The correlation coefficient (r) was 0.97 and the regression equation y = 1.709x - 2.627. The results of the performance characteristics depend on the analyzer used.

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Name	b2M	Ref. male low	*
Abbr. Name	b2M	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	mg/L	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	0,22 mg/L	Panic value low	*
High Conc.	18,00 mg/L	Panic value high	*
Calibrator name	CAL b2M	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 µL		
rerun volume	240 µL		
Sample			
normal volume	3,0 µL		
rerun volume	3,0 µL		
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60 µL		
Predilución	No		
Slope blank	No		
Point one,two	6, 183 sec.		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

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**BIBLIOGRAFIA**

- Bhalla, R.B. et al. Clinical Chemistry 1983; 29: 1560.
- Malaguarnera M et al. Digestive Diseases and Sciences 1997; 42: 762-766
- Chironna et al. Int J Clin Lab Rws 1994; 24: 90-93
- Wibell L et al. Nephron 1973; 10: 320-331
- Berggard B et al. Scand J Clin Lab Invest 1980; 40: 13-25.
- Davey P G et al. Clin Chem 1982; 28/6: 1330-1333
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AAC Press, 1995.

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Ref.: SP1107030

Cont.

R1 Diluente : 2 x 20 mL

R2 Latex : 2 x 5 mL

 β<sub>2</sub>-m CAL : 1 x 1 mL