β₂-m TURBI

β₂-m turbilatex

Latex turbidimetry

Quantitative determination of β₂-microglobulin (β₂-m)

IVD

Store 2 - 8°C.

PRINCIPLE OF THE METHOD

The β₂-m Turbilatex is a quantitative turbidimetric test for the measurement of β₂-microglobulin (β₂-m) in human serum, plasma or urine. Latex particles coated with anti-human β₂-m are agglutinated when mixed with samples containing β₂-m. The agglutination causes an absorbance change, dependent upon the β₂-m contents of the patient sample that can be quantified by comparison from a calibrator of known concentration.

CLINICAL SIGNIFICANCE

β₂-m is a protein located on the surface of human lymphocytes and other nucleated cells. Free β₂-m is filtered by the glomerulus and subsequently reabsorbed in the proximal tubular cells. Increased urinary excretion of β₂-m is a sensitive indicator of renal insufficiency. Also, the β₂-m level in serum is a useful marker of other diseases including carcinomas, lymphoid tumors, rheumatoid arthritis and AIDS.

REAGENTS

β₂-m Diluent (R1)	Tris buffer 20 mmol/L, pH 8.2. Preservative.
β₂-m Latex (R2)	Particles coated with goat IgG anti-human β ₂ -m, pH 8.0. Preservative.
β₂-m CAL	Calibrator. β ₂ -m concentration is stated on the vial label.
Optional:	Ref: 1107034 β ₂ -m Control.

PRECAUTIONS

R1, R2 and β₂-m CAL: H317 - May cause an allergic skin reaction. Contains 2-Methylisothiazol-3(2H)-one (Proclin 950). Follow the precautionary statements given in MSDS and label of the product. Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

CALIBRATION

Use β₂- Microglobulin Calibrator Reference 1107032. The sensitivity of the assay and the target value of the calibrator have been standardized against the 1st International ERM-DA470k/IFCC. Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

PREPARATION

Ready for use.
β₂-m Calibrator:
Serum method: Reconstitute (→) with 1.0 mL of distilled water. Mix gently and bring to room temperature for about 10 minutes before use.
Urine method: Dilute reconstituted calibrator 1/6 with NaCl 9 g/L (50 µL calibrator + 250 µL NaCl 9 g/L).

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date. Do not freeze; frozen Latex or Diluent could change the functionality of the test.
Reagent deterioration: Presence of particles (R1, R2) and turbidity (R1).
β₂-m Calibrator: Stable for 1 month at 2-8°C or 3 months at -20°C.

ADDITIONAL EQUIPMENT

- SPINLAB 180 autoanalyzer
- Laboratory equipment.

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.
Fresh urine. Adjust samples to pH 7-8 by the addition of K₂HPO₄. Stable 2 days at 2-8°C or 2 months at -20°C.
The samples with particles or fibrin should be centrifuged before testing. Do not use hemolyzed or lipemic samples.

QUALITY CONTROL

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT β₂-m Control (Ref: 1107034). Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

Serum: from 1.0 to 3.0 mg/L.
Urine: from 0.1 to 0.3 mg/L.
Each laboratory should establish its own reference range

SPINLAB 180 APPLICATION

Name	b2M	Ref. male low	*
Abbr. Name	b2M	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	mg/L	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	0.22 mg/L	Panic value low	*
High Conc.	18.00 mg/L	Panic value high	*
Calibrator name	CAL b2M	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 µL		
rerun volume	240 µL		
Sample			
normal volume	3.0 µL		
rerun volume	3.0 µL		
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60 µL		
Predilución	No		
Slope blank	No		
Point one,two	6, 183 sec.		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

INTERFERENCES

Serum method: bilirubin (20 mg/L), hemoglobin (10 g/L) and lipids (10 g/L), do not interfere. Rheumatoid factors (150 IU/mL), interfere.
Urine method: urea (urine) (50 g/L), uric ac. (20 g/L) and glucose (100 g/L), do not interfere.
Other substances may interfere⁷.

PERFORMANCE CHARACTERISTICS

- Linearity limit:** Up to 18 mg/L (serum) and 3 mg/L (urine), under the described assay conditions. Samples higher results should be diluted 1/5 in NaCl 9 g/L and retested again. The linearity depends on the sample-reagent ratio, as well as the analyzer used. It will be higher by decreasing sample volume, although the sensitivity of the test will be proportionally decreased.
- Detection limit:** Values less than 0.09 mg/L (serum) and 0.01 mg/L (urine) give non-reproducible results.
- Prozone effect:** No prozone effect was detected up to 100 mg/L (serum) and 20 mg/L (urine).
- Sensitivity:** Δ 0.048 A. mg/L (serum) and Δ 0.228 A. mg/L (urine).
- Precision:** The reagent has been tested for 20 days, using three different β₂-m concentrations in a EP5-based study.

EP5	CV (%)		
	+/- 1 mg/L	+/- 3.2 mg/L	+/- 8.5 mg/L
Total	4.0%	3.4%	1.7%
Within Run	2.8%	2.0%	1.2%
Between Run	1.7%	1.5%	1.2%
Between Day	2.2%	2.4%	0.0%

- Accuracy:** Results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 36 samples of different concentrations of β₂-m were assayed. The correlation coefficient (r) was 0.97 and the regression equation y = 1.709x - 2.627. The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

- Bhalla, R.B. et al. Clinical Chemistry 1983; 29: 1560.
- Malaguarrera M et al. Digestive Diseases and Sciences 1997; 42: 762-766.
- Chironna et al. Int J Clin Lab Rws 1994; 24: 90-93.
- Wibell L et al. Nephron 1973; 10: 320-331.
- Berggard B et al. Scand J Clin Lab Invest 1980; 40: 13-25.
- Davey P G et al. Clin Chem 1982; 28/6: 1330-1333.
- Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Pres, 1995.

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APPLICATION AU SPINLAB 180

Name	b2M	Ref. male low	*
Abbr. Name	b2M	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	mg/L	Ref. Ped. Low	*
Decimals	2	Ref. Ped. High	*
Low Conc.	0.22 mg/L	Panic value low	*
High Conc.	18.00 mg/L	Panic value high	*
Calibrator name	CAL b2M	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000

DUAL MODE			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	240 µL		
rerun volume	240 µL		
Sample			
normal volume	3.0 µL		
rerun volume	3.0 µL		
R2 bottle (mL)	5 mL		
normal volume	60 µL		
rerun volume	60 µL		
Predilución	No		
Slope blank	No		
Point one,two	6, 183 sec.		

Reagent blank	No
Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs
Substr.Depletion	3.000 Abs

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BIBLIOGRAPHIE

- Bhalla, R.B. et al. Clinical Chemistry 1983; 29: 1560.
- Malaguarnera M et al. Digestive Diseases and Sciences 1997; 42: 762-766.
- Chironna et al. Int J Clin Lab Rws 1994; 24: 90-93.
- Wibell L et al. Nephron 1973; 10: 320-331.
- Berggard B et al. Scand J Clin Lab Invest 1980; 40: 13-25.
- Davey P G et al. Clin Chem 1982; 28/6: 1330-1333.
- Young DS. Effects of drugs on clinical laboratory test, 4thed. AACCPres, 1995.

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