**Quantitative determination of Ferritin  
IVD**

Store 2 - 8°C.

**PRINCIPLE OF THE METHOD**

Ferritin-turbilatex is a quantitative turbidimetric test for the measurement of ferritin in human serum or plasma.

Latex particles coated with specific anti-human ferritin are agglutinated when mixed with samples containing ferritin. The agglutination causes an absorbance change, dependent upon the ferritin contents of the sample that can be quantified by comparison from a calibrator of known ferritin concentration.

**CLINICAL SIGNIFICANCE**

Serum ferritin concentration usually reflects body iron stores and is considered one of the most reliable indicators of iron status of patients. Whereas low serum concentrations of ferritin are always indicative of an iron deficiency, elevated concentrations can occur for variety of reasons. Thus, although elevated concentrations often indicate an excessive iron intake, they are also caused by liver disease, chronic inflammation and malignancies. Pregnant women, blood donors, hemodialysis patients, adolescents and children are groups particularly at risk.

**REAGENTS**

<b>Diluent (R1)</b>	Tris buffer 20 mmol/L, pH 8.2. Preservative.
<b>Latex (R2)</b>	Latex particles coated with rabbit IgG anti-human ferritin, pH 8.2. Preservative.
<b>FERR-CAL</b>	Calibrator. Ferritin concentration is stated on the vial label.
<b>Optional</b>	Ref.: 1107044 Ferritin control.

**PRECAUTIONS**

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However, handle cautiously as potentially infectious.

**CALIBRATION**

Use Ferritin Calibrator Reference 1107042.

The sensitivity of the assay and the target value of the calibrator have been standardized against the 3rd International Standard of Ferritin (94/572, 2008 WHO). Recalibrate when control results are out of specified tolerances, when using different lot of reagent and when the instrument is adjusted.

**PREPARATION**

**Ferritin Calibrator:** Reconstitute (→) with 3.0 mL of distilled water. Mix gently and incubate at room temperature for about 10 minutes before testing.

**Calibration curve:** Prepare the following FERR calibrator dilutions in NaCl 9 g/L. Multiply the concentration of the FERR calibrator by the corresponding factor stated in table below to obtain the concentration of each dilution.

Calibrator dilution	1	2	3	4	5	6
Calibrator FERR (µL)	--	25	50	100	200	400
NaCl 9 g/L (µL)	400	375	350	300	200	--
Dilution Factor	0	0,0625	0,125	0,25	0,5	1,0

**STORAGE AND STABILITY**

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Reagents should not be left inside the analyzer after use, they must be stored refrigerated at 2-8°C. Latex may sediment. Mix reagents gently before use. Do not use reagents over the expiration date.

Do not freeze; frozen Latex or Diluent could change the functionality of the test.

**Reagent deterioration:** Presence of particles (R1, R2) and turbidity (R1).

**ADDITIONAL EQUIPMENT**

- SPINLAB 180 autoanalyzer
- Laboratory equipment.

**SAMPLES**

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

**QUALITY CONTROL**

Control Sera are recommended to monitor the performance of manual and automated assay procedures. It should be used the SPINREACT Ferritin Control (Ref.: 1107044).

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

**SPINLAB 180 APPLICATION**

Name	Ferritin	Ref. male low	*
Abbr. Name	FERR	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	µg/L	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	0 µg/L	Panic value low	*
High Conc.	400 µg/L	Panic value high	*
Calibrator name	CAL FERR (6 pts)	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000

**DUAL MODE**

Sample blank	No
R1 bottle (mL)	25 mL
normal volume	180 µL
rerun volume	180 µL
Sample	
normal volume	20.0 µL
5.0 µL	
rerun volume	20.0 µL
3.0 µL	
R2 bottle (mL)	5 mL
normal volume	45 µL
rerun volume	45 µL
Predilución	No
Slope blank	No
Point one,two	6, 236 sec.

Reagent blank	No
Low Absorbance	-0.100 Abs
High Absorbance	3.000 Abs
R. Abs. L. Limit	-0.100 Abs
R. Abs. H. Limit	3.000 Abs
Substr.Depletion	3.000 Abs

**REFERENCE VALUES**

Men: 30 – 220 µg/L.

Women: 20 – 110 µg/L.

Each laboratory should establish its own eference range.

**PERFORMANCE CHARACTERISTICS**

**Measuring range:** Up to 600 µg/L. Samples with higher values should be diluted 1/5 in NaCl 9 g/L and retested. The upper linearity limit increases as the sample volume and the sensitivity decrease.

**Detection limit:** 5,04 µg/L.

**Quantification limit:** Values under 6,6 µg/L may give non-reproducible results.

**Prozone effect:** No prozone effect was detected at least up to 9000 µg/L.

**Precision:** According to the EP5-A2 standards (CLSI), the reagent has been tested for 20 days, measuring each level per duplicate twice a day (n=80):

Mean (µg/L)	Intra-assay (n= 80)			Total (n= 80)		
	33,4	114,5	289,8	33,4	114,5	289,8
SD	1,7	1,4	2,4	2,1	3,4	7,5
CV (%)	5,1	1,2	0,8	6,3	2,9	2,6

**Method comparison:** The reagent was compared to another commercially available Ferritin reagent by testing 144 samples (male and female), with concentrations between 6,97 and 730 µg/L. The coefficient of correlation (r) was 0,988, and the equation  $y = 0,96x + 1,15$

Performance characteristics depend on the analyzer used.

**NOTES**

Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

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**APPLICATION AU SPINLAB 180**

Name	Ferritin	Ref. male low	*
Abbr. Name	FERR	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	µg/L	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	0 µg/L	Panic value low	*
High Conc.	400 µg/L	Panic value high	*
Calibrator name	CAL FERR (6 pts)	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1,000
		Correlat. offset	0,000
<b>DUAL MODE</b>			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	180 µL		
rerun volume	180 µL		
Sample			
normal volume	20,0 µL		
5,0 µL			
rerun volume	20,0 µL		
3,0 µL			
R2 bottle (mL)	5 mL		
normal volume	45 µL		
rerun volume	45 µL		
Predilución	No		
Slope blank	No		
Point one,two	6, 236 sec.		
Reagent blank	No		
Low Absorbance	-0,100 Abs		
High Absorbance	3,000 Abs		
R. Abs. L. Limit	-0,100 Abs		
R. Abs. H. Limit	3,000 Abs		
Substr.Depletion	3,000 Abs		

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Name	Ferritin	Ref. male low	*
Abbr. Name	FERR	Ref. male high	*
Mode	Twopoints	Ref. female low	*
Wavelength	546 nm	Ref. female high	*
Units	µg/L	Ref. Ped. Low	*
Decimals	1	Ref. Ped. High	*
Low Conc.	0 µg/L	Panic value low	*
High Conc.	400 µg/L	Panic value high	*
Calibrator name	CAL FERR (6 pts)	Control 1	*
Prozone check	No	Control 2	*
		Control 3	*
		Correlat. factor	1.000
		Correlat. offset	0.000
<b>DUAL MODE</b>			
Sample blank	No		
R1 bottle (mL)	25 mL		
normal volume	180 µL		
rerun volume	180 µL		
Sample			
normal volume	20.0 µL		
5.0 µL			
rerun volume	20.0 µL		
3.0 µL			
R2 bottle (mL)	5 mL		
normal volume	45 µL		
rerun volume	45 µL		
Predilución	No		
Slope blank	No		
Point one,two	6, 236 sec.		
Reagent blank	No		
Low Absorbance	-0.100 Abs		
High Absorbance	3.000 Abs		
R. Abs. L. Limit	-0.100 Abs		
R. Abs. H. Limit	3.000 Abs		
Substr.Depletion	3.000 Abs		

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