
Quantitative determination of lipase
IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

The pancreatic lipase in presence of colipase, desoxycholate and calcium ions, hydrolyses the substrate 1-2-O-dilauryl-rac-glycero-3-glutaric acid-(6' -methylresorufin)-ester. The sequence of reactions involved in the enzymatic direct lipase determination is the following:

The rate of methylresorufin formation, measured photometrically, is proportional to the catalytic concentration of lipase present in the sample.

CLINICAL SIGNIFICANCE

Lipase (LPS) is a pancreatic enzyme necessary for the absorption and digestion of nutrients that catalyzes the hydrolysis of glycerol esters of fatty acids. Determination of LPS is used for diagnosis of diseases of pancreas such as acute and chronic pancreatitis and obstruction of the pancreatic duct^{1,7,8}. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	TRIS pH 8.3	40 mmol/L
Buffer	Colipase	≥ 1 mg/L
	Desoxycholate	1.8 mmol/L
	Taurodesoxycholate	7.2 mmol/L
R 2	Tartrate pH 4,0	15 mmol/L
Substrate (micro-emulsion)	Lipase Substrate	≥ 0.7 mmol/L
	Calcium chloride (CaCl ₂)	0.1 mmol/L
OPTIONAL	SPINTROL H CAL	

PREPARATION

- R1 – R2 Ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 600 nm ≥ 1.00.
- R 2 is a turbid orange-colored micro-emulsion, discard if turning to red.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment ^(Note 2).

SAMPLES

Serum or plasma with sodium citrate, EDTA or heparin¹.
Avoid repeated frozen and unfrozen. Stability: 2 days at 2-8°C.

REFERENCE VALUES¹

≤ 38 U/L (U/L methylresorufin at 37°C).

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera and calibrators are recommended to monitor the performance of assay procedures: SPINTROL H Calibrator, SPINTROL H Normal and Pathologic (Ref. 1002011, 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINTECH 240

Item Name LIPASE						
<u>DATA INFORMATION</u>				<u>CALIBRATION</u>		
Units	U/L	TYPE	Linear			
Decimals	1					
<u>ANALYSIS</u>				<u>STANDARD</u>		
Type	RATE	#1	*	#4		
		#2		#5		
		#3		#6		
W.Length 1	600	<u>NORMAL RANGE (37°C)</u>		LOW	HIGH	
Method	Kinetic	SERUM				
		FEMALE				
CORR SLOPE	INTER					
1.000 x +	0					
Item Name LIPASE						
<u>ASPIRATION</u>				<u>DATA PROCESS</u>		
KIND	Single <input type="checkbox"/> Double <input checked="" type="checkbox"/>	<u>READ</u>		<u>ABSORBANCE LIMIT</u>		
		START	END	LOW	-3.000	
		MAIN	37	45	HIGH	3.000
		SUB				
SAMPLE	VOLUME			<u>ENDPOINT LIMIT</u>		
REAGENT 1	300 µL			LINEAR CHECK (%) 90		
REAGENT 2	60 µL					
				<u>FACTOR</u>		
Third Mix	<input checked="" type="checkbox"/> OFF <input type="checkbox"/> ON	Blank Correction		1.000		
R1 Blank	Water <input checked="" type="checkbox"/> R1-B					
<u>MONITOR</u>		<u>PROZONE CHECK</u>				
O LEVEL POINT	1	START	END	LIMIT (%)		
SPAN	3.000	FIRST		<input checked="" type="checkbox"/> Low High		
		SECOND		<input checked="" type="checkbox"/> Low High		
		THIRD				

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for 3 days.

PERFORMANCE CHARACTERISTICS

Measuring range: From detection limit of 5 U/L to linearity limit of 250 U/L.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	CV (%)	
Mean (U/L)	40,2	59,35	38,5	58,9
SD	0,410	0,875	1,10	1,25
CV (%)	1,02	1,47	2,86	2,13

Sensitivity: 1 U/L = 0,00059792 (A)

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 101 samples were the following:

Correlation coefficient (r)²: 0,99732.

Regression equation: y = 0,50054x + 3,9443.

The results of the performance characteristics depend on the analyzer used.

NOTES

1. In some storage conditions (i.e. storage at a temperature lower than the one indicated) a precipitate may appear in the vial that will not influence that the reagent performance; however, it is recommended to resuspend the product with a slight rotation.
2. In order to avoid contamination, it is recommended to use disposable material.

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