

Quantitative determination of creatine kinase liquid (CK)
IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Kinetic determination of the creatine kinase based upon IFCC and DGKC recommendations.

Creatine kinase (CK) catalyses the reversible transfer of a phosphate group from phosphocreatine to ADP. This reaction is coupled to those catalysed by hexokinase (HK) and glucose-6-phosphate dehydrogenase (G6P-DH):

 The rate of NADPH formation, measured photometrically, is proportional to the catalytic concentration of CK present in the sample^{1,2}.

CLINICAL SIGNIFICANCE

Creatine kinase is a cellular enzyme with wide tissue distribution in the body. Its physiological role is associated with adenosine triphosphate (ATP) generation for contractile or transport systems.

 Elevated CK values are observed in diseases of skeletal muscle and after myocardial infarction^{1,5,6}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	Imidazol, pH 6.7	125 mmol/L
	D-Glucose	25 mmol/L
	N-Acetyl-L-Cysteine	25 mmol/L
	Magnesium acetate	12.5 mmol/L
	NADP	2.52 mmol/L
	EDTA	2.02 mmol/L
	Hexokinase	≥6 800 U/L
R 2	ADP	15.2 mmol/L
	AMP	25 mmol/L
	di-Adenosine-5- penta-P	103 mmol/L
	Glucose-6-phosphate DH	≥8 800 U/L
	Creatine phosphate	250 mmol/L

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use the tablets if appears broken.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm ≥ 1.60.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

SAMPLES

Serum free of hemolysis or heparin plasma.

Stability 7 days at 2-8°C, protected from light.

The creatin kinase activity decreases 10% after 1 day at 2-5°C or after 1 hour at 15-25°C.

Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1.00	1.56	2.44
30°C	0.64	1.00	1.56
37°C	0.41	0.63	1.00

REFERENCE VALUES¹

	25°C	30°C	37°C
Men, up to	80 U/L	130 U/L	195 U/L
Women, up to	70 U/L	110 U/L	170 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

APPLICATION SPINTECH 240

Item Name CKNAC					
DATA INFORMATION					
Units	U/L			CALIBRATION	
Decimals	0			TYPE	Linear
ANALYSIS					
Type	RATE			STANDARD	
W.Length 1	340			#1	* #4
				#2	#5
				#3	#6
NORMAL RANGE					
Method	Kinetic UV.Liquid			SERUM	MALE LOW HIGH
				FEMALE	24 195 24 170
CORR SLOPE	INTER				
1.000 x +	0				
Item Name CKNAC					
ASPIRATION					
KIND	Single	✓ Double			
		VOLUME			
SAMPLE	5	μL			
REAGENT 1	240	μL			
REAGENT 2	60	μL			
		Third Mix	✓ OFF	ON	
		R1 Blank	✓ Water	R1-B	
MONITOR					
0 LEVEL POINT	1			DATA PROCESS	
SPAN	3.000			START	END
				MAIN	40 52
				SUB	
				ABSORBANCE LIMIT	
					LOW -0.1
					HIGH 3.000
				ENDPOINT LIMIT	
				LINEAR CHECK (%) 90	
FACTOR					
				Blank Correction	1.000
PROZONE CHECK					
				START	END
				FIRST	
				SECOND	✓ Low High
				THIRD	✓ Low High

Blank parameter must be performed in order to get good results in CALIB screen from main menu. This parameter calibration is stable for more than 7days.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

PERFORMANCE CHARACTERISTICS
Measuring range: From detection limit of 2,12 U/L to linearity limit of 2000 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

Mean (U/L)	Intra-assay		Inter-assay	
	147	494	145	485
SD	1,23	3,60	2,91	8,97
CV (%)	0,84	0,73	2,01	1,85

Sensitivity: 1 U/L = 0,00012 ΔA/min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained were the following:

 Correlation coefficient (r)²: 0,9995

Regression equation: y= 1,0846x – 0,3512.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Reccommandation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la créatinine kinase dans le sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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 Ref: TK41250 Cont. R1: 10 x 25 mL, R2: 10 x 7 mL
 Ref: TK41251 Cont. R1: 2 x 25 mL, R2: 2 x 7 mL

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BIBLIOGRAPHIE

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la créatinine kinase dans le sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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BIBLIOGRAFIA

1. Abbot B et al. Creatinine kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-1116.
2. Gerhardt W et al. Creatine kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
4. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.
7. Mathieu M. et coll. Recommendation pour la mesure de la concentration catalytique de la créatinine kinase dans la sérum humain. Ann. Biol. Clin.,40, (1482), 87.

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 Ref: TK41250 R1: 10 x 25 mL, R2: 10 x 7 mL
 Ref: TK41251 Cont. R1: 2 x 25 mL, R2: 2 x 7 mL