
Quantitative determination of lactate dehydrogenase (LDH) IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Lactate dehydrogenase (LDH) catalyses the reduction of pyruvate by NADH, according the following reaction:

 The rate of decrease in concentration of NADPH, measured photometrically, is proportional to the catalytic concentration of LDH present in the sample¹.

CLINICAL SIGNIFICANCE

Lactate dehydrogenase (LDH) is an enzyme with wide tissue distribution in the body.

The higher concentrations of LDH are found in liver, heart, kidney, skeletal muscle and erythrocytes.

 Increased levels of the enzyme are found in serum in liver disease, myocardial infarction, renal disease, muscular dystrophy and anemia^{1,4,5}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

Reagent 1	Imidazol	65 mmol/L
Buffer	Pyruvate	0.6 mmol/L
Reagent 2	NADH	0.18 mmol/L
Substrate		

PREPARATION

All the reagents are ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 340 nm <1.00.

ADDITIONAL EQUIPMENT

- Autoanalyzer Spintech 240.
- General laboratory equipment.

SAMPLES

 Serum¹. Separated from cells as rapidly as possible. Do not use oxalates as anticoagulants since they inhibit the enzyme.

Do not use haemolysed samples. Stability: 2 days at 2-8°C.

Temperature conversion factors

To correct results to other temperatures multiply by:

Assay temperature	Conversion factor to		
	25°C	30°C	37°C
25°C	1.00	1.33	1.92
30°C	0.75	1.00	1.43
37°C	0.52	0.70	1.00

REFERENCE VALUES¹

25°C	30°C	37°C
120-240 U/L	160-320 U/L	230-460 U/L

These values are for orientation purpose; each laboratory should establish its own reference range.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and technique for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

APPLICATION SPINTECH 240

Item Name LDH DATA INFORMATION Units U/L Decimals 0 ANALYSIS Type RATE W.Length 1 340 Method DGKC CORR SLOPE INTER 1.000 x + 0		CALIBRATION TYPE Linear STANDARD #1 * #4 #2 #5 #3 #6 NORMAL RANGE (37°C) SERUM MALE 230 460 FEMALE 230 460	
Item Name LDH ASPIRATION KIND Single <input checked="" type="checkbox"/> Double VOLUME SAMPLE 5 µL REAGENT 1 240 µL REAGENT 2 60 µL Third Mix <input checked="" type="checkbox"/> OFF ON R1 Blank <input checked="" type="checkbox"/> Water R1-B		DATA PROCESS ABSORBANCE LIMIT READ LOW -3.000 START END HIGH 3.000 MAIN 36 48 SUB ENDPOINT LIMIT LINEAR CHECK (%) 90 FACTOR Blank Correction 1.000 PROZONE CHECK START END LIMIT (%) FIRST SECOND <input checked="" type="checkbox"/> Low High THIRD <input checked="" type="checkbox"/> Low High	
MONITOR 0 LEVEL POINT 1 SPAN 3.000			

Blank parameter must be performed in order to get good results in CALIB screen from main menu.

PERFORMANCE CHARACTERISTICS
Measuring range: From detection limit of 3,42 U/L to linearity limit of 1600 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

Precision:

	Intra-assay (n=20)		Inter-assay (n=20)	
	Mean (U/L)	SD	392	773
Mean (U/L)	400	785	6,23	9,93
SD	3,15	10,97	1,59	1,28
CV (%)	0,79	1,40		

Sensitivity: 1 U/L = 0,00009 ΔA/min.

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

 Correlation coefficient (r)²: 0,98382.

Regression equation: y= 0,8988x + 2,583.

The results of the performance characteristics depend on the analyzer used.

BIBLIOGRAPHY

1. Pesce A. Lactate dehydrogenase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1124-117, 438.
2. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
3. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
4. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
5. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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